

Determination of free cyanide by polarography

Summary

This Application Bulletin describes a polarographic method for the determination of cyanide that allows to determine free cyanide fast and accurately. The determination also succeeds in solutions containing sulfides, where other methods fail. Cyanide concentrations in the range $\beta(\text{CN}^-) = 0.01 \dots 10 \text{ mg/L}$ cause no problems. Interference caused by anions and complexed cyanides has been investigated.

Instruments

VA instrument
capable of operating a Multi-Mode Electrode and supporting differential pulse (DP) measuring mode

Electrodes

WE	Multi-Mode Electrode pro	6.1246.120
	Mercury drop capillary	6.1226.030
RE	Ag/AgCl reference electrode	6.0728.x20
	Ag/AgCl/KCl (3 mol/L)	
	Electrolyte vessel	6.1245.010
	Filled with $c(\text{KCl}) = 3 \text{ mol/L}$	
AE	Pt rod electrode	6.0343.x00

Reagents

All of the used reagents must be of purest quality possible (for analysis or for trace analysis*).

- Potassium hydroxide, $w(\text{KOH}) \geq 85\%$, for analysis, CAS 1310-58-3
- Boric acid, H_3BO_3 , for analysis, CAS 10043-35-3
- Potassium cyanide, KCN, for analysis, CAS 151-50-8
- Ultrapure water, resistivity $> 18 \text{ M}\Omega \cdot \text{cm}$ (25 °C), type I grade (ASTM D1193)

Solutions

Supporting electrolyte $c(\text{H}_3\text{BO}_3) = 0.2 \text{ mol/L}$
 $c(\text{KOH}) = 0.17 \text{ mol/L}$
 pH 10.2
 Dissolve 11.2 g KOH in approx. 800 mL ultrapure water. Add

12.4 g boric acid, adjust the pH value to 10.2 if necessary and make up to 1 L with ultrapure water.

Standard solutions

Potassium hydroxide solution $c(\text{KOH}) = 0.01 \text{ mol/L}$
 Dissolve 0.65 g KOH in ultrapure water and make up to 1 L.

Cyanide stock standard solution $\beta(\text{CN}^-) = 1 \text{ g/L}$
 Dissolve 0.2503 g KCN in $c(\text{KOH}) = 0.01 \text{ mol/L}$ and make up to 100 mL.

Cyanide standard solution $\beta(\text{CN}^-) = 100 \text{ mg/L}$
 Additional standards can be prepared by dilution with $c(\text{KOH}) = 0.01 \text{ mol/L}$.

Analysis

Measuring solution

10 mL (diluted) sample solution
 + 10 mL supporting electrolyte

If cyanide is present in concentrations above the linear working range, the sample solution has to be diluted accordingly with ultrapure water.

The concentration is determined by standard addition.


Parameters

Voltammetric


Electrode operating mode	DME
Measuring mode	DP – Differential pulse
Stirring rate	2000 min^{-1}
Equilibration time	5 s
Sweep	
Start potential	0.0 V
End potential	-0.5 V
Potential step	0.008 V
Potential step time	0.8 s

Sweep rate	0.01 V/s
Pulse amplitude	0.05 V
Substance	
Name	Cyanide
Characteristic potential	-0.24 V

Example



CN
 $c = 457.939 \text{ } \mu\text{g/L}$
 $+/- 9.446 \text{ } \mu\text{g/L} (2.06\%)$



Result

Sample	Waste water
Sample size	10.0 mL
$\beta(\text{CN}^-)$	458 $\mu\text{g/L}$

Comments

- Cyanide can be determined in the presence of a 1000-fold excess of phosphate, nitrate and sulfate. A 50'000-

fold excess of chloride does not interfere with the determination either.

- $\text{K}_3[\text{Fe}(\text{CN})_6]$, $\text{K}_4[\text{Fe}(\text{CN})_6]$ and $\text{K}_2[\text{Ni}(\text{CN})_4]$ do not affect the height of the cyanide peak; $\text{K}[\text{Zn}(\text{CN})_3]$ increases the peak height (decomposition of the complex); $\text{K}[\text{Cu}(\text{CN})_2]$ somewhat distorts the shape of the peak when present in tenfold excess.
- If the CN^- concentration in the sample is greater than 10 mg/L, a smaller sample size has to be used. Concentrations as low as $\beta(\text{CN}^-) = 0.01 \text{ mg/L}$ can still be determined with a precision of $\pm 10\%$.
- It is also possible to determine cyanide in the presence of sulfide (e.g. in waste water distillates).
- With the method described only free cyanide can be determined. For the determination of the total cyanide content the cyanide has to be separated by distillation, e.g. according to DIN 38405-13.
- The polarographic method for determining cyanide is quicker, sometimes more sensitive, and, above all, simpler than comparable methods. The times listed below for all the methods refer to a single determination (without standard solution).

Method	Duration
Titration $\beta(\text{CN}^-) = 1 \dots 100 \text{ mg/L}$	5...30 min
Direct potentiometry with the ion-selective electrode $\beta(\text{CN}^-) = 0.26 \dots 26 \text{ mg/L}$	up to 45 min (depending on conc.)
Colorimetry (benzidine/pyridine method) $\beta(\text{CN}^-) = 0.005 \dots 0.1 \text{ mg/L}$	approx. 15 min
Colorimetry (barbituric acid/pyridine method) $\beta(\text{CN}^-) = 0.005 \dots 0.5 \text{ mg/L}$	approx. 40 min
Polarography $\beta(\text{CN}^-) = 0.01 \dots 10 \text{ mg/L}$	approx. 5 min

References

- D. R. Canterford
Simultaneous determination of cyanide and sulfide with rapid direct current polarography
Anal. Chem. 47 (1975) 88–92

Appendix

Report for the example determination of free cyanide in waste water

===== METROHM 797 VA COMPUTRACE (Version 1.0.0.1) (Serial No. 0) =====

Determination : 06261315_waste water.dth
Sample ID : waste water
Creator method : Date : Time:
Creator determ.: Date : 2001-06-26 Time: 13:15:36
Modified by : Date : 2001-06-26 Time: 16:08:59

Method : AB110_Det of CN.mth
Title : Determination of free Cyanide
Remark1 : 10 mL electrolyte + 10 mL sample
Remark2 :

Sample amount : 10.000 mL
Cell volume : 20.000 mL

Substance : CN
Conc. : 228.970 ug/L
Conc.dev. : 4.723 ug/L (2.06%)
Amount : 4.579 ug
Add.amount : 5.000 ug

VR	V	nA	I.mean	Std.Dev.	I.delta	Comments
1 - 1	-0.206	-65.1	-65.1	0.119	0.0	
1 - 2	-0.206	-65.0				
2 - 1	-0.214	-138.0	-138.0	0.085	-73.0	
2 - 2	-0.214	-138.1				
3 - 1	-0.214	-206.1	-204.7	1.966	-66.7	
3 - 2	-0.222	-203.3				

Substance	Calibr.	Y.reg/offset	Slope	Mean deviat.	Corr.Coeff.
CN	std.add.	-6.526e-008	-2.850e-004	1.802e-009	0.99960

Final results	+/-	Res. dev.	%	Comments
---------------	-----	-----------	---	----------

CN:
Cyanide = 457.939 ug/L 9.446 2.063

Method print for the determination of free cyanide

Method parameters

Method : AB110_Det of CN.mth
Title : Determination of free Cyanide
Remark1 : 10 mL electrolyte + 10 mL sample
Remark2 :

Calibration : Standard addition
Technique : Batch

Cell volume (mL) : 20.000
Sample amount (mL) : 10.000
Sample ID : waste water

Voltammetric parameters

Mode : DP - Differential Pulse

Highest current range : 10 mA
Lowest current range : 100 nA

Electrode : DME
Stirrer speed (rpm) : 2000

No. of additions : 2
No. of replications : 2

Measure blank : No
Addition purge time (s) : 10

Initial purge time (s) : 300

Sweep
Equilibration time (s) : 5
Start potential (V) : 0.000
End potential (V) : -0.500
Voltage step (V) : 0.008
Voltage step time (s) : 0.800
Sweep rate (V/s) : 0.010
Pulse amplitude (V) : 0.050
Pulse time (s) : 0.040

Cell off after measurement : Yes

Peak evaluation

Peak evaluation : Height
Minimum peak width (V.steps) : 5
Minimum peak height (A) : 1.000e-010
Smooth factor : 4
Reverse peaks : No
Reverse sweep : No

Substances

CN : -0.240 V +/- 0.050 V
Standard solution : 1 0.100 g/L
Addition volume (mL) : 0.050

Cyanide : Final result (CN) = (1e+006 / 1) * Mass.conc + 0 - 0