

Automated sodium determination in various foods with 859 Titrotherm

Keywords

total sodium / bouillon / gravy / tomato ketchup / corn chips / pretzel sticks / crackers / 859 Titrotherm / automated / 815 Robotic USB Sample Processor

Summary

In an acidic solution (containing $\text{NH}_4\text{F} \cdot \text{HF}$, $\text{Al}(\text{NO}_3)_3$ / KNO_3), sodium forms NaK_2AlF_6 which precipitates in an exothermic solution.

Several foods were analysed, namely bouillon, gravy, tomato ketchup, corn chips, pretzel sticks as well as crackers. The reproducibility of the results was good.

After weighing in and adding solutions, the samples were crushed again with a polytron to ensure homogeneity in the measuring solution.

Relative standard deviations between 0.08% and 3.75% were obtained.

Introduction

In a titration, the titrant reacts with the analyte in the sample either exothermically (gives out heat) or endothermically (takes in heat). The Thermoprobe measures the temperature of the titrating solution. When all of the analyte in the sample has reacted with the titrant, the temperature of the solution will change and the endpoint of the titration is revealed by an inflection in the temperature curve.

The amount of analyte determined is not related to the change in temperature of the solution. Therefore, it is not necessary to use insulated titration vessels.

Theory

IN GENERAL

Thermometric titrations are conducted under conditions of constant titrant addition rate. In this respect they differ from potentiometric titrations, where the titrant addition rate may be varied during the titration according to the electrode response. In thermometric titrations, a constant addition rate of titrant equates to a constant amount of heat being given out or consumed, and hence a more or less constant temperature change up to the endpoint.

TITER DETERMINATION

If a liquid standard is used, pipette it into the titration vessel. Set up a regression plot with the sample size in mL (see appendix) on the x-axis and mL of titrant consumed on the y-axis. The plot will be a linear curve of the form $y = a \cdot x + b$, where the molarity of the titrant is calculated from the slope with the following formula: $a \cdot c(\text{Titrant})$.

BLANK DETERMINATION

The method blank is determined by titrating a number of analyte solutions of different concentrations and plotting the analyte concentration against the titrant consumption. The method blank is determined as the y-intercept from a linear regression of the titration data. Changes in method parameters will require a new determination of method blank.

This parameter is stored along with the other method parameters. For all determinations it is subtracted from the volume of titrant.

Solutions

- *Titrant:* $c(\text{Al}(\text{NO}_3)_3) = 0.5 \text{ mol/L}$ and $c(\text{KNO}_3) = 1.1 \text{ mol/L}$

187.57 g $\text{Al}(\text{NO}_3)_3$ and 111.22 g KNO_3 were weighed into a 1000 mL volumetric flask, dissolved in approx. 800 mL deion. water and filled up to the mark.

- *Complexing agent:* $\beta(\text{NH}_4\text{F} \cdot \text{HF}) = 300 \text{ g/L}$

300 g $\text{NH}_4\text{F} \cdot \text{HF}$ were weighed into a 1000 mL volumetric flask and filled up to the mark with deion. water.

- *Standard:* $c(\text{Na}_2\text{SO}_4) = 0.4 \text{ mol/L}$

28.408 g Na_2SO_4 (dried for 2 h at 105 °C) were weighed into a 500 mL volumetric flask and filled up to the mark with deion. water.

- *Deion. water*

Materials

1x	859 Titrotherm (containing 1 x 800 Dosino and dosing unit 20 mL)	2.856.1010
1x	HF resistant Thermoprobe	6.9011.040
1x	815 USB Robotic Sample Processor XL	2.815.0130
1x	Dispersing aggregate for Polytron, 157 mm	6.9012.010
1x	Swing arm for Polytron and 786 swing head, left swinging	6.1461.250
1x	Robotic arm with holder for titration head, right swinging	6.1462.070
1x	Titration head, 3 x SGJ14	6.1458.040
1x	External Polytron wash- ing station	6.2841.000
1x	Sample Rack 59 x 120 mL	6.2041.850
3x	800 Dosino	2.800.0010
2x	843 Pump Station	2.843.0150
1x	Polytron PT 1300 D	2.136.0100
2x	Dosing unit 5 mL	6.3032.150
1x	Dosing unit 5 mL ETFE	6.1575.150
1x	Dosing unit 50 mL	6.3032.250
1x	Sample beaker 100 x 120 mL	6.1459.300

Samples

- Tomato ketchup
- Bouillon as solution
- Bouillon as powder
- Gravy as powder
- Corn chips
- Pretzel sticks
- Crackers

Procedures*TITER DETERMINATION*

1 to 2.5 mL Na₂SO₄ solution, 5 mL complexing agent as well as 40 mL deion. water were dosed into the titration beaker and titrated with Al(NO₃)₃ / KNO₃ – solution up to the first end point. The consumption of the titrant was regressed against the sample size. For formula see Calculation section.

*BLANK DETERMINATION**Table 1: Sample Sizes for Blank Determination*

Sample	Sample size [g] for 1 st blank determination	Sample size [g] for 2 nd blank determination	Sample size [g] for 3 rd blank determination
Ketchup	1.35	2.21	3.65
Bouillon as solution	0.11 = 5 mL	0.18 = 8 mL	0.21 = 10 mL
Bouillon	0.17	0.21	0.24
Gravy	0.32	0.43	0.50
Cracker	1.65	2.30	3.17
Corn chips	3.06	6.32	9.58
Pretzel sticks	0.91	1.57	2.02

The Sample sizes for the blank determination were chosen depending on the sodium content. A few examples were shown in table 1. The consumption of titrant was regressed against the sample size. The intercept was taken as blank value and was saved automatically as Common variable. This blank value has to be subtracted from each further analyzed sample.

*SAMPLE DETERMINATION**Table 2: Sample Sizes for Sample Determination*

Sample	Sample size
Ketchup	3.5 g
Bouillon as solution	10 mL = 0.2 g
Bouillon	0.3g
Gravy	0.7 g
Cracker	2.7 g
Corn chips	5 g
Pretzel sticks	2 g

Before beginning the determination crush the samples with the mortar. Afterward weigh for each sample type approximately the sample sizes mentioned in Table 2, add 5 mL complexing agent, 40 mL water and titrate with Al(NO₃)₃ / KNO₃ solution up to the first end point. For the calculation of the sodium content the blank value has to be subtracted. The whole

measurement including homogenization with the Polyton takes about 8.5 minutes. The content was calculated in % as well as in mg/100g.

Titration Parameters

	Titer determination	Blank and Sample determination
Stirring rate	10	Depending on sample ¹
Start volume [mL]	0.5	0.5
Pause [s]	30	30
Switch off autom.	no	no
Dosing rate [mL/min]	4	4
Filter factor	45	Depending on sample ²
Damping until [mL]	0.5	0.5
Stop volume [mL]	5.0	5.0
Stop slope	off	off
Add. volume after stop [mL]	2	2
End points [Reaction type]	ex*	ex*
EP criterion [ERC]	-50	Depending on sample ³

* nicht erklärt. Bei end points.

1	Sample	Stir rate
	- Ketchup	10
	- Bouillon	11
	- All others	12

2	Sample	Filter factor
	- Ketchup	45
	- Bouillon as solution	45
	- Gravy	50
	- Pretzel sticks	50
	- Cracker	60
	- Corn chips	60
	- Bouillon	65

3	Sample	EP criterion
	- Ketchup	- 45
	- Bouillon as solution	- 45
	- Gravy	-40
	- Cracker	-40
	- Corn chips	-40
	- Pretzel sticks	-40
	- Bouillon	-20

Calculations

TITER DETERMINATION

Assignment	RS name	formula
RS01	EP1	'TET.EP{1}.VOL'
RS02	Slope	'RS.EP1.SLO'
RS03	Intercept	'RS.EP1.ITS'
RS04	Correlation (R ²)	'RS.EP.COR' · 'RS.EP.COR'
RS05	Molarity[mol/L]	2 · 'Add.CONC' / 'RS.Slope'
RS06	Titer	(Add.CONC · 2) / ('RS.Slope' · 'TET.CONC')

BLANK DETERMINATION

Assignment	RS name	Formula
RS01	EP1	'TET.EP{1}.VOL'
RS02	Slope	'RS.EP1.SLO'
RS03	Intercept	'RS.EP1.ITS'
RS04	Correlation (R ²)	'RS.EP.COR' · 'RS.EP.COR'

SAMPLE DETERMINATION

Assignment	RS name	Formula
RS01	EP1	'TET.EP{1}.VOL'
RS02	Content [%]	'TET.TITER' · 'TET.CONC' · ('TET.EP{1}.VOL' - 'CV.Blank') · M(Na) · 100 / (1000 · 'MV.Sample Size')
RS03	Content [mg/100g]	'TET.TITER' · 'TET.CONC' · ('TET.EP{1}.VOL' - 'CV.Blank') · M(Na) · 100 / 'MV.Sample Size'

Legend formula

TET.EP1.VOL	= thermometric end point [mL]
.SLO	= slope of the linear regression
.ITS	= intercept of the linear regression
.COR	= correlation factor
Add.CONC	= Concentration of Na ₂ SO ₄ solution = 0.4 mol/L
TET.CONC	= Concentration of Al(NO ₃) ₃ / KNO ₃ - solution = 0.5 mol/L
TET.TITER	= Titer of the Al(NO ₃) ₃ / KNO ₃ - solution
Blank	= intercept of the linear regression at the blank determination
M(Na)	= 22.990 g/mol
MV.Sample Size	= Sample size in g
2	= two Na ⁺ -Ions in Na ₂ SO ₄
100	= for % and for mg/100g
1000	= conversion factor from g to mg

Results*TITER DETERMINATION*

	Values
Slope	1.5819 mL/mL
Intercept	0.1976 mL
Molarity	0.7909 mol/L
Titer	0.9887
Corr. coefficient	1.0000

BLANK DETERMINATION

Sample	Slope [g/mL]	Intercept [mL]	Corr. Coefficient
Ketchup	1.0968	0.2143	0.99997
Bouillon as Solution	13.8983	0.2181	0.99986
Bouillon	16.3180	-0.2337	0.99975
Gravy	5.6772	0.2847	0.99813
Cracker	1.0288	0.1948	1.00000
Corn Chips	0.4203	0.0410	0.99945
Pretzel sticks	1.5793	0.2577	0.99827

*SAMPLE DETERMINATION**Ketchup*

	Calculation in [%]	Calculation in [mg/100g]
Mean	1.257	1256.8
abs. standard deviation	0.0024	2.427
rel. standard deviation	0.19 %	0.19 %
n =	6	6

Bouillon as solution

	Calculation in [%]	Calculation in [mg/100g]
Mean	16.04	16041.4
abs. standard deviation	0.013	14.101
rel. standard deviation	0.08 %	0.09 %
n =	6	6

Bouillon weighed in directly

	Calculation in [%]	Calculation in [mg/100g]
Mean	18.02	18019.9
abs. standard deviation	0.34	336.6
rel. standard deviation	1.89 %	1.87 %
n =	6	6

Gravy

	Calculation in [%]	Calculation in [mg/100g]
Mean	6.34	6335.8
abs. standard deviation	0.06	57.45
rel. standard deviation	0.95 %	0.91 %
n =	6	6

Crackers

	Calculation in [%]	Calculation in [mg/100g]
Mean	1.169	1168.9
abs. standard deviation	0.0058	5.777
rel. standard deviation	0.5 %	0.49 %
n =	6	6

Corn chips




	Calculation in [%]	Calculation in [mg/100g]
Mean	0.507	506.6
abs. standard deviation	0.019	18.76
rel. standard deviation	3.75 %	3.7 %
n =	6	6






Pretzel sticks

	Calculation in [%]	Calculation in [mg/100g]
Mean	1.809	1809.4
abs. standard deviation	0.0185	18.726
rel. standard deviation	1.02 %	1.03 %
n =	6	6

Comparison labeled / found results

Sample	Labeled [g/100g]	Found [g/100g]	Difference [%]
Ketchup	3.1 g NaCl = 1.22 g Na ⁺	1.257 g Na ⁺	+ 2.94%
Bouillon as solution	41 g NaCl =16.13 g Na ⁺	16.04 g Na ⁺	- 0.56 %
Bouillon	41 g NaCl =16.13 g Na ⁺	18.02 Na ⁺	+ 10.49 %
Gravy	6.67 g Na ⁺	6.34 g Na ⁺	- 5.21 %
Crackers	2.5 g NaCl =0.984 g Na ⁺	1.169 g Na ⁺	+ 15.83 %
Corn chips	0.5 g Na ⁺	0.507 g Na ⁺	+ 1.38 %
Pretzel sticks	4.6 g NaCl =1.810 g Na ⁺	1.809 g Na ⁺	- 0.055 %

Curves**TITER DETERMINATION****SAMPLE DETERMINATION****Ketchup****Bouillon as solution**

Bouillon weighed in directly**Corn chips****Gravy****Pretzel sticks****Crackers****Comments**

- For the cleaning of the Polytron and the electrodes, a dip rinsing was done in a beaker filled up with hexane. The further cleaning was carried out with deion. water.
- The speed of the Polytron was set to 20000 rpm.
- After crushing with the Polytron a milky suspension was obtained.
- For the dosing of the $\text{NH}_4\text{F} \cdot \text{HF}$ – solution an ETFE dosing unit was used.

Date

16.2.2010 – revised 12.4.2012


Name

Metrohm AG, CH-9100 Herisau
Iris Kalkman, CC Titration

Appendix

METHODS

Titer Determination



Blank / Sample Determination

