

## IC Application Note No. C-22

Title:	Sodium, potassium, calcium and magnesium in drip feeding formula using Dialysis as
	Sample Preparation

Summary: Determination of sodium, potassium, calcium and magnesium in drip feeding formula using cation chromatography with direct conductivity detection and dialysis as sample preparation.

Sample: Drip feeding formula


1 g/100 mL of 1 mmol/L HNO<sub>3</sub>, injection using 754 Dialysis
Unit (acceptor solution: 1 mmol/L HNO<sub>3</sub>)

Column: 6.1010.000 Metrosep Cation 1-2

Eluent: 4 mmol/L tartaric acid, 1 mmol/L dipicolinic acid

Flow: 1.0 mL/min

*Injection Volume:* 20 µL after dialysis



Results:	Na <sup>+</sup>	K <sup>+</sup>	Ca <sup>2+</sup>	Mg <sup>2+</sup>
	μg/L	μg/L	µg/L	μg/L
	682	1288	762	192

Peak No. 2: Ammonium, not quantified