


Cations and amines in the water-steam circuit



Water in steel-based cooling systems requires a pH value slightly above 7 to prevent corrosion. Often ammonium or organic amines are applied for pH adjustment. This application shows the separation of typical amines besides inorganic cations. Sample preconcentration applies combined Inline Preconcentration and Matrix Elimination (MiPCT-ME).

Results

	Concentration [mg/L]		Concentration [mg/L]
Sodium	1.0	Nickel	1.0
Ammonium	1.0	3-methoxypropylamine	1.0
Ethanolamine	1.0	Zinc	1.0
Potassium	1.0	Magnesium	1.0
Dimethylamine	1.0	Calcium	1.0
Morpholine	1.0		

Sample

Spiked water-steam circuit water

Sample preparation

None

Columns

Metrosep C 4 - 250/4.0	6.1050.430
Metrosep C 4 Guard/4.0	6.1050.500
Metrosep C PCC 1 HC/4.0	6.1010.310

Solutions

Eluent (inline eluent preparation)	2.5 mmol/L nitric acid 0.5 mmol/l <u>oxalic acid</u>
Liquid handling	Ultrapure water

Parameters

Flow rate	0.9 mL/min
Injection volume	100 µL
P _{max}	20 MPa
Recording time	30 min
Column temperature	32 °C

Analysis

Direct conductivity detection

Instrumentation


850 Professional IC – Cation	2.850.1010
IC Conductivity Detector	2.850.9010
858 Professional Sample Processor	2.858.0010
2 × 800 Dosino	2.800.0010
849 Level Control for Inline Eluent Preparation	2.849.1030

Calibration MiPCT-ME

Calibration range	Factor of 100
Standard solution:	
All cations	10 mg/L
1. Level	5 µL = 0.1 mg/L
2. Level	10 µL = 0.2 mg/L
3. Level	25 µL = 0.5 mg/L
4. Level	50 µL = 1.0 mg/L
5. Level	100 µL = 2.0 mg/L
6. Level	250 µL = 5.0 mg/L
7. Level	500 µL = 10.0 mg/L



Flow chart



www.metrohm.com

 **Metrohm**