




Thermo. Titr. Application Note No. 71

Title: Determination of Ammonium Ion by Titration with Hypochlorite

Scope: Determination of ammonium ion in ammonium salts and mixtures containing ammonium ion

Principle: Hypochlorite ions react with bromide ions to form hypobromite ions, which in turn rapidly oxidize ammonium ions to nitrogen. Hypobromite reacts more rapidly with ammonium than hypochlorite, and is formed *in situ* (Vogel, 1961). The titration is carried out with in a solution containing bromide and bicarbonate:



Reagents:

~0.25mol/L NaOCl solution. Prepare from commercial concentrated sodium hypochlorite solution. Sodium hypochlorite solution from a swimming pool supply shop is satisfactory. Dilute 185-200mL of concentrated solution to 1000mL. Store in a brown glass bottle.

Bromide/bicarbonate solution. Dissolve 200g KBr and 120g KHCO₃ in DI water and make to 1000mL with DI water. Store in a brown glass bottle.

1mol/L Na₂S₂O₃ solution. Dissolve 248.2 A.R. Na₂S₂O₃·5H₂O in DI water and make to 1000mL. Store in a brown glass bottle.

0.04mol/L KIO₃ solution. Dry at least 10g A.R. KIO₃ at 120°C for 2 hours. Cool in a dessicator, and weigh accurately approximately 8.56g. Make to 1000mL with DI water in a volumetric flask.

4mol/L KI solution. Dissolve 166g A.R. KI in DI water, and make to 250mL with DI water. Store in a brown glass bottle.

Glacial acetic acid.

Sulfuric acid, 25% w/v.

Method:*Basic Experimental Parameters:*

Titrant delivery rate (mL/min.)	4
Titration blank (mL)	0.050
Bromide/bicarbonate solution addition (mL):	10
No. of exothermic endpoints	1
Delay start of titration (secs.)	5
Data smoothing factor (DSF)	60
Stirring speed (802 stirrer)	8

Basic Method: The sample added to the titration vessel should ideally contain approximately 6 - 8mmole NH_4^+ .

In the examples using reagent grade ammonium salts illustrated here, solutions were made to approximately 0.2mol/L with respect to NH_4^+ , and added using a 800 Dosino. This was done to reduce errors in dispensing aliquots.

For the examples of commercial fertilizers illustrated here, weigh accurately approximately 1g of finely ground MAP, and transfer to a 250mL volumetric flask. Make to volume and mix well. Add a small magnetic spin bar, and stir contents for 10 minutes. Allow to stand for a few minutes to allow coarser solids to settle. The solution can be filtered if desired, but this is not necessary. Take a 25mL aliquot for analysis. This will contain the equivalent of approximately 0.1g of original material.

Standardization of NaOCl titrant.





It is necessary to trace the calibration of the NaOCl solution to a recognized primary standard. NaOCl solutions are conveniently calibrated against standard $\text{Na}_2\text{S}_2\text{O}_3$ solution, but $\text{Na}_2\text{S}_2\text{O}_3$ is not a recognized primary standard. However, $\text{Na}_2\text{S}_2\text{O}_3$ solutions may be standardized against solutions of KIO_3 , which is a primary standard.

Standardization of $\text{Na}_2\text{S}_2\text{O}_3$ solution. Prepare a series of solutions in titration vessels as follows:

mL KIO ₃ solution	mL DI water	mL 25% H ₂ SO ₄
5	25	5
10	20	5
15	15	5
20	10	5
25	5	5

Prepare a titration program with a Na₂S₂O₃ titrant dose rate of 4mL/min. and a digital filtration value of 60. Immediately prior to the titration of each prepared solution, and while the solution is stirring, add 10mL 4mol/L KI solution. This can be added conveniently from a Dosino

Prepare a regression analysis, plotting mmol KIO₃ (x-axis) against mL Na₂S₂O₃ (y-axis), and calculate the gradient to the line of best fit. The molarity is calculated as 6/gradient. This is derived from:









So that 1 mol IO₃⁻ ≡ 6 mol S₂O₃²⁻

Standardization of NaOCl against Na₂S₂O₃. Because NaOCl solutions do not possess good long term stability, it is necessary to check the molarity of the NaOCl on a regular basis, say every few days. It is thus convenient to create a special standardization program in Titrotherm, where the NaOCl solution is dispensed as a pre-dose which is titrated with Na₂S₂O₃. Prepare solutions of 25mL DI water to which has been added 5mL 4mol/L KI and 5mL glacial acetic acid. In successive titrations, program the addition of 2, 4, 6, 8 and 10mL of NaOCl prior to the titration.

Prepare a regression analysis, plotting mmol Na₂S₂O₃ (from the volume of Na₂S₂O₃ titrated) on the

x-axis, and mL of NaOCl dispensed on the y-axis. The molarity is calculated as 0.5/gradient. This is derived from:



So that 2 mol $\text{Na}_2\text{S}_2\text{O}_3 \equiv 1$ mol NaOCl

Sample Preparation.

MAP Samples. Weigh accurately approximately 1g of finely ground MAP, and transfer to a 250mL volumetric flask. Make to volume and mix well. Add a small magnetic spin bar, and stir contents for 10 minutes. Allow to stand for a few minutes to allow coarser solids to settle. The solution can be filtered if desired, but this is not necessary.

A 25mL aliquot is taken for titration. This will contain the equivalent of approximately 0.1g of original material.

“Ammoniated” TSP sample. As for MAP samples, but weigh out accurately approximately 2g of material.

Titration. Set up a titration program with a NaOCl titration rate of 4mL/min, a digital filter value of 60, and a 10mL pre-dose of bromide/bicarbonate reagent

Examples:

Salts analyzed were reagent grade, but had been opened for some years (except for the $(\text{NH}_4)_2\text{SO}_4$, which was freshly opened)

Salt	NH_4^+ content, %	
	Theoretical	Analyzed
$(\text{NH}_4)_2\text{SO}_4$	27.32	27.26±0.007 (n=5)
$\text{NH}_4\text{H}_2\text{PO}_4$	15.68	15.71±0.005 (n=5)
$(\text{NH}_4)_2\text{HPO}_4$	27.31	27.22±0.010 (n=5)
$\text{AlNH}_4(\text{SO}_4)_4 \cdot 12\text{H}_2\text{O}$	3.98	3.99±0.001 (n=5)

Fertilizers were supplied by a customer

Sample	Total Nitrogen	
	As %N	As % NH_4^+
MAP #1	11.09±0.03 (n=9)	14.28±0.04
MAP #2	10.94±0.01 (n=8)	14.09±0.01
“ammoniated” TSP	3.08±0.01 (n=8)	3.97±0.01


Calculation:

$$\% \text{NH}_4^+ = \frac{((\text{Titre, mL} - \text{blank, mL}) \times M \text{NaOCl} \times 18.034 \times 100 \times 2)}{(\text{sample mass, g} \times 1000 \times 3)}$$

Titration Plot:

Legend:

Red = solution temperature curve
Black = second derivative curve (for endpoints)



Titrant calibration:

*Gradient = 1.94674,
molarity = 0.5/gradient =
0.2568 mol/L*

