
Merkblatt 6.2802.200 Wartungsset zu Arbeitselektroden 6.1257.2XX Leaflet 6.2802.200 Maintenance set for 6.1257.2XX working electrodes

Teile / Parts

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EN

Leaflet 6.2802.200 Maintenance set for 6.1257.2XX working electrodes

General

The maintenance set is exclusively intended for use as a manually operated tool for regenerating heavily eroded working electrodes (6.1257.2XX). It allows users to grind the affected surface to a planar parallel state again.

To ensure that the measuring cell remains tight even with a ground working electrode, the working electrode requires a subsequent finish grinding step.

If the working electrode is used in the DC mode, then it can be additionally polished with the polishing set (6.2802.000) to ensure optimum performance.

Contents of the maintenance set

- The maintenance set includes:
- Base with guide rod (6)

- Threaded sleeve (4) to secure the grinding disk in place
- Electrode holder (1)
- Bag containing grinding disks (5), 3 disks (coarse-grit) for main grinding, 5 disks (fine-grit) for finish grinding (6.2802.210)

Application

Note

We recommend that the working electrode should always be ground with a moist grinding disk. This prevents the formation of grinding dust floating around.

Note

When the electrode is ground to the maximum (approx. 1.5 mm), the electrode holder rests on the edge of the threaded sleeve, thus preventing further grinding of the electrode surface. If this is the case, then the electrode can no longer be regenerated and must be replaced.

Regeneration grinding

1. Place a coarse-grit grinding disk onto the base's guide rod and screw on the threaded sleeve. This secures the grinding disk in place so that it can no longer slip.
2. Evenly wet the surface of the grinding disk with 2 - 3 drops of water.
3. Insert the shaft of the working electrode (3) into the electrode holder and secure it in place with the knurled screw (2).
4. Put the electrode holder with the working electrode in place onto the guide rod.
5. Rotate the electrode holder against the base without applying too much pressure. The working electrode is thus ground.
6. Take off the electrode holder from time to time to check the progress.

7. When the surface has a smooth and clean appearance, the working electrode requires finish grinding.

Finish grinding

1. Replace the coarse-grit grinding disk with a fine-grit grinding disk.
2. Evenly wet the surface of the grinding disk with 2 - 3 drops of water.
3. Carry out finish grinding for the working electrode while applying almost no pressure at all.

Notes

- Use a separate grinding disk for each electrode material.
- Grinding disks may be only slightly wet. Grinding disks that are too wet tend to form waves and folds.
- To prevent annoying squeaking and premature wear, you may lubricate the guide rod with a little petroleum jelly.
- After several years of use, the tool might lose its accuracy (the backlash between electrode holder and guide rod might increase too much). This increases the risk of the electrodes being ground with a slope or too round, possibly resulting in leakage. If this is the case, the tool must be replaced.
- In isolated cases, an electrode might become leaky due to grinding.
 - If liquid leaks between the electrode and the cell, try to grind the electrode again to a flat surface.
 - If liquid leaks at the electrode plug, the working electrode can no longer be regenerated and must be replaced.