

# Monograph



## Complexometric (Chelometric) Titrations

Peter Bruttel, revised by Iris Kalkman and Lucia Meier



# Table of contents

<b>Preface</b>	<b>06</b>
<b>Theoretical fundamentals</b>	<b>07–23</b>
<b>Indication of the equivalence point of the titration</b>	<b>09–23</b>
Visual and photometric	09
Potentiometric	10–11
Commonly used complexing agents/titrants	12
Buffer solutions	13
Masking	14
Color indicators	15–21
Preparation of the standard solutions	22
Titer determinations	22–23
<b>Practical examples</b>	<b>24–39</b>
<b>Water hardness, calcium, and magnesium</b>	<b>24–29</b>
Potentiometric titration with the Ca-ISE	25–26
Potentiometric titration with the Cu-ISE	26–27
Photometric titration with the Optrode at 610 nm	27–29
<b>Aluminum and magnesium, alloys and antacids</b>	<b>29–31</b>
Potentiometric titration with the Cu-ISE	31–33
<b>Magnesium and zinc</b>	<b>31–33</b>
Potentiometric titration with the Cu-ISE	31–33
<b>Aluminum, calcium, iron, and magnesium in cement and clinker</b>	<b>33–37</b>
Photometric titrations with the Optrode at 610 nm	33–37
<b>Sulfate</b>	<b>38–39</b>
Potentiometric with the Ca-ISE	38–39
<b>Single determinations</b>	<b>39–92</b>
<b>Al – Aluminum</b>	<b>39–42</b>
General	39–40
Potentiometric with the Cu-ISE	40–41
Photometric with the Optrode at 610 nm	41–42
<b>Ba – Barium and Sr – Strontium</b>	<b>43–45</b>
General	43
Potentiometric with the Cu-ISE	44
Photometric with the Optrode at 574 nm	45
<b>Bi – Bismuth</b>	<b>46–48</b>
General	46
Potentiometric with the Cu-ISE	46–47
Photometric with the Optrode at 520 nm	47–48
<b>Ca – Calcium</b>	<b>48–51</b>
General	48–49
Potentiometric with the Ca-ISE	49–50
Potentiometric with the Cu-ISE	50
Photometric with the Optrode at 610 nm	51

<b>Cd – Cadmium</b>	<b>52–54</b>
General	52
Potentiometric with the Cu-ISE	52–53
Photometric with the Optrode at 610 nm	53–54
<b>Co – Cobalt</b>	<b>54–56</b>
General	54
Potentiometric with the Cu-ISE	55
Photometric with the Optrode at 574 nm	56
<b>Cu – Copper</b>	<b>57–59</b>
General	57
Potentiometric with the Cu-ISE	58
Photometric with the Optrode at 520 nm	59
<b>Fe – Iron</b>	<b>60–62</b>
General	60
Potentiometric with the Cu-ISE	61
Photometric with the Optrode at 610 nm	62
<b>Ga – Gallium and In – Indium</b>	<b>63–66</b>
General	63
Potentiometric with the Cu-ISE	63–64
Photometric with the Optrode at 610 nm	65
<b>Hg – Mercury</b>	<b>66–67</b>
General	66
Photometric with the Optrode at 502 nm	67
<b>Mg – Magnesium</b>	<b>68–70</b>
General	68
Potentiometric with the Cu-ISE	68–69
Photometric with the Optrode at 610 nm	69–70
<b>Mn – Manganese</b>	<b>70–73</b>
General	70–71
Potentiometric with the Cu-ISE	71–72
Photometric with the Optrode at 610 nm	72–73
<b>Ni – Nickel</b>	<b>73–76</b>
General	73–74
Potentiometric with the Cu-ISE	74–75
Photometric with the Optrode at 574 nm	75–76
<b>Pb – Lead</b>	<b>76–79</b>
General	76–77
Potentiometric with the Cu-ISE	77–78
Photometric with the Optrode at 574 nm	78–79
<b>Pd – Palladium</b>	<b>79–80</b>
General	79
Photometric with the Optrode at 610 nm	80

<b>Sn – Tin</b>	<b>81–82</b>
General	81
Photometric with the Optrode at 574 nm	81–82
<b>Th – Thorium</b>	<b>82–84</b>
General	82–83
Photometric with the Optrode at 574 nm	83–84
<b>Tl – Thallium</b>	<b>84–86</b>
General	84
Potentiometric with the Cu-ISE	85–86
Photometric with the Optrode at 610 nm	86
<b>Zn – Zinc</b>	<b>87–89</b>
General	87
Potentiometric with the Cu-ISE	87–88
Photometric with the Optrode at 610 nm	88–89
<b>Zr – Zirconium and Hf – Hafnium</b>	<b>89–92</b>
General	89–90
Potentiometric with the Cu-ISE	90–91
Photometric with the Optrode at 520 nm	91–92
<b>Appendix</b>	<b>93–95</b>
<b>Index</b>	<b>93–94</b>
<b>Photometric indication of the equivalence point of the titration</b>	<b>95</b>
<b>Table of figures</b>	<b>96–98</b>

# Preface

06

In 1945, Gerold Schwarzenbach discovered that aminocarboxylic acids (e.g., NTA, EDTA) form stable complexes with metal ions. On the basis of this discovery, he developed the complexometric titration technique, which became increasingly popular from 1950 onwards, particularly for determining water hardness. Soon, in addition to calcium and magnesium, other metal ions could also

be titrated. The use of new metal indicators and masking agents also enabled determinations of metal ion mixtures. Nowadays, together with acid/base, redox, and precipitation titrations, complexometric titrations are the most frequently used volumetric methods, and are recommended in many international standards and directives.



# Theoretical fundamentals

The term *complex* refers to species composed of individual ions or molecules with the formula  $[ML_n]^{±z}$ .

The *coordination number* specifies how many monodentate ligands (or complexing agents) are bound. Cations are usually present in aqueous solutions as *aqua complexes* with the formula  $[M(OH_2)_4]^{2+}$  or  $[M(OH_2)_6]^{2+}$ .

The term *denticity* refers to the number of coordination sites of the respective ligand. Complexes with multidentate ligands are known as *chelates*.

The stability of a complex is determined by thermodynamic and kinetic factors. A measure of the thermodynamic stability is the *complex formation constant*  $K_f$  or the *dissociation constant*  $K = 1/K_f$ .

The stability of metal complexes can be described, for example, by the following simplified equation, which uses EDTA as ligand:

$$K_f = \frac{[MeEDTA^{2-}] \times [H_3O^+]^2}{[Me^{2+}] \times [H_2EDTA^{2-}]}$$



The complex formation constant is usually given as  $\log K_f$  in data tables. The stability of the complex increases with increasing magnitude of  $\log K_f$ .

The effective complex formation constants are generally lowered in the presence of protons (acids) because these compete with the metal ions in the reaction with the complexing agent, for example, by protonating the carboxyl group(s).

As a general rule:

- Complex formation constant  $< 10$ ; alkaline titration (e.g.,  $Ca^{2+}$ ,  $Mg^{2+}$ )
- Complex formation constant  $> 15$ ; slightly acidic titration (e.g.,  $Al^{3+}$ ,  $Pb^{2+}$ ) ( $Fe^{3+}$  and  $Bi^{3+}$  can even be titrated at  $pH = 2$ ).

Titration must be carried out in buffered solutions because the complexation reaction always liberates protons. This is illustrated briefly by the following examples:



Apart from a few exceptions, the chelating agents used for titration form 1:1 complexes with metal ions, regardless of whether the metal ions have a charge  $z = 2+$  or  $z = 3+$ .

The following table contains a few selected complex formation constants ( $\log K_f$ ) for common titrants and metal ions:

Metal ion	EDTA <sup>1</sup>	EGTA <sup>2</sup>	DCTA <sup>2</sup>	DTPA <sup>2</sup>	NTA <sup>2</sup>
Al(III)	16.4	13.9	18.6	18.4	9.5
Ba(II)	7.9	8.4	8.6	8.6	4.8
Bi(III)	27.8 <sup>a</sup>	23.8	31.2	29.7	–
Ca(II)	10.7	11.0	12.5	10.7	6.4
Cd(II)	16.5	16.7	19.2	19.3	9.5
Co(II)	16.5	12.5	18.9	18.4	10.4
Co(III)	41.4	–	–	–	–
Cr(II)	13.6 <sup>a</sup>	–	–	–	–
Cr(III)	23.4 <sup>a</sup>	2.5	–	–	> 10.0
Cu(II)	18.8	17.8	21.3	21.5	13.0
Fe(II)	14.3	11.9	16.3	16.6	8.8
Fe(III)	25.1	20.5	28.1	28.6	15.9
Ga(III)	21.7	–	22.9	23.0	13.6
Hf(IV)	29.5	–	–	35.4	20.3
Hg(II)	21.5	23.1	24.3	27.0	14.6
In(III)	24.9	–	28.8	29.0	16.9
Mg(II)	8.8	5.2	10.3	9.3	5.5
Mn(II)	13.9	12.3	16.8	15.6	7.4
Ni(II)	18.4	13.6	19.4	20.3	11.5
Pb(II)	18.0	14.7	19.7	18.8	11.4
Pd(II)	25.6 <sup>a</sup>	–	–	–	–
Sn(II)	18.3 <sup>b</sup>	23.9	–	–	–
Sr(II)	8.7	8.5	10.5	9.7	5.0
Th(IV)	23.2	–	29.3	28.8	12.4
Tl(I)	6.4	–	5.3	6.0	4.8
Tl(III)	35.3	–	38.3	48.0	18.0
Zn(II)	16.5	14.5	18.7	18.8	10.7
Zr(IV)	29.3	–	20.7	36.9	20.85

DCTA	trans-Diaminocyclohexanetetraacetic acid
DTPA	Diethylenetriaminepentaacetic acid
EDTA	Ethylenediaminetetraacetic acid
EGTA	Ethylene glycol bis-(2-aminoethyl)tetraacetic acid
NTA	Nitrilotriacetic acid

The following conclusions can be drawn from the above table:

- Metal ions with a triple positive charge form stronger complexes than metal ions with a double positive charge (cf.  $\text{Fe}^{2+}$  –  $\text{Fe}^{3+}$ ).
- Ligands with a higher denticity form stronger complexes (cf. NTA – DTPA).
- The weakest complex is formed between Ba and NTA, and the strongest between  $\text{Tl}^{3+}$  and DTPA.

<sup>1</sup> If not otherwise indicated, the values are based on measurements at 25 °C and an ionic strength of 0.1 mol/L.

<sup>a</sup> Values measured at 20 °C and an ionic strength of 0.1 mol/L.

<sup>b</sup> Values measured at 20 °C and an ionic strength of 1 mol/L.

Harris, Daniel (2007): Quantitative Chemical Analysis. 7<sup>th</sup> ed., New York.

<sup>2</sup> Values were taken from Dojindo Molecular Technologies (publisher) <http://www.dojindo.com> (Issued 28.9.2015).



# Indication of the equivalence point of the titration

09

## Visual or photometric

Visual indication with color indicators is the oldest method of determining the equivalence point of a titration, and it is still frequently used. It is inexpensive and requires little equipment. The particular drawback of this method is that it cannot be automated and is difficult to validate. The color sensitivity of individuals differs, and also depends on the lighting conditions. Furthermore, the endpoint is even more difficult to observe in colored and/or turbid solutions.

The endpoint detection can be improved by the use of photometric indication. In this case, the subjective human eye is replaced by an objective sensor. This kind of indication can be automated and validated – provided that the correct color indicator is chosen and the solutions are not too turbid before or during the titration.

Color indicators do not have a precise *color change point* – rather, they have a *color change range*. This has a considerable influence on the accuracy of the results of titrations where visual indication of the equivalence point is employed. As a general rule: the human eye can perceive a color change when the concentration ratio changes from roughly 1/10 to 10/1. It is thus expedient to select color indicators that have a sharp color change with very different colors.

The color indicators react with the analyte metal ions to form a color complex whose color suddenly changes when all the metal ions have been (re)complexed or titrated by the titrant.

Photometric titration is carried out with fiber-optic sensors. These comprise a light source, two glass fiber (optical) conductors, and an amplifier. Monochromatic light exits the fiber-optic conductor and enters the solution where its color is partially absorbed. Depending on the design of the photometric sensor, the optical signal is sent to a photodiode directly or via a concave mirror and a second fiber-optic conductor. This generates a voltage signal in the amplifier that is sent to the titrator via its electrode input.

During photometric titrations, attention should be paid to the fact that air bubbles in the sample solution may have an unfavorable effect on the shape of the titration curve, and may also influence the results, leading to spikes or equivalence point artifacts.

The following procedure generally applies:

- Stirring is adjusted so that the solution is thoroughly mixed, but no air bubbles are entrained.
- The solution is degassed under vacuum or in an ultrasonic bath before the analysis. If necessary, the sample solution is diluted with dist. H<sub>2</sub>O and pre-neutralized.

### Potentiometric

The measuring setup for potentiometric measurements always consists of two electrodes – an indicator electrode (measuring electrode) and a reference electrode. These do not measure potentials (also known as Galvanic potentials), but potential differences (volts).

The indicator electrode (ion-selective electrode, ISE) provides an electrode potential that depends on the composition of the measuring solution (specifically, the concentration of the measuring ion).

The reference electrode (usually Ag/AgCl) provides an electrode signal (reference potential) that is as independent as possible from the measuring solution.

The potential is measured under almost zero-current conditions using a voltmeter (e.g., the titrator) with a high-resistance measuring input. This is important because it avoids voltage drops.

If a suitable measuring setup is used, the potential  $U$  measured between the two electrodes depends only on the measuring solution, or to be more precise, on the activity  $a_i$  of the measuring ion, since only dissociated ions are measured. In case of ISE measurements, this relationship is described by the Nernst equation:

$$U = U_0 + \frac{2.303 \times R \times T}{z \times F} \times \log a_i = U_0 + U_N \times \log a_i$$

U: Measured potential between measuring and reference electrodes

U<sub>0</sub>: Standard potential of the measuring chain (depends on its configuration)

2.303: Conversion factor from the natural to the decadic logarithm

R: Gas constant (8.31441 J/K/mol)

T: Absolute temperature in Kelvin (273.15 + x °C)

z: Charge of the measuring ion, including algebraic sign (e.g., +2 for Ca<sup>2+</sup>)

F: Faraday constant (96484.56 C/mol)

a<sub>i</sub>: Activity of the measuring ion

U<sub>N</sub>: Nernst slope (e.g., 29.58 mV at 25 °C and z = +2)

The Nernst slope  $U_N$  is the theoretical electrode slope. It corresponds to a potential change induced by a change in  $a_i$  by one order of magnitude. It depends on the temperature of the measuring solution and on the charge  $z$  of the measuring ion. The following table illustrates these dependencies:

Temperature / °C	U / mV, z = +1	U / mV, z = +2
0	54.20	27.10
10	56.18	28.09
20	58.17	29.09
25	59.16	29.58
30	60.15	30.08
40	62.14	31.07
50	64.12	32.06

### Ion-selective electrodes

As their name indicates, ion-selective electrodes respond more or less selectively to the ion for whose measurement they are designed. This suggests that a corresponding ISE is required for each measuring ion to be determined, which would make indication of the equivalence point using ion-selective electrodes extremely unattractive. Furthermore, such a concept would be much too expensive and also difficult to implement with certain measuring ions. How to manage with only one ion-selective electrode is explained in the chapter Cu-ISE.

#### a) Ca-ISE

The Ca-ISE is a polymer membrane electrode and, as its name indicates, it is particularly responsive to Ca ions. It is therefore mainly used in the titration of calcium. However, there is a trick that allows simultaneous titration of Ca and Mg (e.g., to determine the water hardness).


This involves treating the measuring solution with an auxiliary complexing agent (acetylacetonate in TRIS buffer). If the mixture is then titrated with EDTA, Ca is complexed first. When there is no Ca left, the auxiliary complexing agent comes into play. Mg reacts with CaEDTA to form MgEDTA. The corresponding quantity of Ca is liberated and can then be titrated with EDTA.

The determination of sulfate exploits the cross-sensitivity of Ca ISE with respect to Ba ions. The sulfate is precipitated with  $\text{BaCl}_2$ , and the excess Ba is then back-titrated with EDTA. In the presence of Ca, the titration has two equivalence points. The first equivalence point corresponds to that of Ca and the second to that of Ba.

More details of these two strategies are given in the practical examples (determination of water hardness and of sulfate).

### b) Cu-ISE

The Cu-ISE is a crystal membrane electrode. In theory, it responds only to Cu ions. However, once again a trick allows it to be used not only for titrating Cu ions. To do this, a small amount of CuEDTA is added to the sample solution. CuEDTA is a very stable complex, and liberates only a very small fraction of the Cu(II) ions. However, in alkaline solutions and in the presence of Ni(II) ions, for example, all Cu(II) ions are liberated from the complex, and the NiEDTA complex forms instead. If the mixture is then titrated with EDTA, the concentration of the Cu(II) ions decreases continuously until it reaches its very small initial value at the equivalence point. The following equations illustrate the three steps:



#### Commonly used complexing agents/titrants

The first complexing agent used for titrations was NTA or its trisodium salt. NTA is cheap, but it has the disadvantage that its complex formation constants are relatively small. It was therefore soon replaced by EDTA, which is still the most common agent used today, and is the specified titrant in many national and international standards.

Additional complexing agents were developed for special purposes (e.g., better selectivity, higher complex formation constants, better separation of mixtures); however, most of them are relatively expensive.

Commercially available complexing agents include:

#### NTA

- Nitritriacetic acid  $\text{C}_6\text{H}_9\text{NO}_6$ ;  
 $M_R = 191.14 \text{ g/mol}$
- Nitritriacetic acid trisodium salt monohydrate  $\text{C}_6\text{H}_6\text{NNa}_3\text{O}_6 \cdot \text{H}_2\text{O}$ ;  
 $M_R = 275.12 \text{ g/mol}$

#### EDTA

- Ethylenediaminetetraacetic acid  $\text{C}_{10}\text{H}_{16}\text{N}_2\text{O}_8$ ;  $M_R = 292.25 \text{ g/mol}$
- Ethylenediaminetetraacetic acid disodium salt dihydrate  $\text{C}_{10}\text{H}_{14}\text{N}_2\text{Na}_2\text{O}_8 \cdot 2\text{H}_2\text{O}$ ;  
 $M_R = 372.24 \text{ g/mol}$

#### EGTA

- Ethylene glycol-bis-(2-aminoethyl)tetraacetic acid  $\text{C}_{14}\text{H}_{24}\text{N}_2\text{O}_{10}$ ;  
 $M_R = 380.35 \text{ g/mol}$

#### DCTA


- trans-1,2-Diaminocyclohexanetetraacetic acid monohydrate  $\text{C}_{14}\text{H}_{22}\text{N}_2\text{O}_6 \cdot \text{H}_2\text{O}$ ;  
 $M_R = 364.36 \text{ g/mol}$

#### DTPA

- Diethylenetriaminepentaacetic acid  $\text{C}_{14}\text{H}_{23}\text{N}_3\text{O}_{10}$ ;  $M_R = 393.35 \text{ g/mol}$

### Buffer solutions

As already mentioned in the section on theoretical fundamentals, two protons are liberated per metal ion during complexation. The sample solution must be buffered to prevent the formation of metal complexes being hindered at lower pH. This is particularly important for metal ions that must be titrated in alkaline solutions.



#### a) Alkaline buffer solutions

The most commonly buffer is probably a mixture of ammonia and ammonium chloride – its pH value is approx. 10.

The solution normally contains 5–8 mol/L  $\text{NH}_3$  and 1 mol/L  $\text{NH}_4\text{Cl}$ .

In special cases, the following alternatives can be used:

- 1 mol/L each of NaOH and  $\text{H}_3\text{BO}_3$
- 0.2 mol/L tris(hydroxymethyl)amino-methane (TRIS)
- 1 mol/L triethanolamine, adjusted to pH 8 with HCl

#### b) Acidic buffer solutions

The most commonly used solution in this case is the acetic acid/ammonium acetate buffer: a mixture of 1.4 mol/L  $\text{CH}_3\text{COOH}$  and 1.5 mol/L  $\text{NH}_4\text{CH}_3\text{COO}$  has a pH value of approx. 4.7. Buffer solutions can be prepared in the pH range between 4 and 5, mainly by varying the amount of acetic acid.

Highly acidic buffer solutions (pH 2–3) are prepared by adding HCl to 0.5 mol/L glycine.

For all buffer solutions, it is important that their constituents do not form a stronger complex with the analyte metal ion than does the titrant.

## Masking

Masking agents can be added to the sample solution to prevent interfering metal ions from reacting with the titrant. These form stronger complexes or bonds with the interfering ions than does the titrant. In this case, it is important that the metal ion being determined is not masked itself.

### a) Masking by precipitation

Masking by precipitation means that the interfering metal ions are precipitated from the solution; however, in contrast to a physical separation, the precipitates are left in the solution and are not removed by filtration.

This method was used at an early stage in the photometric determination of water hardness. The sum of calcium and magnesium (total hardness) is determined in a first titration at  $\text{pH} = 10$ . The

second titration is then carried out at  $\text{pH} = 12$ . Magnesium is precipitated as  $\text{Mg}(\text{OH})_2$  and can thus no longer interfere in the titration – only the calcium is determined. However, the results for calcium are generally too low because some of it adsorbs onto the precipitate. This is a general problem in precipitation masking. A further disadvantage of this method is that the color indicator is often adsorbed, or its color change is more difficult to detect.

### b) Masking by complexation

The masking agent used for masking by complexation must be a selective complexing agent that binds the analyte metal ion as weakly as possible and the interfering ion as strongly as possible.

A small selection of commonly used masking agents

Masking agent	Element
Ammonium fluoride	Al, Ti, Be, Ca, Mg, Sr, Ba
Potassium cyanide	Zn, Cd, Hg, Cu, Ag, Ni, Co
Acetylacetone	Fe, Al, Pd, $\text{UO}_2$
Tiron (dihydroxybenzenedisulfonic acid)	Al, Fe, Ti
Triethanolamine	Fe, Al
2,3-Dimercaptopropanol (BAL)	Zn, Cd, Hg, As, Sb, Sn, Pb, Bi

### Color indicators

A large number of color indicators have already been described, and all have advantages and disadvantages. Some of them have become well established and

can be used for most complexometric titrations. A small overview of such metallochromic indicators is given below:

Indicator	Element	Color with metal ion	Color without metal ion
Pyrocatechol violet	$\text{Bi}^{3+}$ , $\text{Cd}^{2+}$ , $\text{Co}^{2+}$ , $\text{In}^{3+}$ , $\text{Mn}^{2+}$ , $\text{Zn}^{2+}$	blue	yellow
Dithizone	$\text{Zn}^{2+}$	red	green-violet
Eriochrome Black T	$\text{Ca}^{2+}$ , $\text{Cd}^{2+}$ , $\text{Hg}^{2+}$ , $\text{Mg}^{2+}$ , $\text{Pb}^{2+}$ , $\text{Zn}^{2+}$	red	blue
Calconcarboxylic acid	$\text{Ca}^{2+}$ in the presence of significant $\text{Mg}^{2+}$	red	blue
Hydroxynaphthol blue	$\text{Ca}^{2+}$ , $\text{Mg}^{2+}$	red	blue
Murexide	$\text{Co}^{2+}$ , $\text{Cu}^{2+}$ , $\text{Ni}^{2+}$	yellow	violet
PAN	$\text{Cd}^{2+}$ , $\text{Cu}^{2+}$ , $\text{Zn}^{2+}$	red	yellow
Phthalein purple	$\text{Ba}^{2+}$ , $\text{Ca}^{2+}$ , $\text{Sr}^{2+}$	purple	colorless
Tiron	$\text{Fe}^{3+}$ , $\text{SO}_4^{2-}$	blue	yellow
Xylenol orange	$\text{Al}^{3+}$ , $\text{Bi}^{3+}$ , $\text{La}^{3+}$ , $\text{Sc}^{3+}$ , $\text{Th}^{4+}$ , $\text{Zr}^{4+}$	red	yellow


Most metallochromic indicators are unstable in solution. They are thus often added in a solid form (as tablets or triturations with NaCl) to the solution being titrated. This circumstance, however, makes it poorly suited or even unsuitable for automation purposes. This can be overcome by preparing alcoholic or

stabilized aqueous solutions. The shelf life of the resulting indicator solution must be determined on a case-by-case basis.


The visible spectra of some common indicators are shown below.

### a) Pyrocatechol violet

Pyrocatecholsulfonephthalein;  $C_{19}H_{14}O_7S$   
0.1% in water



**Figure 1** Overlaid Vis spectra of complexed and uncomplexed pyrocatechol violet




**Figure 2** Vis absorption difference spectrum of complexed and uncomplexed pyrocatechol violet


### b) Dithizone

1,5-Diphenylthiocarbazone;  $C_{13}H_{12}N_4S$

Approx. 30 mg in 100 mL ethanol (shelf life of 1 month in a refrigerator)



**Figure 3** Overlaid Vis spectra of complexed and uncomplexed dithizone




**Figure 4** Vis absorption difference spectrum of complexed and uncomplexed dithizone




### c) Eriochrome Black T

1-(1-Hydroxy-2-naphthylazo)-6-nitro-2-naphthol-4-sulfonic acid sodium salt;  
 $C_{20}H_{12}N_3NaO_7S$

1% triturated in NaCl or 100 mg in 100 mL ethanol or 100 mg in 100 mL water  
 and stabilized with 100 mg ascorbic acid



**Figure 5** Overlaid Vis spectra of complexed and uncomplexed Eriochrome Black T




**Figure 6** Vis absorption difference spectrum of complexed and uncomplexed Eriochrome Black T


### d) Calconcarboxylic acid (HHSNN)

2-Hydroxy-1-(2-hydroxy-4-sulfo-1-naphthylazo)-3-naphthoic acid;  $C_{21}H_{14}N_2O_7S$

1% triturated in NaCl or 50 mg in 100 mL  $c(NaOH) = 0.1$  mol/L (shelf life of one day)




**Figure 7** Overlaid Vis spectra of complexed and uncomplexed calconcarboxylic acid




**Figure 8** Vis absorption difference spectrum of complexed and uncomplexed calconcarboxylic acid

### e) Hydroxynaphthol blue

2,2'-Dihydroxy-1,1'-azonaphthalene-3',4,6'-trisulfonic acid trisodium salt;  
 $C_{20}H_{11}N_2Na_3O_{11}S_3$   
 10% triturated in NaCl




**Figure 9** Overlaid Vis spectra of complexed and uncomplexed hydroxynaphthol blue




**Figure 10** Vis absorption difference spectrum of complexed and uncomplexed hydroxynaphthol blue

### f) Methylthymol blue

3,3'-Bis[N,N-di(carboxymethyl)aminomethyl]thymolsulfonephthalein tetrasodium salt;  $C_{37}H_{40}N_2Na_4O_{13}S$   
 100 mg in 100 mL ethanol



**Figure 11** Overlaid Vis spectra of complexed and uncomplexed methylthymol blue




**Figure 12** Vis absorption difference spectrum of complexed and uncomplexed methylthymol blue


### g) Murexide

Ammonium purpurate;  $C_8H_8N_6O_6$

1% triturated in NaCl or 50 mg in 100 mL dist.  $H_2O$



**Figure 13** Overlaid Vis spectra of complexed and uncomplexed murexide




**Figure 14** Vis absorption difference spectrum of complexed and uncomplexed murexide


### h) PAN

1-(2-Pyridylazo)-2-naphthol;  $C_{15}H_{11}NO_3$

100 mg in 100 mL ethanol




**Figure 15** Overlaid Vis spectra of complexed and uncomplexed PAN




**Figure 16** Vis absorption difference spectrum of complexed and uncomplexed PAN

### i) Phthalein purple

3,3'-Bis[bis-(carboxymethyl)-aminomethyl]-cresolphthalein;  $C_{32}H_{32}N_2O_{12} \cdot H_2O$   
 100 g in 100 mL dist.  $H_2O$  containing 1 mL w/w( $NH_3$ ) = 25% (not stable)




**Figure 17** Overlaid Vis spectra of complexed and uncomplexed phthalein purple




**Figure 18** Vis absorption difference spectrum of complexed and uncomplexed phthalein purple

### j) Tiron

Pyrocatechol-3,5-disulfonic acid disodium salt;  $C_6H_4Na_2O_2S_2$   
 Added as a solid or 2% in dist.  $H_2O$



**Figure 19** Overlaid Vis spectra of complexed and uncomplexed Tiron




**Figure 20** Vis absorption difference spectrum of complexed and uncomplexed Tiron

### k) Xylenol orange


3,3'-Bis[N,N-bis(carboxymethyl)aminomethyl]-cresolsulfonephthalein tetrasodium salt;

$C_{31}H_{28}N_2Na_4O_{13}S$

100 mg are dissolved in 100 mL deion.  $H_2O$



**Figure 21** Overlaid Vis spectra of complexed and uncomplexed xylenol orange




**Figure 22** Vis absorption difference spectrum of complexed and uncomplexed xylenol orange


### l) Eriochrome Cyanine R

3,3'-(1,1-Dioxido-3H-2,1-benzoxathiol-3-ylidene)bis[6-hydroxy-5-methylbenzoic acid] trisodium salt;  $C_{23}H_{15}Na_3O_9S$

Dissolve 40 mg in 100 mL dist.  $H_2O$



**Figure 23** Overlaid Vis spectra of complexed and uncomplexed Eriochrome Cyanine R



**Figure 24** Vis absorption difference spectrum of complexed and uncomplexed Eriochrome Cyanine R

### Preparation of standard solutions

Many standard solutions are commercially available as ready-to-use products – and these should be given precedence. The manufacturer has adjusted their titer to 1.0000 at 20 °C.

#### a) Direct titration

Most complexing agents are available in a pure form as an acid. These have very poor solubility in water and must be converted into their soluble di- or trisodium salts during preparation. The substance concentration  $c$  is usually 0.1 mol/L. The solutions have a stable titer.

The preparation of  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  is described as an example:

Weigh 29.5 g EDTA ( $\text{C}_{10}\text{H}_{16}\text{N}_2\text{O}_8$ ) into a glass beaker and, while stirring, slurry in approx. 300 mL dist.  $\text{H}_2\text{O}$ . Add  $c(\text{NaOH}) = 5 \text{ mol/L}$  until everything has dissolved. Allow to cool, mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .

#### b) Back titration

$\text{Cu}^{2+}$  or  $\text{Zn}^{2+}$  solutions are usually used for back titrations. These solutions also have a stable titer and are generally prepared with a concentration of  $c(\text{Me}) = 0.1 \text{ mol/L}$ .

$c(\text{CuSO}_4) = 0.1 \text{ mol/L}$

Dissolve 25.2 g  $\text{CuSO}_4 \cdot 5 \text{ H}_2\text{O}$  (99%) in approx. 500 mL dist.  $\text{H}_2\text{O}$ . Add 0.5 mL  $w/w(\text{H}_2\text{SO}_4) = 96\%$  and then make up to 1 L with dist.  $\text{H}_2\text{O}$  and mix.

$c(\text{ZnSO}_4) = 0.1 \text{ mol/L}$

Dissolve 28.9 g  $\text{ZnSO}_4 \cdot 7 \text{ H}_2\text{O}$  (99.5%) in approx. 500 mL dist.  $\text{H}_2\text{O}$ . Add 0.5 mL  $w/w(\text{H}_2\text{SO}_4) = 96\%$  and then make up to 1 L with dist.  $\text{H}_2\text{O}$  and mix.

### Titer determinations

Calcium carbonate is the primary standard for complexing agents (CA in the following).  $\text{CaCO}_3$  is dried overnight in a drying oven at 140 °C and then left to cool in a desiccator for at least 2 h.

The titer is usually determined in triplicate, and the mean value is used as titer.

#### Titer determination for direct titration:

Weigh approx. 100 mg  $\text{CaCO}_3$  with an accuracy of 0.1 mg into a titration beaker and add 20 mL dist.  $\text{H}_2\text{O}$ . Add  $c(\text{HCl}) = 3 \text{ mol/L}$  while stirring until the  $\text{CaCO}_3$  has completely dissolved. Add approx. 80 mL dist.  $\text{H}_2\text{O}$  and 10 mL ammoniacal buffer solution, and titrate with  $c(\text{CA}) = 0.1 \text{ mol/L}$ .

### Calculation

$$\text{Titer} = \frac{m_S}{V_{EP1} \times c_{KB} \times M_{CaCO_3}}$$

$m_S$ : Weight of  $CaCO_3$  in mg

$V_{EP1}$ : Volume of CA consumed up to the first equivalence point in mL

$c_{CA}$ : Concentration of the complexing agent in mol/L

$M_{CaCO_3}$ : Molecular mass of calcium carbonate; 100.09 g/mol

### Titer determination for back titration:

To a titration beaker, add approx. 50 mL dist.  $H_2O$ , 10 mL  $c(CA) = 0.1$  mol/L, and 10 mL buffer solution (alkaline or acidic). Then titrate with  $c(Cu^{2+}) = 0.1$  mol/L or  $c(Zn^{2+}) = 0.1$  mol/L.

### Calculation


$$\text{Titer} = \frac{V_S \times f_{KB}}{V_{EP1}}$$

$V_S$ : Added volume of  $c(CA) = 0.1$  mol/L in mL


$f_{CA}$ : Titer of the CA solution

$V_{EP1}$ : Volume of  $Cu^{2+}$  or  $Zn^{2+}$  solution consumed up to the first equivalence point in mL

### Examples: Titer determination of EDTA and EGTA with $CaCO_3$



**Figure 25** Titration curve of the potentiometric titer determination of EDTA with the Ca-ISE



**Figure 26** Titration curve of the potentiometric titer determination of EGTA with the Ca-ISE

## Practical examples

### Water hardness, calcium, and magnesium

#### General

Hardness constituents in water are usually calcium and magnesium ions. They are mainly present as hydrogen carbonates and sulfates or, in rare cases, as chlorides. Barium and strontium, which are also determined in the total hardness determination, are present only in very small quantities, if at all, and thus do not interfere (0.1 mg/L  $\text{Ba}^{2+}$  would consume only 0.7  $\mu\text{L}$   $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  for a sample size of 100 mL).

Water hardness is given in mmol/L. Terms such as "soft", "hard", " $^\circ\text{dH}$ ", " $^\circ\text{fH}$ " etc. are also used for information purposes, for example, on detergent packaging. The relationships are shown in the following table:

mmol/L	mg/L $\text{CaCO}_3$	$^\circ\text{fH}$	$^\circ\text{dH}$	$^\circ\text{USH}$	Hardness
0–0.7	0–70	0–7	0–3.92	0–3.77	very soft
0.7–1.5	70–150	7–15	3.92–8.4	3.77–8.07	soft
1.5–2.5	150–250	15–25	8.4–14.0	8.07–9.50	moderately hard
2.5–3.2	250–320	25–32	14.0–17.92	9.50–17.22	fairly hard
3.2–4.2	320–420	32–42	17.92–23.5	17.22–22.6	hard
> 4.2	> 420	> 42	> 23.5	> 22.6	very hard

#### Conversion factors

mmol/L  $\times$  100  $\rightarrow$  mg/L  $\text{CaCO}_3$

mmol/L  $\times$  10  $\rightarrow$   $^\circ\text{fH}$  (French degree of hardness)

mmol/L  $\times$  5.6  $\rightarrow$   $^\circ\text{dH}$  (German degree of hardness)

mmol/L  $\times$  5.38  $\rightarrow$   $^\circ\text{USH}$  (US degree of hardness)

If the water being investigated has an unfavorable  $\text{Ca}^{2+}/\text{Mg}^{2+}$  ratio (small amount of  $\text{Mg}^{2+}$  and a large amount of  $\text{Ca}^{2+}$ ), this

leads to errors in all titrimetric methods. In such cases, the recommendation is to use a buffer solution with an accurately added amount of  $\text{Mg}^{2+}$  standard (e.g., 2.4305 mg  $\text{Mg}^{2+}$ , which corresponds to 1 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$ ) and then to subtract this value during the calculation.

In the case of soft water, the best method is to titrate with  $c(\text{Na}_2\text{EDTA}) = 0.05 \text{ mol/L}$ .



## Potentiometric titration with the Ca ISE

### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$   
 $c(\text{KOH}) = 0.1 \text{ mol/L}$
- Auxiliary complex solution: weigh 24.3 g tris(hydroxymethyl)amino-methane (TRIS) into a 1 L volumetric flask and dissolve in approx. 500 mL dist.  $\text{H}_2\text{O}$ . Add 10 mL acetylacetone, mix, and make up to the mark with dist.  $\text{H}_2\text{O}$ . Allow the solution to stand for 24 h before use (max. shelf life is 1 week).

### Analyses

Measure a 100 mL water sample into a titration beaker and add 15 mL auxiliary complex solution. Using the "DET U" mode, titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the second equivalence point ( $\text{EP}_2$ ).  $\text{EP}_1$  corresponds to  $\text{Ca}^{2+}$  and the difference  $\text{EP}_2 - \text{EP}_1$  corresponds to  $\text{Mg}^{2+}$ .

### Calculation

$$\text{Calcium hardness} = \frac{V_{\text{EP1}} \times c_{\text{EDTA}} \times f \times 1000}{V_S}$$

$$\text{Magnesium hardness} = \frac{(V_{\text{EP2}} - V_{\text{EP1}}) \times c_{\text{EDTA}} \times f \times 1000}{V_S}$$


$$\text{Total hardness} = \frac{V_{\text{EP2}} \times c_{\text{EDTA}} \times f \times 1000}{V_S}$$

$$\text{Calcium content} = \frac{V_{\text{EP1}} \times c_{\text{EDTA}} \times f \times 1000 \times M_{\text{Ca}}}{V_S}$$

$$\text{Magnesium cont.} = \frac{(V_{\text{EP2}} - V_{\text{EP1}}) \times c_{\text{EDTA}} \times f \times 1000 \times M_{\text{Mg}}}{V_S}$$

- Calcium hardness: Calcium content in mmol/L
- Magnesium hardn.: Magnesium content in mmol/L
- Total hardness: Total hardness in mmol/L
- Calcium content: Calcium content in mg/L
- Magnesium cont.: Magnesium content in mg/L

- $V_{\text{EP1}}$ : Volume of EDTA consumed up to the first equivalence point in mL
- $V_{\text{EP2}}$ : Volume of EDTA consumed up to the second equival. point in mL
- $c_{\text{EDTA}}$ : Concentration of EDTA in mol/L
- $f$ : Titer of the EDTA solution
- 1000: Conversion factor
- $V_S$ : Sample size in mL
- $M_{\text{Ca}}$ : Molar mass of calcium; 40.078 g/mol
- $M_{\text{Mg}}$ : Molar mass of magnesium; 24.305 g/mol



**Figure 27** Titration curve of tap water from Herisau (Switzerland), measured with the Ca-ISE

**Further literature on this topic can be accessed using the Application Finder under**

[www.metrohm.com/applications](http://www.metrohm.com/applications)

- *Application Bulletin AB-125 – Simultaneous determination of calcium, magnesium, and alkalinity by complexometric titration with potentiometric or photometric indication in water and beverage samples*

### Potentiometric with the Cu ISE Reagents

- Titrant 1:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Titrant 2:  $c(\text{Na}_2\text{EGTA}) = 0.1 \text{ mol/L}$
- Buffer solution: dissolve 54 g  $\text{NH}_4\text{Cl}$  and 350 mL w/w( $\text{NH}_3$ ) = 25% in dist.  $\text{H}_2\text{O}$ , mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Cu complex solution:  $c[\text{Cu}(\text{NH}_4)_2\text{EDTA}] = 0.1 \text{ mol/L}$  (e.g., Merck No. 105217)
- CuEGTA: dissolve 2.497 g  $\text{CuSO}_4 \cdot 5 \text{ H}_2\text{O}$  in dist.  $\text{H}_2\text{O}$ , make up to 100 mL and mix well. Add 25.0 mL of this solution to 25.0 mL  $c(\text{Na}_2\text{EGTA}) = 0.1 \text{ mol/L}$ .

### Analysis

#### Total hardness (sum of $\text{Ca}^{2+}$ & $\text{Mg}^{2+}$ )

Measure a 100 mL water sample into a titration beaker and add 5 mL buffer solution. Add 0.5 mL Cu complex solution, wait for 20 s while stirring, and titrate in the "MET U" mode with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$ . There may be two equivalence points, whereby  $\text{EP}_2$  corresponds to the sum  $\text{Ca}^{2+} + \text{Mg}^{2+}$  ( $\text{EP}_1$  cannot be used for the calculations

because it would give inaccurate results). Since only the largest equivalence point is required, select "largest" as the EP criterion so that only one EP is shown.

#### Calcium (hardness)

Measure a 100 mL water sample into a titration beaker and add 5 mL buffer solution. Add 0.5 mL CuEGTA, wait for 20 s while stirring, and titrate in the "MET U" mode with  $c(\text{Na}_2\text{EGTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.

#### Magnesium (hardness)

This is calculated by subtracting the calcium hardness from the total hardness.

#### Calculation

$$\text{Total hardness} = \frac{V_{\text{EP}1.1} \times c_{\text{EDTA}} \times f_{\text{EDTA}} \times 1000}{V_S}$$

$$\text{Calcium hardn.} = \frac{V_{\text{EP}1.2} \times c_{\text{EGTA}} \times f_{\text{EGTA}} \times 1000}{V_S}$$

$$\text{Magnesium hardn.} = \text{total hardness} - \text{calcium hardness}$$

Total hardness: Total hardness in mmol/L

Calcium hardness: Calcium content in mmol/L

Magnesium hardn.: Magnesium content in mmol/L

$V_{\text{EP}1.1}$ : Volume of EDTA consumed up to the largest equivalence point of the total hardness titration in mL


$c_{\text{EDTA}}$ : Concentration of EDTA in mol/L  
 $f_{\text{EDTA}}$ : Titer of the EDTA solution

- 1000: Conversion factor  
 $V_S$ : Sample size in mL  
 $V_{EP1,2}$ : Volume of EGTA consumed up to the first equivalence point of the calcium hardness titration in mL  
 $c_{EGTA}$ : Concentration of EGTA in mol/L  
 $f_{EGTA}$ : Titer of the EGTA solution


### Further literature on this topic is given in

- *Application Bulletin AB-101 – Complexometric titrations with the Cu-ISE*
- *Titration Application Note AN-T-131 – Automatic determination of calcium, magnesium, and total hardness in water with the Cu-ISE*

### Examples: Titration curves of Herisau tap water with the Cu ISE



**Figure 28** Titration curve of the total hardness determination



**Figure 29** Titration curve of the calcium hardness determination

### Photometric with the Optrode at 610 nm

#### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Buffer solution: Dissolve 54 g  $\text{NH}_4\text{Cl}$  in dist. water and add 350 mL  $w/w(\text{NH}_3) = 25\%$ . Mix the solution and make up to 1 L with dist. water.
- Sodium hydroxide:  $c(\text{NaOH}) = 2 \text{ mol/L}$
- $\text{Na}_2\text{MgEDTA} \cdot \text{H}_2\text{O}$
- Indicator Solution 1: Eriochrome Black T; dissolve 100 mg Eriochrome Black T in 100 mL ethanol.
- Indicator Solution 2: HHSNN; dissolve 50 mg HHSNN in 100 mL  $c(\text{NaOH}) = 0.1 \text{ mol/L}$ .

## General remarks

There are a few special requirements that must be taken into account when carrying out photometric titration with optical sensors:

- No interfering gas bubbles should be present in the light path. For this reason, we recommend degassing the water samples before adding the buffer and adjusting the stirrer speed so that no gas bubbles are entrained.
- During determination of the calcium hardness, precipitation of magnesium hydroxide leads to ever increasing turbidity. The mixture must therefore be titrated immediately after adding NaOH and preferably with a start volume. Although the resulting titration curves are not very smooth and look unattractive, this does not affect the accuracy of the results.

## Analysis

### Total hardness (sum of $\text{Ca}^{2+}$ & $\text{Mg}^{2+}$ )

Measure a 100 mL water sample into a titration beaker, add 0.1 g  $\text{Na}_2\text{MgEDTA} \cdot \text{H}_2\text{O}$ , 10 mL buffer solution, and 0.25 mL Indicator Solution 1, and titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  in the "MET U" mode past the color change.

### Calcium (hardness)

Measure a 100 mL water sample into the titration beaker, add 2 mL  $c(\text{NaOH}) = 2 \text{ mol/L}$  ( $\text{pH} > 12$ ) and 1.5 mL Indicator Solution 2, and titrate immediately with

$c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  in the "MET U" mode up to the color change.

## Magnesium (hardness)

This is calculated by subtracting the calcium hardness from the total hardness.

## Calculation

$$\text{Total hardness} = \frac{V_{\text{EP1.1}} \times c_{\text{EDTA}} \times f \times 1000}{V_S}$$

$$\text{Calcium hardn.} = \frac{V_{\text{EP1.2}} \times c_{\text{EDTA}} \times f \times 1000}{V_S}$$

$$\text{Magnesium hardn.} = \text{total hardness} - \text{calcium hardness}$$

Total hardness: Total hardness in mmol/L

Calcium hardness: Calcium content in mmol/L

Magnesium hardn.: Magnesium content in mmol/L

$V_{\text{EP1.1}}$ : Volume of EDTA consumed up to the first equivalence point of the total hardness titration in mL

$c_{\text{EDTA}}$ : Concentration of EDTA in mol/L


f: Titer of the EDTA solution

1000: Conversion factor


$V_S$ : Sample size in mL

$V_{\text{EP1.2}}$ : Volume of EDTA consumed up to the first equivalence point of the calcium hardness titration in mL

## Examples: Photometric titration curves of Herisau tap water



**Figure 30** Titration curve of the total hardness determination with the Optrode at 610 nm



**Figure 31** Titration curve of the calcium hardness determination with the Optrode at 610 nm

### Further literature on this topic is given in

- *Application Bulletin AB-125 – Simultaneous determination of calcium, magnesium, and alkalinity by complexometric titration with potentiometric or photometric indication in water and beverage samples*
- *Titration Application Note AN-T-084 – Fully automatic determination of the total, calcium and magnesium hardness of water samples using photometric titration*

## Aluminum and magnesium, alloys and antacids

### Potentiometric with the Cu ISE

#### Reagents

- Titrant 1:  $c(\text{DCTA}) = 0.1 \text{ mol/L}$ ; dissolve 36.463 g *trans*-1,2-diaminocyclohexanetetraacetic acid monohydrate in 400 mL  $c(\text{NaOH}) = 0.5 \text{ mol/L}$ . Allow to cool, then mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Titrant 2:  $c(\text{CuSO}_4) = 0.1 \text{ mol/L}$ ; dissolve 24.968 g  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$  in approx. 500 mL dist.  $\text{H}_2\text{O}$ .
- Add 0.5 mL  $w/w(\text{H}_2\text{SO}_4) = 96\%$ , mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Buffer solution pH = 4.7: Dissolve 123 g sodium acetate in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 86 mL glacial acetic acid. Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Ammonia solution:  $w/w(\text{NH}_3) = 10\%$
- Sodium hydroxide:  $c(\text{NaOH}) = 2 \text{ mol/L}$

### Sample preparation

Preneutralize the acidic sample solution, which should not contain more than 12 mg Al<sup>3+</sup> and/or 20 mg Mg<sup>2+</sup>, with c(NaOH) = 2 mol/L to pH 2–3, and then dilute with dist. H<sub>2</sub>O to approx. 50 mL.

### Analysis

To the prepared sample solution, add 5 mL buffer solution pH 4.7 and 6 mL c(DCTA) = 0.1 mol/L. Allow the mixture to react for 1 min while stirring. Next, titrate with

c(CuSO<sub>4</sub>) = 0.1 mol/L past the first equivalence point → aluminum content. Store the value of the excess titrant of Titrant 2 as a "common variable".

To the completely titrated sample solution, add 20 mL ammonia solution and a further 6 mL c(DCTA) = 0.1 mol/L. Next, titrate with c(CuSO<sub>4</sub>) = 0.1 mol/L past the first equivalence point → magnesium content.

### Calculation

$$\beta_{\text{Al}} = \frac{(V_{\text{DCTA}} \times f_1 \times c_{\text{DCTA}} - V_{\text{EP1.1}} \times f_2 \times c_{\text{CuSO}_4}) \times M_{\text{Al}}}{V_{\text{S}}}$$

$$V_{\text{Ex}} = V_{\text{End}} - V_{\text{EP1.1}}$$

$$\beta_{\text{Mg}} = \frac{(V_{\text{DCTA}} \times f_1 \times c_{\text{DCTA}} - (V_{\text{EP1.2}} + V_{\text{Ex}}) \times f_2 \times c_{\text{CuSO}_4}) \times M_{\text{Mg}}}{V_{\text{S}}}$$

$\beta_{\text{Al}}$ : Aluminum content of the sample in g/L

$\beta_{\text{Mg}}$ : Magnesium content of the sample in g/L

$V_{\text{DCTA}}$ : Added quantity of DCTA standard solution in mL

$V_{\text{EP1.1}}$ : Volume of CuSO<sub>4</sub> standard solution consumed up to the first equivalence point of the aluminum titration in mL

$V_{\text{EP1.2}}$ : Volume of CuSO<sub>4</sub> standard solution consumed up to the first equivalence point of the magnesium titration in mL

$V_{\text{Ex}}$ : Excess quantity of titrant after the aluminum titration in mL

$V_{\text{End}}$ : End volume of the aluminum titration in mL

$f_1$ : Titer of the DCTA standard solution

$f_2$ : Titer of the CuSO<sub>4</sub> standard solution

$c_{\text{DCTA}}$ : Concentration of the DCTA standard solution in mol/L

$c_{\text{CuSO}_4}$ : Concentration of the CuSO<sub>4</sub> standard solution in mol/L

$M_{\text{Al}}$ : Molar mass of aluminum; 26.982 g/mol


$M_{\text{Mg}}$ : Molar mass of magnesium; 24.305 g/mol

$V_{\text{S}}$ : Sample size in mL


### Remarks

The surface of the Cu ISE must be polished occasionally with aluminum oxide power (e.g., polishing set 6.2802.000).

### Example: $\text{Al}^{3+}$ and $\text{Mg}^{2+}$ in aqueous solution



**Figure 32** Titration curve of the potentiometric aluminum determination with the Cu-ISE



**Figure 33** Titration curve of the potentiometric magnesium determination with the Cu-ISE

#### Further literature on this topic is given in

- *Application Bulletin AB-181 – Automated potentiometric titration of aluminum and magnesium in the same solution*
- *Titration Application Note No. AN-T-117 – Automatic determination of aluminum and magnesium mixtures with ion-selective copper electrode (Cu-ISE)*

## Magnesium and zinc

### Potentiometric with the Cu ISE

#### Reagents

- Titrant 1:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Titrant 2:  $c(\text{CuSO}_4) = 0.1 \text{ mol/L}$ ; dissolve 24.968 g  $\text{CuSO}_4 \cdot 5 \text{ H}_2\text{O}$  in approx. 500 mL dist.  $\text{H}_2\text{O}$ . Add 0.5 mL  $w/w(\text{H}_2\text{SO}_4) = 96\%$ , mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Acetate buffer,  $\text{pH} = 4.7$ : Dissolve 123 g sodium acetate in approx. 500 mL dist.  $\text{H}_2\text{O}$  and add 86 mL glacial acetic acid. Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Ammonia buffer,  $\text{pH} = 10$ : Dissolve 54 g  $\text{NH}_4\text{Cl}$  in approx. 500 mL dist.  $\text{H}_2\text{O}$  and add 350 mL  $w/w(\text{NH}_3) = 25\%$ . Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Sodium hydroxide:  $c(\text{NaOH}) = 1 \text{ mol/L}$

### Sample preparation

If necessary, dilute the sample solution, which should not contain more than 50 mg Zn<sup>2+</sup> and/or 20 mg Mg<sup>2+</sup>, with dist. H<sub>2</sub>O to approx. 50 mL and adjust to pH 4–4.5 with c(NaOH) = 1 mol/L.

### Analysis

To the prepared sample solution, add 5 mL acetate buffer and 10 mL c(DCTA) = 0.1 mol/L, and allow the mixture to react for 1 min while stirring. Next, titrate with c(CuSO<sub>4</sub>) = 0.1 mol/L past the first equivalence

point → zinc content. Store the value of the excess titrant of Titrant 2 as a "common variable".

To the completely titrated sample solution, add 10 mL ammonia buffer and a further 10 mL c(DCTA) = 0.1 mol/L. Next, titrate with c(CuSO<sub>4</sub>) = 0.1 mol/L past the first equivalence point → magnesium content.

### Calculation

$$\beta_{\text{Zn}} = \frac{(V_{\text{DCTA}} \times f_1 \times c_{\text{DCTA}} - V_{\text{EP1.1}} \times f_2 \times c_{\text{CuSO}_4}) \times M_{\text{Zn}}}{V_{\text{S}}}$$

$$V_{\text{Ex}} = V_{\text{End}} - V_{\text{EP1.1}}$$

$$\beta_{\text{Mg}} = \frac{(V_{\text{DCTA}} \times f_1 \times c_{\text{DCTA}} - (V_{\text{EP1.2}} + V_{\text{Ex}}) \times f_2 \times c_{\text{CuSO}_4}) \times M_{\text{Mg}}}{V_{\text{S}}}$$

$\beta_{\text{Zn}}$ : Zinc content of the sample in g/L

$\beta_{\text{Mg}}$ : Magnesium content of the sample in g/L

$V_{\text{DCTA}}$ : Added quantity of DCTA standard solution in mL

$V_{\text{EP1.1}}$ : Volume of CuSO<sub>4</sub> standard solution consumed up to the first equivalence point of the zinc titration in mL

$V_{\text{EP1.2}}$ : Volume of CuSO<sub>4</sub> standard solution consumed up to the first equivalence point of the magnesium titration in mL

$V_{\text{Ex}}$ : Excess quantity of titrant after the zinc titration in mL

$V_{\text{End}}$ : End volume of the zinc titration in mL

$f_1$ : Titer of the DCTA standard solution

$f_2$ : Titer of the CuSO<sub>4</sub> standard solution

$c_{\text{DCTA}}$ : Concentration of the DCTA standard solution in mol/L

$c_{\text{CuSO}_4}$ : Concentration of the CuSO<sub>4</sub> standard solution in mol/L

$M_{\text{Zn}}$ : Molar mass of zinc;  
65.409 g/mol

$M_{\text{Mg}}$ : Molar mass of magnesium;  
24.305 g/mol


$V_{\text{S}}$ : Sample size in mL

### Remarks


The surface of the Cu-ISE must be polished occasionally with aluminum oxide powder (e.g., polishing set 6.2802.000).



### Example: $\text{Zn}^{2+}$ and $\text{Mg}^{2+}$ in aqueous solution



**Figure 34** Titration curve of the potentiometric zinc determination with the Cu-ISE



**Figure 35** Titration curve of the potentiometric magnesium determination with the Cu-ISE

#### Further literature on this topic is given in

- *Titration Application Note AN-T-118 – Automatic determination of zinc and magnesium mixtures with ion-selective copper electrode (Cu-ISE)*

## Aluminum, calcium, iron, and magnesium in cement and clinker

### Photometric with the Optrode at 610 nm

#### Reagents

- Hydrochloric acid: w/w(HCl) = 36%
- Nitric acid: w/w(HNO<sub>3</sub>) = approx. 65%
- Sodium hydroxide: w/w(NaOH) > 99%
- Ammonia: w/w(NH<sub>3</sub>) = 25%
- Ammonium chloride: NH<sub>4</sub>Cl p.a.
- Titrant 1 for Ca<sup>2+</sup> and Mg<sup>2+</sup>: c(Na<sub>2</sub>EDTA) = 0.1 mol/L
- Titrant 2 for Fe<sup>3+</sup>: c(Na<sub>2</sub>EDTA) = 0.025 mol/L
- Titrant 3 for Al<sup>3+</sup>: c(Bi(NO<sub>3</sub>)<sub>3</sub>) = 0.05 mol/L; dissolve 24.25 g Bi(NO<sub>3</sub>)<sub>3</sub>·5H<sub>2</sub>O in approx. 500 mL HNO<sub>3</sub> (2 mol/L) and transfer into a 1 L volumetric flask. Next, make up to 1 L with dist. H<sub>2</sub>O.

- Ammonia:  $c(\text{NH}_3) = 2 \text{ mol/L}$ ; add 144 mL  $w/w(\text{NH}_3) = 25\%$  to a 1 L volumetric flask and make up to the mark with dist.  $\text{H}_2\text{O}$ .
- Sodium hydroxide:  $c(\text{NaOH}) = 2 \text{ mol/L}$ ; dissolve 80 g NaOH in approx. 600 mL dist.  $\text{H}_2\text{O}$ , transfer into a 1 L volumetric flask, and make up to the mark with dist.  $\text{H}_2\text{O}$ .
- Hydrochloric acid:  $c(\text{HCl}) = 6 \text{ mol/L}$ ; add 590 mL  $w/w(\text{HCl}) = 37\%$  to a 1 L volumetric flask that already contains approx. 200 mL dist.  $\text{H}_2\text{O}$ . After the solution has cooled, make up to the mark with dist.  $\text{H}_2\text{O}$ .
- Nitric acid:  $c(\text{HNO}_3) = 2 \text{ mol/L}$ ; add 192 mL  $w/w(\text{HNO}_3) = 65\%$  to a 1 L volumetric flask that already contains approx. 500 mL dist.  $\text{H}_2\text{O}$ . After the solution has cooled, make up to the mark with dist.  $\text{H}_2\text{O}$ .
- Acetate buffer: Dissolve 60 g ammonium acetate in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 400 mL glacial acetic acid. Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Ammonia buffer: Dissolve 54 g  $\text{NH}_4\text{Cl}$  in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 350 mL  $w/w(\text{NH}_3) = 25\%$ . Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Indicator for  $\text{Ca}^{2+}$ , murexide: triturate 100 mg murexide with 10 g NaCl.
- Indicator for  $\text{Mg}^{2+}$ , methylthymol blue: triturate 100 mg methylthymol blue with 10 g NaCl.
- Indicator solution for  $\text{Fe}^{3+}$ , sulfosalicylic acid: dissolve 4 g sulfosalicylic acid in 100 mL dist.  $\text{H}_2\text{O}$ .

- Indicator solution for  $\text{Al}^{3+}$ , xylenol orange: dissolve 100 mg xylenol orange disodium salt in 100 mL dist.  $\text{H}_2\text{O}$ .

### Sample preparation

Mix approx. 4 g cement with 4 g  $\text{NH}_4\text{Cl}$ , add 48 mL  $c(\text{HCl}) = 6 \text{ mol/L}$  and 3 mL  $c(\text{HNO}_3) = 2 \text{ mol/L}$ . Heat the mixture and boil for a while. Stir the mixture with a magnetic stirrer for approx. 60 min. Take up the mixture in 50 mL hot dist.  $\text{H}_2\text{O}$  and filter it through a black band filter into a 500 mL volumetric flask. Rinse the filter thoroughly with hot dist.  $\text{H}_2\text{O}$ . Allow to cool, then mix the combined filtrates and make up to the mark with dist.  $\text{H}_2\text{O}$ .

### Analyses

Carry out all analyses with the Optrode at a wavelength of 610 nm.

### Calcium

To a 2.5 mL sample solution, add approx. 70 mL dist.  $\text{H}_2\text{O}$  and adjust to pH 12 with  $c(\text{NaOH}) = 2 \text{ mol/L}$ . Add a spatula tip of murexide indicator, then titrate the solution with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the second breakpoint. Store the titrant consumption up to  $\text{BP}_2$  as a "common variable" because it will be used in the determination of the  $\text{Mg}^{2+}$  content.

### Calculation

$$\text{CaO \%} = \frac{V_{\text{BP2}} \times C_{\text{EDTA}} \times f_{\text{EDTA}} \times 40 \times M_{\text{CaO}}}{m_s}$$

- CaO%: Calcium content as calcium oxide in %  
 $V_{\text{BP2}}$ : Volume of EDTA consumed up to the second breakpoint in mL  
 $C_{\text{EDTA}}$ : Concentration of EDTA in mol/L;  $c(\text{EDTA}) = 0.1 \text{ mol/L}$   
 $f_{\text{EDTA}}$ : Titer of the EDTA solution  
 40: Conversion factor  
 $M_{\text{CaO}}$ : Molar mass of calcium oxide; 56.08 g/mol  
 $m_s$ : Sample size in g

### Magnesium

This titration is used to determine the sum of  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$ , and to calculate the magnesium content by subtracting the titrant consumption for the  $\text{Ca}^{2+}$  determination from the total consumption.

Pipette 2.5 mL of the digested solution into a titration beaker and dilute to approx. 80 mL with dist.  $\text{H}_2\text{O}$ . Adjust the solution to pH 10 with  $c(\text{NaOH}) = 2 \text{ mol/L}$ , and add a spatula tip of methylthymol blue. Next, titrate the solution with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.

### Calculation

$$\text{MgO \%} = \frac{(V_{\text{EP1}} - V_{\text{Ca, BP2}}) \times C_{\text{EDTA}} \times f_{\text{EDTA}} \times 40 \times M_{\text{MgO}}}{m_s}$$

- MgO%: Magnesium content as magnesium oxide in %  
 $V_{\text{EP1}}$ : Volume of EDTA consumed up to the first equivalence point in mL  
 $V_{\text{Ca, BP2}}$ : Volume of EDTA consumed in the calcium determination in mL  
 $C_{\text{EDTA}}$ : Concentration of EDTA in mol/L;  $c(\text{EDTA}) = 0.1 \text{ mol/L}$   
 $f_{\text{EDTA}}$ : Titer of the EDTA solution  
 40: Conversion factor  
 $M_{\text{MgO}}$ : Molar mass of magnesium oxide; 40.32 g/mol  
 $m_s$ : Sample size in g

### Iron

Pipette 10 mL of the digested solution into a titration beaker, dilute with dist.  $\text{H}_2\text{O}$ , and add 1 mL sulfosalicylic acid indicator. Adjust the mixture to pH 1.5–2 with  $c(\text{NH}_3) = 2 \text{ mol/L}$  or with  $c(\text{HCl}) = 6 \text{ mol/L}$ , as appropriate. Next, titrate with  $c(\text{Na}_2\text{EDTA}) = 0.025 \text{ mol/L}$  past the first breakpoint. (Use small volume increments – e.g., 25  $\mu\text{L}$ ).

Store the titrant consumption up to  $\text{BP}_1$  as a "common variable" because it will be used in the determination of  $\text{Al}^{3+}$ .

### Calculation

$$\text{Fe}_2\text{O}_3 \% = \frac{V_{\text{EP1}} \times c_{\text{EDTA}} \times f_{\text{EDTA}} \times 10 \times M_{\text{Fe}_2\text{O}_3}}{m_s \times 2}$$

- Fe<sub>2</sub>O<sub>3</sub> %: Iron content as iron oxide in %  
 V<sub>EP1</sub>: Volume of EDTA consumed up to the first equivalence point in mL  
 c<sub>EDTA</sub>: Concentration of EDTA in mol/L; c(EDTA) = 0.025 mol/L  
 f<sub>EDTA</sub>: Titer of the EDTA solution  
 10: Conversion factor  
 M<sub>Fe<sub>2</sub>O<sub>3</sub></sub>: Molar mass of iron oxide; 159.69 g/mol  
 m<sub>s</sub>: Sample size in g  
 2: Stoichiometric factor

### Aluminum

As in the Mg<sup>2+</sup> determination, this is also a cumulative determination method. Al<sup>3+</sup> and Fe<sup>3+</sup> are detected together. The aluminum content is then calculated by subtracting the titrant consumption for the Fe<sup>3+</sup> determination from the total consumption.


Pipette 10 mL sample into a glass beaker and add 70 mL dist. H<sub>2</sub>O. Add 10 mL acetate buffer, adjust to pH 3.5 with c(HCl) = 6 mol/L, and add 1.25 mL c(Na<sub>2</sub>EDTA) = 0.1 mol/L and 1 mL xylenol orange indicator solution. Titrate the solution with c(Bi(NO<sub>3</sub>)<sub>3</sub>) = 0.05 mol/L past the first equivalence point.

### Calculation


$$\text{Al}_2\text{O}_3 \% = \frac{((V_{\text{EDTA1}} \times c_{\text{EDTA1}} \times f_{\text{EDTA1}}) - (V_{\text{EDTA2}} \times c_{\text{EDTA2}} \times f_{\text{EDTA2}}) - (V_{\text{EP1}} \times c_{\text{Bi}} \times f_{\text{Bi}})) \times M_{\text{Al}_2\text{O}_3} \times 10}{m_s \times 2}$$

- |  |   |
|--|---|
| Al <sub>2</sub> O <sub>3</sub> %: Aluminum content as aluminum oxide in %                                      | f <sub>EDTA2</sub> : Titer of the Na <sub>2</sub> EDTA solution<br>c(Na <sub>2</sub> EDTA) = 0.025 mol/L  |
| V <sub>EDTA1</sub> : Initial volume of c(Na <sub>2</sub> EDTA) = 0.1 mol/L in mL                               | V <sub>EP1</sub> : Volume of c(Bi(NO <sub>3</sub> ) <sub>3</sub> ) = 0.05 mol/L consumed up to the first equivalence point in mL                  |
| c <sub>EDTA1</sub> : Concentration of the Na <sub>2</sub> EDTA solution in mol/L; c(EDTA) = 0.1 mol/L          | c <sub>Bi</sub> : Concentration of the Bi(NO <sub>3</sub> ) <sub>3</sub> solution in mol/L;<br>c(Bi(NO <sub>3</sub> ) <sub>3</sub> ) = 0.05 mol/L |
| f <sub>EDTA1</sub> : Titer of the Na <sub>2</sub> EDTA solution<br>c(Na <sub>2</sub> EDTA) = 0.1 mol/L         | f <sub>Bi</sub> : Titer of the Bi(NO <sub>3</sub> ) <sub>3</sub> solution;<br>c(Bi(NO <sub>3</sub> ) <sub>3</sub> ) = 0.05 mol/L                  |
| V <sub>EDTA2</sub> : Volume of c(Na <sub>2</sub> EDTA) = 0.025 mol/L consumed for the iron determination in mL | 10: Conversion factor   |
| c <sub>EDTA2</sub> : Concentration of the Na <sub>2</sub> EDTA solution in mol/L;<br>c(EDTA) = 0.025 mol/L     | M <sub>Al<sub>2</sub>O<sub>3</sub></sub> : Molar mass of aluminum oxide; 101.96 g/mol   |
|  | m <sub>s</sub> : Sample size in g   |
|  | 2: Stoichiometric factor  |


## Examples: Titration curves of a cement sample




**Figure 36** Titration curve of the photometric determination of calcium in cement with the Optrode at 610 nm



**Figure 37** Titration curve of the photometric determination of magnesium in cement with the Optrode at 610 nm



**Figure 38** Titration curve of the photometric determination of iron in cement with the Optrode at 610 nm



**Figure 39** Titration curve of the photometric determination of aluminum in cement with the Optrode at 610 nm

### Further literature on this topic is given in

- Application Bulletin AB-063 – Determination of silicon, calcium, magnesium, iron, and aluminum in digested cement samples by photometric titrations
- Titration Application Note AN-T-078 – Determination of aluminum in cement using photometric titration
- Titration Application Note AN-T-079 – Determination of calcium in cement by photometric titration of the solubilized product according to EN 196-2
- Titration Application Note AN-T-080 – Determination of iron in cement using photometric titration
- Titration Application Note AN-T-081 – Determination of magnesium in cement using photometric titration

# Sulfate

38

## Potentiometric with the Ca ISE

### Reagents

- Titrant:  $c(\text{Na}_2\text{EGTA}) = 0.05 \text{ mol/L}$ ; slurry 19.4 g ethylene glycol-bis(2-aminoethyl-ether)tetraacetic acid ( $w = 98\%$ ) in approx. 200 mL dist.  $\text{H}_2\text{O}$ . Add  $c(\text{NaOH}) = 10 \text{ mol/L}$  while stirring until everything has dissolved. Allow to cool, mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Barium chloride solution:  $c(\text{BaCl}_2) = 0.05 \text{ mol/L}$ ; dissolve 12.34 g  $\text{BaCl}_2 \cdot 2\text{H}_2\text{O}$  ( $w = 99\%$ ) in  $c(\text{HCl}) = 0.01 \text{ mol/L}$ , mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Buffer solution pH = 10: dissolve 9 g  $\text{NH}_4\text{Cl}$  and 60 mL  $w/w(\text{NH}_3) = 25\%$  in dist.  $\text{H}_2\text{O}$ , mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- $\text{Ca}^{2+}$  standard:  $c(\text{Ca}^{2+}) = 0.1 \text{ mol/L}$ ; e.g., Metrohm No. 6.2301.070
- Nitric acid:  $c(\text{HNO}_3) = 1 \text{ mol/L}$

### Blank value determination

To a glass beaker, add approx. 50 mL dist.  $\text{H}_2\text{O}$ , 1 mL  $c(\text{HNO}_3) = 1 \text{ mol/L}$ , 0.5 mL  $\text{Ca}^{2+}$  standard ( $c(\text{Ca}^{2+}) = 0.1 \text{ mol/L}$ ), and 7.5 mL  $c(\text{BaCl}_2) = 0.05 \text{ mol/L}$ . Allow to react for 3 min while stirring. Add 5 mL buffer solution pH 10, allow to react for a further 30 s, and titrate with  $c(\text{Na}_2\text{EGTA}) = 0.05 \text{ mol/L}$  past the second equivalence point.  $\text{EP}_1$  corresponds to the calcium content and the difference  $\text{EP}_2 - \text{EP}_1$  to the barium content. Store  $\text{EP}_2 - \text{EP}_1$  (mL) as blank value (common variable).

### Analysis

If necessary, dilute the sample solution, which should not contain more than 20 mg sulfate, with dist.  $\text{H}_2\text{O}$  to approx. 50 mL and adjust to  $\text{pH} < 4$  with  $c(\text{HNO}_3) = 1 \text{ mol/L}$ . Add 0.5 mL  $\text{Ca}^{2+}$  standard ( $c(\text{Ca}^{2+}) = 0.1 \text{ mol/L}$ ) and 7.5 mL  $c(\text{BaCl}_2) = 0.05 \text{ mol/L}$ , and allow to react for 3 min while stirring. Add 5 mL buffer solution pH 10, allow to react for a further 30 s, and titrate with  $c(\text{Na}_2\text{EGTA}) = 0.05 \text{ mol/L}$  past the second equivalence point.

### Calculation

$$\beta_{\text{SO}_4^{2-}} = \frac{(\text{Blank} - (V_{\text{EP}_2} - V_{\text{EP}_1})) \times c_{\text{EGTA}} \times f \times M_{\text{SO}_4^{2-}}}{V_S}$$

$\beta_{\text{SO}_4^{2-}}$ : Sulfate content of the sample in g/L

Blank: Blank value in mL

$V_{\text{EP}_2}$ : Volume of EGTA standard solution consumed up to the second equivalence point in mL

$V_{\text{EP}_1}$ : Volume of EGTA standard solution consumed up to the first equivalence point in mL


$c_{\text{EGTA}}$ : Concentration of the EGTA standard solution in mol/L

f: Titer of the EGTA standard solution


$M_{\text{SO}_4^{2-}}$ : Molar mass of sulfate; 96.063 g/mol

$V_S$ : Sample size in mL

### Example: Sulfate in aqueous solution



**Figure 40** Titration curve of the blank value determination with the Ca ISE



**Figure 41** Titration curve of the sulfate determination with the Ca ISE

#### Further literature on this topic is given in

- *Application Bulletin AB-140 – Titrimetric sulfate determination*
- *Titration Application Note AN-T-116 – Automatic sulfate determination in aqueous solution using a combined ion-selective calcium electrode (Ca-ISE)*

## Single determinations

### Al – Aluminum

#### General

The complexometric titration of aluminum is associated with a particularities that should be taken into consideration when determining this metal.

#### Reaction rate

Complexation is slow at room temperature. For this reason, many titration

procedures were formerly carried out in a hot solution. However, if the titrant excess is sufficiently large (> 50%), the reaction is fast and the titration can be performed at room temperature. Hence,  $\text{Al}^{3+}$  can be determined by back titration.

### Formation of hydroxo complexes

Al forms hydroxo complexes that react very sluggishly (if at all) with the complexing agent. Al must thus be titrated in acidic solutions. If the preneutralization is carried out by adding hydroxide ions, for example, this must be taken into account and the solution should only be neutralized to  $\text{pH} \approx 4$ . Otherwise, these hydroxo complexes may form (also locally during

the addition), thus falsifying the result. Al is best titrated at pH values between 4 and 5.

Since the titration is carried out in slightly acidic solution,  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  do not interfere, if present.  $\text{Fe}^{3+}$  is also detected quantitatively and must be determined separately.

### Table of equivalent masses

1 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L} = x \text{ mg equivalent}$ :

Formula	Molar mass in g/mol	Equivalent in mg
Al	26.982	2.698
$\text{Al}(\text{OH})_3$	78.004	7.800
$\text{Al}_2\text{O}_3$	101.961	5.098
$\text{AlCl}_3$	133.341	13.334
$\text{Al}(\text{NO}_3)_3$	212.996	21.300
$\text{Al}(\text{NO}_3)_3 \cdot 9 \text{ H}_2\text{O}$	375.134	37.513
$\text{Al}_2(\text{SO}_4)_3$	342.151	17.108
$\text{Al}_2(\text{SO}_4)_3 \cdot 18 \text{ H}_2\text{O}$	666.426	33.321

### Potentiometric with the Cu-ISE

#### Reagents

- Titrant 1:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Titrant 2:  $c(\text{CuSO}_4) = 0.1 \text{ mol/L}$ ; dissolve 24.97 g  $\text{CuSO}_4 \cdot 5 \text{ H}_2\text{O}$  in approx. 500 mL dist.  $\text{H}_2\text{O}$ , add 0.5 mL  $w/w(\text{H}_2\text{SO}_4) = 96\%$ , mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Acetate buffer,  $\text{pH} = 4.7$ : Dissolve 116 g ammonium acetate in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 86 mL glacial acetic acid. Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .

#### Analysis

If necessary, dilute the sample solution, which should not contain more than 15 mg  $\text{Al}^{3+}$ , with dist.  $\text{H}_2\text{O}$  to approx. 50 mL. Add 5 mL acetate buffer and 10 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$ , allow to react for 1–3 min while stirring, and titrate with  $c(\text{CuSO}_4) = 0.1 \text{ mol/L}$  past the first equivalence point.



### Calculation

$$\beta_{\text{Al}} = \frac{(V_{\text{EDTA}} \times f_1 \times C_{\text{EDTA}} - V_{\text{EP1}} \times f_2 \times C_{\text{CuSO}_4}) \times M_{\text{Al}}}{V_S}$$

$\beta_{\text{Al}}$ : Aluminum content of the sample in g/L

$V_{\text{EDTA}}$ : Added volume of EDTA standard solution in mL

$f_1$ : Titer of the EDTA standard solution

$C_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L


$V_{\text{EP1}}$ : Volume of  $\text{CuSO}_4$  standard solution consumed up to the first equivalence point in mL

$f_2$ : Titer of the  $\text{CuSO}_4$  standard solution

$C_{\text{CuSO}_4}$ : Concentration of the  $\text{CuSO}_4$  standard solution in mol/L

$M_{\text{Al}}$ : Molar mass of aluminum; 26.982 g/mol

$V_S$ : Sample size in mL



**Figure 42** Example curve of a potentiometric aluminum determination

### Further literature on this topic is given in

- *Application Bulletin AB-101 – Complexometric titrations with the Cu-ISE*
- *Titration Application Note AN-T-103 – Fully automated determination of aluminum in aqueous solution*

### Photometric with the Optrode at 610 nm

#### Reagents

- Titrant:  $c(\text{Bi}(\text{NO}_3)_3) = 0.05 \text{ mol/L}$ ; dissolve 24.25 g  $\text{Bi}(\text{NO}_3)_3 \cdot 5 \text{ H}_2\text{O}$  in approx. 500 mL  $\text{HNO}_3$  (2 mol/L) and transfer into a 1 L volumetric flask. Next, make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Auxiliary solution:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Hydrochloric acid:  $c(\text{HCl}) = 6 \text{ mol/L}$ ; add 590 mL w/w(HCl) = 37% to a 1 L volumetric flask that already contains approx. 200 mL dist.  $\text{H}_2\text{O}$ . After the

solution has cooled, make up to the mark with dist.  $\text{H}_2\text{O}$ .

- Acetate buffer: Dissolve 60 g ammonium acetate in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 400 mL glacial acetic acid. Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Indicator solution: Xylenol orange; dissolve 100 mg xylenol orange disodium salt in 100 mL dist.  $\text{H}_2\text{O}$ .

## Analysis

Pipette 10 mL sample solution into a glass beaker and add 70 mL dist. H<sub>2</sub>O. Add 10 mL acetate buffer, adjust to pH 3.5 with c(HCl) = 6 mol/L, and then add 1.25 mL c(Na<sub>2</sub>EDTA) = 0.1 mol/L and 1 mL indica-

tor solution. Titrate the solution with c(Bi(NO<sub>3</sub>)<sub>3</sub>) = 0.05 mol/L past the first equivalence point.

## Calculation

$$\text{Al}_2\text{O}_3 \% = \frac{((V_{\text{EDTA}} \times c_{\text{EDTA}} \times f_{\text{EDTA}}) - (V_{\text{EP1}} \times c_{\text{Bi}} \times f_{\text{Bi}})) \times M_{\text{Al}_2\text{O}_3} \times 10}{m_s \times 2}$$

Al<sub>2</sub>O<sub>3</sub> %: Aluminum content as aluminum oxide in %

V<sub>EDTA</sub>: Added volume of EDTA standard solution in mL

c<sub>EDTA</sub>: Concentration of the EDTA standard solution in mol/L

f<sub>EDTA</sub>: Titer of the EDTA standard solution

V<sub>EP1</sub>: Volume of c(Bi(NO<sub>3</sub>)<sub>3</sub>) = 0.05 mol/L consumed up to the first equivalence point in mL

c<sub>Bi</sub>: Concentration of the Bi(NO<sub>3</sub>)<sub>3</sub> solution in mol/L


f<sub>Bi</sub>: Titer of the Bi(NO<sub>3</sub>)<sub>3</sub> solution; c(Bi(NO<sub>3</sub>)<sub>3</sub>) = 0.05 mol/L

10: Conversion factor

M<sub>Al<sub>2</sub>O<sub>3</sub></sub>: Molar mass of aluminum oxide; 101.96 g/mol

m<sub>s</sub>: Sample size in g

2: Stoichiometric factor



**Figure 43** Example of a photometric titration curve of aluminum at 610 nm

### Further literature on this topic is given in

- Application Bulletin AB-063 – Determination of silicon, calcium, magnesium, iron, and aluminum in digested cement samples by photometric titrations

# Ba – Barium and Sr – Strontium

## General

Barium and strontium can be titrated either directly or indirectly in alkaline solution.  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  are also determined.  $\text{Ba}^{2+}$  can be separated by precipitation as  $\text{BaSO}_4$  (Caution:  $\text{BaSO}_4$  dissolves in excess  $\text{Na}_2\text{EDTA}$ ). Heavy metal ions are masked by adding KCN. Since the titration can only be performed in alkaline solution, in which both barium and strontium form sparingly soluble carbonates, the results of the determination may be too low. This is why the determinations are often carried out using back titration.

The following table contains some solubility products as  $\text{pK}_{\text{sp}}$  values as a guide

$$K_{\text{sp}} = a[\text{A}] \cdot a[\text{B}]$$

$\text{pK}_{\text{sp}}$  is the negative logarithm of  $K_{\text{sp}}$ . As  $\text{pK}_{\text{sp}}$  increases, the solubility of the respective compound decreases.

Compound	$\text{pK}_{\text{sp}}$
$\text{BaCO}_3$	8.2
$\text{BaSO}_4$	10.0
$\text{CaCO}_3$	7.9
$\text{CaSO}_4$	4.3
$\text{MgCO}_3$	3.8
$\text{SrCO}_3$	8.8
$\text{SrSO}_4$	6.6

## Table of equivalent masses

1 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L} = x \text{ mg equivalent}$

Formula	Molar mass in g/mol	Equivalent in mg
Ba	137.327	13.733
$\text{BaCO}_3$	197.336	19.734
$\text{BaCl}_2$	208.233	20.823
$\text{Ba}(\text{NO}_3)_2$	261.337	26.134
BaO	153.326	15.333
$\text{Ba}(\text{OH})_2$	171.342	17.134
$\text{BaSO}_4$	233.390	23.339
Sr	87.620	8.762
$\text{SrCO}_3$	147.629	14.763
$\text{SrCl}_2$	158.526	15.853
$\text{Sr}(\text{NO}_3)_2$	211.630	21.163
SrO	103.619	10.362
$\text{Sr}(\text{OH})_2$	121.635	12.164
$\text{SrSO}_4$	183.683	18.368

## Potentiometric with the Cu-ISE

### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Cu complex solution:  $c[\text{Cu}(\text{NH}_4)_2\text{EDTA}] = 0.1 \text{ mol/L}$  (e.g., Merck No. 105217)
- Ammonia buffer,  $\text{pH} = 10$ : Dissolve 54 g  $\text{NH}_4\text{Cl}$  in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 350 mL  $w/w(\text{NH}_3) = 25\%$ . Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Sodium hydroxide:  $c(\text{NaOH}) = 1 \text{ mol/L}$

### Analysis

Preneutralize very acidic sample solutions to approx.  $\text{pH} 4$  with  $c(\text{NaOH}) = 1 \text{ mol/L}$ . If necessary, dilute the sample solution with dist.  $\text{H}_2\text{O}$  to approx. 50 mL. Add 1 mL Cu complex solution and 5 mL ammonia buffer, wait for 10–20 s while stirring, and then titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.

### Calculation

$$\beta_{\text{Ba/Sr}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Ba/Sr}}}{V_{\text{S}}}$$

$\beta_{\text{Ba/Sr}}$ : Barium or strontium content of the sample in g/L


$V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL

$f$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L

$M_{\text{Ba/Sr}}$ : Molar mass of barium and strontium; 137.327 g/mol and 87.620 g/mol, respectively

$V_{\text{S}}$ : Sample size in mL



**Figure 44** Example of a potentiometric titration curve of barium with the Cu-ISE

### Further literature on this topic is given in

- *Application Bulletin AB-101 – Complexometric titrations with the copper ion-selective electrode*
- *Titration Application Note AN-T-104 – Fully automated determination of barium in aqueous solution*

## Photometric with the Optrode at 574 nm

### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Indicator solution: Phthalein purple; dissolve 0.1 g phthalein purple in approx. 75 mL dist.  $\text{H}_2\text{O}$  and 1 mL  $w/w(\text{NH}_3) = 25\%$ . Mix and make up to 100 mL with dist.  $\text{H}_2\text{O}$ .
- Ammonia:  $w/w(\text{NH}_3) = 25\%$
- Sodium hydroxide:  $c(\text{NaOH}) = 1 \text{ mol/L}$

### Analysis

Preneutralize very acidic samples to pH 4–5 with  $c(\text{NaOH}) = 1 \text{ mol/L}$ . To 50 mL sample solution, which should not contain more than 100 mg  $\text{Ba}^{2+}$  or 80 mg  $\text{Sr}^{2+}$ , add 100 mL dist.  $\text{H}_2\text{O}$ . Add 5 mL  $w/w(\text{NH}_3) = 25\%$  and 0.5 mL indicator solution, and titrate *immediately* with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point (colorless).

### Calculation

$$\beta_{\text{Ba/Sr}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Ba/Sr}}}{V_{\text{S}}}$$

$\beta_{\text{Ba/Sr}}$ : Barium or strontium content of the sample in g/L


$V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL

$f$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L

$M_{\text{Ba/Sr}}$ : Molar mass of barium and strontium; 137.327 g/mol and 87.620 g/mol, respectively

$V_{\text{S}}$ : Sample size in mL



**Figure 45** Photometric titration curve of barium with the Optrode at 574 nm

### Further literature on this topic is given in

- *Titration Application Note AN-T-142 – Determination of barium using automatic photometric titration*

# Bi – Bismuth

46

## General

Similar to  $Al^{3+}$ ,  $Bi^{3+}$  also forms hydroxo complexes in aqueous solutions.  $Bi^{3+}$  reacts sluggishly with the complexing agent but, on the other hand, Bi complexes have a very high complex formation constant. This enables  $Bi^{3+}$  to be titrated without problems even in very acidic solutions (pH 1–2). Complexo-

metric determination of Bi is thus very selective. Interference can only be expected from metal ions whose complex formation constant is  $> 20$  ( $Fe^{3+}$ ,  $Hg^{2+}$ ,  $Sn^{2+}$ ,  $Sn^{4+}$ ). Determination of Bi in some pharmaceuticals (e.g., Bi subgallate, Bi subsalicylate) is carried out after digestion of the sample.

## Table of equivalent masses

1 mL  $c(Na_2EDTA) = 0.1$  mol/L = x mg equivalent

Formula	Molar mass in g/mol	Equivalent in mg
Bi	208.980	20.898
$Bi(OH)_3$	261.013	26.101
$Bi_2O_3$	465.959	23.298
$(BiO)_2CO_3$	509.969	25.498
$BiONO_3 \cdot H_2O$	305.000	30.500
$Bi(NO_3)_3$	394.995	39.500

## Potentiometric with the Cu ISE

### Reagents

- Titrant 1:  $c(Na_2EDTA) = 0.1$  mol/L
- Titrant 2:  $c(CuSO_4) = 0.1$  mol/L; dissolve 24.97 g  $CuSO_4 \cdot 5 H_2O$  in approx. 500 mL dist.  $H_2O$ . Add 0.5 mL  $w/w(H_2SO_4) = 96\%$ , mix, and make up to 1 L with dist.  $H_2O$ .
- Acetate buffer, pH = 4.7: Dissolve 116 g ammonium acetate in approx. 200 mL dist.  $H_2O$  and add 86 mL glacial acetic acid. Mix, and make up to 1 L with dist.  $H_2O$ .
- Ammonia:  $w/w(NH_3) = 25\%$

### Analysis

If necessary, dilute the acidic sample solution, which should not contain more than 160 mg  $Bi^{3+}$ , with dist.  $H_2O$  to approx. 50 mL. Add 5 mL acetate buffer and 10 mL  $c(Na_2EDTA) = 0.1$  mol/L, and, if necessary, adjust to pH 4.7 with  $w/w(NH_3) = 25\%$ . Allow to react for approx. 1 min while stirring, and titrate with  $c(CuSO_4) = 0.1$  mol/L past the first equivalence point.

### Calculation

$$\beta_{\text{Bi}} = \frac{(V_{\text{EDTA}} \times f_1 \times c_{\text{EDTA}} - V_{\text{EP1}} \times f_2 \times c_{\text{CuSO}_4}) \times M_{\text{Bi}}}{V_{\text{S}}}$$

- $\beta_{\text{Bi}}$ : Bismuth content of the sample in g/L
- $V_{\text{EDTA}}$ : Added volume of EDTA standard solution in mL
- $f_1$ : Titer of the EDTA standard solution
- $c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L
- $V_{\text{EP1}}$ : Volume of  $\text{CuSO}_4$  standard solution consumed up to the first equivalence point in mL
- $f_2$ : Titer of the  $\text{CuSO}_4$  standard solution
- $c_{\text{CuSO}_4}$ : Concentration of the  $\text{CuSO}_4$  standard solution in mol/L
- $M_{\text{Bi}}$ : Molar mass of bismuth; 208.980 g/mol
- $V_{\text{S}}$ : Sample size in mL


### Photometric with the Optrode at 520 nm

#### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Nitric acid:  $w/w(\text{HNO}_3) = 65\%$
- Indicator solution: Xylenol orange; dissolve 100 mg xylenol orange disodium salt in 100 mL dist.  $\text{H}_2\text{O}$ .

#### Analysis

Dissolve the sample in a small amount of  $w/w(\text{HNO}_3) = 65\%$ , and dilute to approx. 100 mL with dist.  $\text{H}_2\text{O}$ . Add 0.5 mL indicator solution, and titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.



**Figure 46** Potentiometric titration curve of bismuth with the Cu-ISE


#### Further literature on this topic is given in

- *Application Bulletin AB-101 – Complexometric titrations with the Cu-ISE*
- *Titration Application Note AN-T-105 – Fully automated determination of bismuth(III) in aqueous solution*

### Calculation

$$\beta_{\text{Bi}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Bi}}}{V_{\text{S}}}$$

- $\beta_{\text{Bi}}$ : Bismuth content of the sample in g/L
- $V_{\text{EP1}}$ : Volume of titrant consumed up to the first equivalence point in mL
- $f$ : Titer of the EDTA standard solution
- $c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L
- $M_{\text{Bi}}$ : Molar mass of bismuth; 208.980 g/mol
- $V_{\text{S}}$ : Sample size in mL



**Figure 47** Example of a photometric titration curve of bismuth at 520 nm

### Further literature on this topic is given in

- *Titration Application Note AN-T-088 – Photometric EDTA titration of bismuth nitrate according to Ph. Eur. and USP*

## Ca – Calcium

### General

$\text{Ca}^{2+}$  is one of the first metals that was determined using complexometric titration. Since then, hundreds, if not thousands, of publications on this topic have appeared. In this chapter, we will only discuss the determination of  $\text{Ca}^{2+}$ . If the sample contains both  $\text{Ca}^{2+}$  with  $\text{Mg}^{2+}$ , please refer to the section on "Water hardness".

Since  $\text{Ca}^{2+}$  is always titrated in alkaline solution, there is always a risk that it may absorb  $\text{CO}_2$  from the air and precipitate

as  $\text{CaCO}_3$ . If the titration is carried out too quickly, this could lead to incorrect results that are too low. ( $\text{CaCO}_3$  reacts slowly with the titrant, which is acting as a complexing agent, and the resulting turbidity interferes, particularly in the case of photometric titration).

If  $\text{Ca}^{2+}$  is to be determined in the presence of  $\text{Mg}^{2+}$ , the EDTA can be replaced by EGTA, which is a more selective titrant (or by the Cu-EGTA complex if a Cu ISE is used).



**Table of equivalent masses**

1 mL  $c(\text{CA}) = 0.1 \text{ mol/L} = x \text{ mg equivalent (CA = complexing agent)}$

Formula	Molar mass in g/mol	Equivalent in mg
Ca	40.078	4.008
CaC <sub>2</sub> O <sub>4</sub>	128.097	12.810
CaCO <sub>3</sub>	100.087	10.009
CaCl <sub>2</sub>	110.984	11.098
Ca(NO <sub>3</sub> ) <sub>2</sub>	164.088	16.409
CaO	56.077	5.608
Ca(OH) <sub>2</sub>	74.093	7.409
CaSO <sub>3</sub>	120.141	12.014
CaSO <sub>4</sub>	136.141	13.614

**Potentiometric with the Ca-ISE****Reagents**

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Buffer solution: Dissolve 54 g  $\text{NH}_4\text{Cl}$  in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 350 mL  $w/w(\text{NH}_3) = 25\%$ . Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Sodium hydroxide:  $c(\text{NaOH}) = 1 \text{ mol/L}$

**Analysis**

Preneutralize acidic sample solutions with  $c(\text{NaOH}) = 1 \text{ mol/L}$  to pH 5–7 and dilute, if necessary, with dist.  $\text{H}_2\text{O}$  to approx. 50 mL. Add 5 mL buffer solution and titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.

**Calculation**

$$\beta_{\text{Ca}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Ca}}}{V_{\text{S}}}$$

$\beta_{\text{Ca}}$ : Calcium content of the sample in g/L


$V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL

$f$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L

$M_{\text{Ca}}$ : Molar mass of calcium; 40.078 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 48** Example curve of a potentiometric calcium determination with the Ca-ISE

### Further literature on this topic is given in

- *Application Bulletin AB-125 – Simultaneous determination of calcium, magnesium, and alkalinity by complexometric titration with potentiometric or photometric indication in water and beverage samples*

### Potentiometric with the Cu ISE Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Buffer solution: Dissolve 54 g  $\text{NH}_4\text{Cl}$  in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 350 mL  $w/w(\text{NH}_3) = 25\%$ . Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Cu complex solution:  $c[\text{Cu}(\text{NH}_4)_2\text{EDTA}] = 0.1 \text{ mol/L}$  (e.g., Merck No. 105217)
- Sodium hydroxide:  $c(\text{NaOH}) = 1 \text{ mol/L}$

### Analysis

Preneutralize acidic sample solutions with  $c(\text{NaOH}) = 1 \text{ mol/L}$  to pH 5–7 and dilute, if necessary, with dist.  $\text{H}_2\text{O}$  to approx. 50 mL. Add 5 mL buffer solution and 0.5 mL Cu complex solution, wait 20 s while stirring, and titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.


### Further literature on this topic is given in

- *Application Bulletin AB-101 – Complexometric titrations with the Cu-ISE*
- *Titration Application Note AN-T-106 – Fully automated determination of calcium in milk*

### Calculation

$$\beta_{\text{Ca}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Ca}}}{V_{\text{S}}}$$

- $\beta_{\text{Ca}}$ : Calcium content of the sample in g/L
- $V_{\text{EP1}}$ : Volume in mL of EDTA standard solution consumed up to the first equivalence point
- f: Titer of the EDTA standard solution
- $c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L
- $M_{\text{Ca}}$ : Molar mass of calcium; 40.078 g/mol
- $V_{\text{S}}$ : Sample size in mL



**Figure 49** Potentiometric titration curve of the determination of calcium in milk with the Cu-ISE

## Photometric with the Optrode at 610 nm

### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Sodium hydroxide:  $c(\text{NaOH}) = 1 \text{ mol/L}$
- Indicator solution: HHSNN; dissolve 20 mg HHSNN in 50 mL ethanol.


### Analysis

Preneutralize acidic sample solutions with  $c(\text{NaOH}) = 1 \text{ mol/L}$  to pH 5–7 and dilute, if necessary, with dist.  $\text{H}_2\text{O}$  to approx. 50 mL. Add 4 mL  $c(\text{NaOH}) = 1 \text{ mol/L}$  and 1.5 mL indicator solution, and titrate *immediately* with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.

### Calculation

$$\beta_{\text{Ca}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Ca}}}{V_{\text{S}}}$$

- $\beta_{\text{Ca}}$ : Calcium content of the sample in g/L
- $V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL
- f: Titer of the EDTA standard solution
- $c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L
- $M_{\text{Ca}}$ : Molar mass of calcium; 40.078 g/mol
- $V_{\text{S}}$ : Sample size in mL



**Figure 50** Example of a photometric titration curve of the calcium determination with the Optrode at 610 nm

### Further literature on this topic is given in

- *Application Bulletin AB-125 – Simultaneous determination of calcium, magnesium, and alkalinity by complexometric titration with potentiometric or photometric indication in water and beverage samples*

# Cd – Cadmium

52

## General

The complexation behavior of  $\text{Cd}^{2+}$  is very similar to that of zinc. It can be titrated in alkaline solution without

problems. However,  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  interfere if present because they are also determined.

## Table of equivalent masses

1 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L} = x \text{ mg equivalent}$

Formula	Molar mass in g/mol	Equivalent in mg
Cd	112.411	11.241
$\text{CdCO}_3$	172.420	17.242
$\text{CdCl}_2$	183.317	18.332
$\text{Cd}(\text{NO}_3)_2$	236.421	23.642
$\text{Cd}(\text{NO}_3)_2 \cdot 2\text{H}_2\text{O}$	272.451	27.245
CdO	128.410	12.841
$\text{Cd}(\text{OH})_2$	146.426	14.643
$\text{CdSO}_4$	208.474	20.847

## Potentiometric with the Cu-ISE

### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Ammonia buffer: Dissolve 54 g  $\text{NH}_4\text{Cl}$  in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 350 mL  $w/w(\text{NH}_3) = 25\%$ . Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Cu complex solution:  $c[\text{Cu}(\text{NH}_4)_2\text{EDTA}] = 0.1 \text{ mol/L}$  (e.g., Merck No. 105217)
- Sodium hydroxide:  $c(\text{NaOH}) = 1 \text{ mol/L}$

### Analysis

Preneutralize acidic sample solutions with  $c(\text{NaOH}) = 1 \text{ mol/L}$  to pH 4–5 and dilute, if necessary, with dist.  $\text{H}_2\text{O}$  to approx. 50 mL. Add 5 mL ammonia buffer and 0.5 mL Cu complex solution, wait 20 s while stirring, and titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.

### Calculation

$$\beta_{\text{Cd}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Cd}}}{V_{\text{S}}}$$

$\beta_{\text{Cd}}$ : Cadmium content of the sample in g/L


$V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL

$f$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L

$M_{\text{Cd}}$ : Molar mass of cadmium; 112.411 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 51** Potentiometric titration curve of cadmium with the Cu-ISE

## Photometric with the Optrode at 610 nm

### Reagents

- Titrant 1:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Titrant 2:  $c(\text{ZnSO}_4) = 0.1 \text{ mol/L}$ ;  
dissolve 28.8 g  $\text{ZnSO}_4 \cdot 7 \text{ H}_2\text{O}$  in approx. 500 mL dist.  $\text{H}_2\text{O}$ , add 0.5 mL  $w/w(\text{H}_2\text{SO}_4) = 96\%$ , mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Ammonia buffer,  $\text{pH} = 10$ : Dissolve 54 g  $\text{NH}_4\text{Cl}$  in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 350 mL  $w/w(\text{NH}_3) = 25\%$ . Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Sodium hydroxide:  $c(\text{NaOH}) = 1 \text{ mol/L}$
- Indicator solution: Eriochrome Black T; dissolve 100 mg Eriochrome Black T in 100 mL ethanol.

### Analysis

Preneutralize acidic sample solutions with  $c(\text{NaOH}) = 1 \text{ mol/L}$  to approx.  $\text{pH} 5$  and dilute, if necessary, with dist.  $\text{H}_2\text{O}$  to approx. 50 mL. Add 5 mL ammonia buffer, 10 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$ , and 0.25 mL indicator solution. Then titrate with  $c(\text{ZnSO}_4) = 0.1 \text{ mol/L}$  past the first equivalence point.

### Further literature on this topic is given in

- *Application Bulletin AB-101 – Complexometric titrations with the Cu-ISE*
- *Titration Application Note AN-T-125 – Automatic determination of cadmium in aqueous solution with the Cu-ISE*

### Calculation

$$\beta_{\text{Cd}} = \frac{(V_{\text{EDTA}} \times f_1 \times c_{\text{EDTA}} - V_{\text{EP1}} \times f_2 \times c_{\text{ZnSO}_4}) \times M_{\text{Cd}}}{V_{\text{S}}}$$

$\beta_{\text{Cd}}$ : Cadmium content of the sample in g/L

$V_{\text{EDTA}}$ : Added volume of EDTA standard solution in mL

$f_1$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L


$V_{\text{EP1}}$ : Volume of  $\text{ZnSO}_4$  standard solution consumed up to the first equivalence point in mL

$f_2$ : Titer of the  $\text{ZnSO}_4$  standard solution

$c_{\text{ZnSO}_4}$ : Concentration of the  $\text{ZnSO}_4$  standard solution in mol/L

$M_{\text{Cd}}$ : Molar mass of cadmium; 112.411 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 52** Photometric titration curve of the cadmium determination with the Optrode at 610 nm

### Further literature on this topic is given in

- *Titration Application Note AN-T-153 – Determination of cadmium using automatic photometric titration*

## Co – Cobalt

### General

The complexation behavior of cobalt is similar to that of nickel, and it has almost the same complex formation constants. In alkaline solution and in the presence of excess EDTA as well as added  $\text{H}_2\text{O}_2$ ,  $\text{Co}^{2+}$  forms extremely robust and strongly colored  $\text{Co}^{3+}$ -EDTA complexes that are also stable in acidic media. This situation was

used by various authors for the complexometric determination of cobalt with enhanced selectivity. These methods, however, did not gain acceptance because the quantity of added  $\text{H}_2\text{O}_2$  is critical and, in the case of photometric titration, the intrinsic color of the  $\text{Co}^{3+}$  complex interferes with the color change of the indicator.

### Table of equivalent masses

1 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L} = x \text{ mg equivalent}$

Formula	Molar mass in g/mol	Equivalent in mg
Co	58.933	5.893
$\text{Co}(\text{CH}_3\text{COO})_2 \cdot 4 \text{H}_2\text{O}$	249.083	24.908
$\text{CoCO}_3$	118.942	11.894
$\text{CoCl}_2$	129.839	12.984
$\text{CoCl}_2 \cdot 6\text{H}_2\text{O}$	237.930	23.793
$\text{Co}(\text{NO}_3)_2$	182.943	18.294
$\text{Co}(\text{NO}_3)_2 \cdot 6 \text{H}_2\text{O}$	291.034	29.103
CoO	74.933	7.493
$\text{Co}_2\text{O}_3$	165.865	8.293
$\text{CoSO}_4$	154.996	15.500
$\text{CoSO}_4 \cdot 7\text{H}_2\text{O}$	281.103	28.110

## Potentiometric with the Cu-ISE

### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Cu complex solution:  $c[\text{Cu}(\text{NH}_4)_2\text{EDTA}] = 0.1 \text{ mol/L}$  (e.g., Merck No. 105217)
- Ammonia buffer: Dissolve 54 g  $\text{NH}_4\text{Cl}$  in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 350 mL  $w/w(\text{NH}_3) = 25\%$ . Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Sodium hydroxide:  $c(\text{NaOH}) = 1 \text{ mol/L}$

### Analysis

If necessary, dilute the sample solution with dist.  $\text{H}_2\text{O}$  to approx. 50 mL and pre-neutralize, if necessary, to pH 4–5 with  $c(\text{NaOH}) = 1 \text{ mol/L}$ . Add 5 mL ammonia buffer and 0.5 mL Cu complex solution, wait for 10–20 s while stirring, and titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.

### Calculation

$$\beta_{\text{Co}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Co}}}{V_{\text{S}}}$$

$\beta_{\text{Co}}$ : Cobalt content of the sample in g/L


$V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL

$f$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L

$M_{\text{Co}}$ : Molar mass of cobalt; 58.933 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 53** Example of a potentiometric titration curve of cobalt with the Cu-ISE

### Further literature on this topic is given in

- *Application Bulletin AB-101 – Complexometric titrations with the Cu-ISE*
- *Titration Application Note AN-T-126 – Automatic determination of cobalt in aqueous solution with the Cu-ISE*

## Photometric with the Optrode at 574 nm

### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Indicator solution: murexide; dissolve 0.2 g murexide (trituated 1:100 with NaCl) in 50 mL dist.  $\text{H}_2\text{O}$ .
- Ammonia buffer: Dissolve 54 g  $\text{NH}_4\text{Cl}$  in approx. 500 mL dist.  $\text{H}_2\text{O}$  and add 350 mL  $w/w(\text{NH}_3) = 25\%$ . Adjust to pH 9 with  $c(\text{HCl}) = 6 \text{ mol/L}$ , mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .

### Analysis

If necessary, dilute the acidic sample solution, which should not contain more than 40 mg  $\text{Co}^{2+}$ , with dist.  $\text{H}_2\text{O}$  to 100 mL and degas. Add 10 mL ammonia buffer solution and 5 mL indicator solution, and then titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.

### Calculation

$$\beta_{\text{Co}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Co}}}{V_{\text{S}}}$$

$\beta_{\text{Co}}$ : Cobalt content of the sample in g/L


$V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL

$f$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L

$M_{\text{Co}}$ : Molar mass of cobalt; 58.933 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 54** Photometric titration curve of cobalt with the Optrode at 574 nm

### Further literature on this topic is given in

- *Titration Application Note AN-T-143 – Determination of cobalt using automatic photometric titration*



# Cu – Copper

57

## General

Owing to its relatively high complex formation constant, copper can be titrated in both alkaline and slightly acidic solutions. Only the  $\text{Cu}^{2+}$  ion can be determined. In case of photometric titration in an ammoniacal solution, the ammonia excess should not be too large and the titrations should always be carried out in dilute solutions. This is because the  $\text{Cu}^{2+}$ -amine complex is relatively stable and strongly colored and thus competes with the color indicator (e.g., murexide).

Keep adding  $\text{NH}_3$  to the sample solution until  $\text{Cu}(\text{OH})_2$  precipitates. Then add further  $\text{NH}_3$  until everything has re-dissolved. Dilute with dist.  $\text{H}_2\text{O}$  to at least 1:1. Add murexide, and titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$ .

Better results are obtained by photometric titration at  $\text{pH} = 5.0$  with PAN as indicator. To avoid having to titrate at high temperatures, dilute the sample solution with ethanol in an approx. 1:1 ratio.

## Table of equivalent masses

1 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L} = x \text{ mg equivalent}$

Formula	Molar mass in g/mol	Equivalent in mg
Cu	63.546	6.355
$\text{Cu}(\text{CH}_3\text{COO})_2 \cdot \text{H}_2\text{O}$	199.650	19.965
CuCN	89.564	8.956
$\text{Cu}(\text{CN})_2$	115.581	11.558
$\text{CuCl}_2$	134.452	13.445
$\text{CuCl}_2 \cdot 2 \text{H}_2\text{O}$	170.842	17.084
$\text{Cu}(\text{NO}_3)_2$	187.556	18.756
$\text{Cu}(\text{NO}_3)_2 \cdot 3 \text{H}_2\text{O}$	241.601	24.160
CuO	79.544	7.954
$\text{Cu}(\text{OH})_2$	97.561	9.756
$\text{CuSO}_4$	159.600	15.960
$\text{CuSO}_4 \cdot 5 \text{H}_2\text{O}$	249.680	24.968

## Potentiometric with the Cu-ISE

### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Ammonia buffer,  $\text{pH} = 10$ : Dissolve 54 g  $\text{NH}_4\text{Cl}$  in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 350 mL  $w/w(\text{NH}_3) = 25\%$ . Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Sodium hydroxide:  $c(\text{NaOH}) = 1 \text{ mol/L}$

### Analysis

If necessary, dilute acidic sample solutions with dist.  $\text{H}_2\text{O}$  to approx. 50 mL and preneutralize to  $\text{pH} 4\text{--}5$  with  $c(\text{NaOH}) = 1 \text{ mol/L}$ . Add 5 mL ammonia buffer, and titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.

### Calculation

$$\beta_{\text{Cu}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Cu}}}{V_{\text{S}}}$$

$\beta_{\text{Cu}}$ : Copper content of the sample in g/L


$V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL

$f$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L

$M_{\text{Cu}}$ : Molar mass of copper; 63.546 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 55** Potentiometric titration curve of copper with the Cu-ISE

### Further literature on this topic is given in

- *Application Bulletin AB-101 – Complexometric titrations with the Cu-ISE*
- *Titration Application Note AN-T-127 – Automatic determination of copper in aqueous solution with the Cu-ISE*

## Photometric with the Optrode at 520 nm

### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Acetate buffer,  $\text{pH} = 5.0$ : dissolve 116 g ammonium acetate in approx. 500 mL dist.  $\text{H}_2\text{O}$ . Adjust to  $\text{pH} 5.0$  with glacial acetic acid, mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Sodium hydroxide:  $c(\text{NaOH}) = 1 \text{ mol/L}$
- Ethanol:  $w/w(\text{ethanol}) = 96\text{--}98\%$
- Indicator solution: PAN; dissolve 100 mg PAN in 100 mL ethanol.

### Analysis

If necessary, dilute acidic sample solutions with dist.  $\text{H}_2\text{O}$  to approx. 50 mL and then preneutralize to approx.  $\text{pH} 4\text{--}5$  with  $c(\text{NaOH}) = 1 \text{ mol/L}$ . Add 50 mL ethanol and degas. Add 5 mL acetate buffer and 0.5 mL indicator solution, and titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.

### Calculation

$$\beta_{\text{Cu}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Cu}}}{V_{\text{S}}}$$

$\beta_{\text{Cu}}$ : Copper content of the sample in g/L


$V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL

$f$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L

$M_{\text{Cu}}$ : Molar mass of copper; 63.546 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 56** Photometric titration curve of copper with the Optrode at 520 nm

### Further literature on this topic is given in

- *Titration Application Note AN-T-124 – Photometric copper determination in aqueous solution*

# Fe – Iron

60

## General

Iron forms complexes in both its divalent and trivalent states. The complex formation constants, however, are very different (approx. 12 orders of magnitude).  $\text{Fe}^{3+}$  forms very strong complexes that are also stable in acidic solution. This means that Fe can be titrated in the presence of other metal ions almost without interference. In the presence of, for example, EDTA,  $\text{Fe}^{2+}$  becomes a very strong reducing agent (standard potential  $-0.1\text{ V}$ ), which can even reduce  $\text{Ag}^+$  to its metal.

As a general rule, titration of  $\text{Fe}^{2+}$  ions is not advisable. Owing to the relatively low complex formation constant, it must be titrated in alkaline solution under inert gas to prevent oxidation to  $\text{Fe}^{3+}$ . Conversely, interfering  $\text{Fe}^{3+}$  can be masked by reducing it, for example, with ascorbic acid. Thus, only the  $\text{Fe}^{3+}$  ion is titrated in practice.  $\text{Fe}^{2+}$  can be oxidized using  $\text{H}_2\text{O}_2$  or peroxydisulfate.

## Table of equivalent masses

1 mL  $c(\text{Na}_2\text{EDTA}) = 0.1\text{ mol/L} = x\text{ mg equivalent}$

Formula	Molar mass in g/mol	Equivalent in mg
Fe	55.845	5.585
$\text{FeCO}_3$	115.856	11.586
$\text{FeCl}_2$	126.753	12.675
$\text{FeCl}_2 \cdot 4\text{H}_2\text{O}$	198.813	19.881
$\text{FeCl}_3$	162.206	16.221
$\text{FeCl}_3 \cdot 6\text{H}_2\text{O}$	270.297	27.030
$\text{Fe}(\text{NO}_3)_3$	241.862	24.186
$\text{Fe}(\text{NO}_3)_3 \cdot 6\text{H}_2\text{O}$	349.953	34.995
FeO	71.846	7.185
$\text{Fe}_2\text{O}_3$	159.692	7.985
$\text{FeSO}_4$	151.910	15.191
$\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$	278.010	27.801
$\text{Fe}_2(\text{SO}_4)_3$	399.870	19.994
$\text{Fe}_2(\text{SO}_4)_3 \cdot 9\text{H}_2\text{O}$	562.010	28.101
$\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}$	392.130	39.213
$\text{Fe}(\text{NH}_4)(\text{SO}_4)_2 \cdot 12\text{H}_2\text{O}$	482.180	48.218

## Potentiometric with the Cu-ISE

### Reagents

- Titrant 1:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Titrant 2:  $c(\text{CuSO}_4) = 0.1 \text{ mol/L}$ ;  
dissolve 24.97 g  $\text{CuSO}_4 \cdot 5 \text{ H}_2\text{O}$  in  
approx. 500 mL dist.  $\text{H}_2\text{O}$ , add  
0.5 mL w/w( $\text{H}_2\text{SO}_4$ ) = 96%, mix, and  
make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Acetate buffer, pH = 4.7: Dissolve  
116 g ammonium acetate in approx.  
200 mL dist.  $\text{H}_2\text{O}$  and add 86 mL gla-  
cial acetic acid. Mix and make up to  
1 L with dist.  $\text{H}_2\text{O}$ .
- Sodium hydroxide:  $c(\text{NaOH}) =$   
1 mol/L


### Analysis

If necessary, dilute the acidic sample solution, which should not contain more than 50 mg  $\text{Fe}^{3+}$ , with dist.  $\text{H}_2\text{O}$  to approx. 50 mL, add 5 mL acetate buffer, and adjust to pH 4.7 with  $c(\text{NaOH}) = 1 \text{ mol/L}$ . Add 10 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$ , allow to react for 1–2 min while stirring, and titrate with  $c(\text{CuSO}_4) = 0.1 \text{ mol/L}$  past the first equivalence point.

### Calculation

$$\beta_{\text{Fe}} = \frac{(V_{\text{EDTA}} \times f_1 \times c_{\text{EDTA}} - V_{\text{EP1}} \times f_2 \times c_{\text{CuSO}_4}) \times M_{\text{Fe}}}{V_{\text{S}}}$$

- $\beta_{\text{Fe}}$ : Iron content of the sample in g/L  
 $V_{\text{EDTA}}$ : Added volume of EDTA standard solution in mL  
 $f_1$ : Titer of the EDTA standard solution  
 $c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L  
 $V_{\text{EP1}}$ : Volume of  $\text{CuSO}_4$  standard solution consumed up to the first equivalence point in mL  
 $f_2$ : Titer of the  $\text{CuSO}_4$  standard solution  
 $c_{\text{CuSO}_4}$ : Concentration of the  $\text{CuSO}_4$  standard solution in mol/L  
 $M_{\text{Fe}}$ : Molar mass of iron; 55.845 g/mol  
 $V_{\text{S}}$ : Sample size in mL



**Figure 57** Example of a potentiometric titration curve of iron with the Cu-ISE

### Further literature on this topic is given in

- *Application Bulletin AB-101 – Complexometric titrations with the Cu-ISE*
- *Titration Application Note AN-T-107 – Fully automated determination of total iron in cement*

## Photometric with the Optrode at 610 nm

### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Ammonia:  $w/w(\text{NH}_3) = 25\%$
- Indicator solution: sulfosalicylic acid; dissolve 4 g sulfosalicylic acid in 100 mL dist.  $\text{H}_2\text{O}$ .
- Glycine

### Analysis

If necessary, dilute the acidic sample solution with dist.  $\text{H}_2\text{O}$  to approx. 100 mL. Add 0.5 g glycine and 1 mL indicator solution, and adjust to pH 1.5–2 with  $w/w(\text{NH}_3) = 25\%$ . Heat the mixture to 50 °C to dissolve the glycine. Allow to cool, then titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.

### Calculation

$$\beta_{\text{Fe}} = \frac{V_{\text{BP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Fe}}}{V_{\text{S}}}$$

$\beta_{\text{Fe}}$ : Iron content of the sample in g/L


$V_{\text{BP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL

$f$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L

$M_{\text{A}}$ : Molar mass of iron; 55.847 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 58** Photometric titration curve of iron with the Optrode at 610 nm

### Further literature on this topic is given in

- *Application Bulletin AB-063 – Determination of silicon, calcium, magnesium, iron, and aluminum in digested cement samples by photometric titrations.*
- *Titration Application Note AN-T-080 – Determination of iron in cement by photometric titration*

# Ga – Gallium and In – Indium

## General

Gallium and indium form complexes with quite high formation constants (e.g., InEDTA;  $\log K_f = 24.9$ ). However, they also tend to form hydroxo complexes that cannot be completely titrated, if at all. Therefore acetate ions are added as auxiliary complexing agents, and the metals are determined by back titration in slightly acidic solution. This avoids possible interference by  $\text{Ca}^{2+}$ - and/or  $\text{Mg}^{2+}$  ions.

Photometric back titration with a Mn standard solution with Eriochrome Black T indicator would also be possible in alkaline solution. However, the usual ammonia buffer pH 10 cannot be used for pH adjustment in this case, because this would lead to precipitation of Ga and In as  $\text{Ga}(\text{OH})_3$  and  $\text{In}(\text{OH})_3$ , respectively. To avoid this, tartrate is added as auxiliary complexing agent before adjusting the pH to 8.5–9 using ammonia. If present,  $\text{Zn}^{2+}$ ,  $\text{Cd}^{2+}$ ,  $\text{Hg}^{2+}$ ,  $\text{Cu}^{2+}$ ,  $\text{Ni}^{2+}$ , and  $\text{Co}^{2+}$  can be masked with KCN.

## Table of equivalent masses

1 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L} = x \text{ mg equivalent}$

Formula	Molar mass in g/mol	Equivalent in mg
Ga	69.723	6.972
$\text{Ga}(\text{NO}_3)_3 \cdot 8\text{H}_2\text{O}$	399.860	39.986
$\text{Ga}_2\text{O}_3$	187.444	9.372
$\text{Ga}_2(\text{SO}_4)_3$	427.634	21.382
$\text{Ga}_2(\text{SO}_4)_3 \cdot 18\text{H}_2\text{O}$	751.909	37.595
In	114.818	11.482
$\text{InCl}_3$	221.177	22.118
$\text{In}_2\text{O}_3$	277.634	13.882
$\text{In}_2(\text{SO}_4)_3$	517.824	25.891
$\text{In}_2(\text{SO}_4)_3 \cdot 9\text{H}_2\text{O}$	679.961	33.998

## Potentiometric with the Cu-ISE

### Reagents

- Titrant 1:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Titrant 2:  $c(\text{CuSO}_4) = 0.1 \text{ mol/L}$ ; dissolve 24.97 g  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$  in approx. 500 mL dist.  $\text{H}_2\text{O}$ , add 0.5 mL  $w/w(\text{H}_2\text{SO}_4) = 96\%$ , mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Acetate buffer, pH = 4.7: Dissolve 116 g ammonium acetate in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 86 mL glacial acetic acid. Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Ammonia:  $w/w(\text{NH}_3) = 25\%$


## Analysis

If necessary, dilute the acidic sample solution, which should not contain more than 60 mg Ga<sup>3+</sup> and/or 100 mg In<sup>3+</sup>, with dist. H<sub>2</sub>O to approx. 50 mL. Add 5 mL acetate buffer and 10 mL c(Na<sub>2</sub>EDTA) = 0.1 mol/L, and, if necessary, adjust to pH 4.7 with w/w(NH<sub>3</sub>) = 25%. Allow to react for approx. 30 s while stirring, and then titrate with c(CuSO<sub>4</sub>) = 0.1 mol/L past the first equivalence point.

## Calculation

$$\beta_{\text{Ga/In}} = \frac{(V_{\text{EDTA}} \times f_1 \times c_{\text{EDTA}} - V_{\text{EP1}} \times f_2 \times c_{\text{CuSO}_4}) \times M_{\text{Ga/In}}}{V_S}$$

- $\beta_{\text{Ga/In}}$ : Gallium or indium content of the sample in g/L
- $V_{\text{EDTA}}$ : Added volume of EDTA standard solution in mL
- $f_1$ : Titer of the EDTA standard solution
- $c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L
- $V_{\text{EP1}}$ : Volume of CuSO<sub>4</sub> standard solution consumed up to the first equivalence point in mL
- $f_2$ : Titer of the CuSO<sub>4</sub> standard solution
- $c_{\text{CuSO}_4}$ : Concentration of the CuSO<sub>4</sub> standard solution in mol/L
- $M_{\text{Ga/In}}$ : Molar mass of gallium and indium; 69.723 g/mol and 114.818 g/mol, respectively
- $V_S$ : Sample size in mL



**Figure 59** Potentiometric titration curve of indium with the Cu-ISE

## Further literature on this topic is given in

- *Application Bulletin AB-101 – Complexometric titrations with the Cu-ISE*
- *Titration Application Note AN-T-121 – Automatic indium determination in aqueous solution using the ion-selective copper electrode (Cu-ISE)*



## Photometric with the Optrode at 610 nm

### Reagents

- Titrant 1:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Titrant 2:  $c(\text{ZnSO}_4) = 0.1 \text{ mol/L}$ ; dissolve 28.8 g  $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$  in approx. 500 mL dist.  $\text{H}_2\text{O}$ , add 0.5 mL w/w( $\text{H}_2\text{SO}_4$ ) = 96%, mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Acetate buffer, pH = 4.7: Dissolve 116 g ammonium acetate in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 86 mL glacial acetic acid. Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Indicator solution: xylenol orange; dissolve 25 mg xylenol orange disodium salt in 100 mL dist.  $\text{H}_2\text{O}$ .

### Analysis

If necessary, dilute the acidic sample solution, which should not contain more than 60 mg  $\text{Ga}^{3+}$  and/or 100 mg  $\text{In}^{3+}$ , with dist.  $\text{H}_2\text{O}$  to approx. 100 mL and degas. Add 5 mL acetate buffer, 10 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$ , and 2 mL indicator solution and titrate with  $c(\text{ZnSO}_4) = 0.1 \text{ mol/L}$  past the first equivalence point.

### Calculation

$$\beta_{\text{Ga/In}} = \frac{(V_{\text{EDTA}} \times f_1 \times c_{\text{EDTA}} - V_{\text{EP1}} \times f_2 \times c_{\text{ZnSO}_4}) \times M_{\text{Ga/In}}}{V_S}$$

$\beta_{\text{Ga/In}}$ : Gallium or indium content of the sample in g/L

$V_{\text{EDTA}}$ : Added volume of EDTA standard solution in mL

$f_1$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L


$V_{\text{EP1}}$ : Volume of  $\text{ZnSO}_4$  standard solution consumed up to the first equivalence point in mL

$f_2$ : Titer of the  $\text{ZnSO}_4$  standard solution

$c_{\text{ZnSO}_4}$ : Concentration of the  $\text{ZnSO}_4$  standard solution in mol/L

$M_{\text{Ga/In}}$ : Molar mass of gallium and indium; 69.723 g/mol and 114.818 g/mol, respectively

$V_S$ : Sample size in mL



**Figure 60** Example of a photometric titration curve of gallium with the Optrode at 610 nm

### Further literature on this topic is given in

- *Titration Application Note AN-T-158 – Determination of gallium using automatic photometric titration*

# Hg – Mercury

66

## General

Mercury forms relatively strong complexes and is thus suitable for titration. However, it must be present as  $\text{Hg}^{2+}$  because  $\text{Hg}^+$  disproportionates into  $\text{Hg}^{2+}$  and metallic Hg in the presence of, e.g., EDTA.  $\text{Hg}^{2+}$  can be masked relatively easily. It can be precipitated as  $\text{HgI}_2$  by adding potassium iodide solution, which liberates the corresponding amount of complexing agent (e.g., EDTA). This allows consecutive titrations, for example, in the presence of  $\text{Cu}^{2+}$ .

In some pharmaceutical formulations,  $\text{Hg}^{2+}$  can be titrated directly without having to separate it from the organic matrix beforehand.

The Cu-ISE cannot be used for potentiometric indication because  $\text{Hg}^{2+}$  poisons the crystal membrane, making the electrode inoperable.

Hg complexes are serious environmental pollutants! The solutions used during the titration should thus be collected and detoxified. We recommend the following procedure for this:

The collected solutions are made alkaline, for example, by adding NaOH. Excess sodium sulfide solution is then added to the stirred mixture. The precipitated (very sparingly soluble)  $\text{Hg}_2\text{S}$  can be collected and disposed of separately.

## Table of equivalent masses

1 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L} = x \text{ mg equivalent}$

Formula	Molar mass in g/mol	Equivalent in mg
Hg	200.590	20.059
$\text{HgCl}_2$	271.496	27.150
$\text{Hg}(\text{NO}_3)_2$	324.600	32.460
$\text{Hg}(\text{NO}_3)_2 \cdot \text{H}_2\text{O}$	342.615	34.262
HgO	216.589	21.659
$\text{HgSO}_4$	296.653	29.665

## Photometric with the Optrode at 502 nm

### Reagents

- Titrant 1:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Titrant 2:  $c(\text{ZnSO}_4) = 0.1 \text{ mol/L}$ ;  
dissolve 28.8 g  $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$  in dist.  $\text{H}_2\text{O}$ , mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Ammonia buffer, pH = 10: dissolve 54 g  $\text{NH}_4\text{Cl}$  in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 350 mL w/w( $\text{NH}_3$ ) = 25%. Mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Sodium hydroxide:  $c(\text{NaOH}) = 1 \text{ mol/L}$
- Indicator solution: Eriochrome Black T; dissolve 100 mg Eriochrome Black T in 100 mL ethanol.

### Analysis

If necessary, dilute the acidic sample solution, which should not contain more than 75 mg  $\text{Hg}^{2+}$ , with dist.  $\text{H}_2\text{O}$  to approx. 50 mL. Add 5.00 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  and neutralize the solution with  $c(\text{NaOH}) = 1 \text{ mol/L}$  to approx. pH 5–7. Add 5 mL ammonia buffer and 1 mL indicator solution, and titrate with  $c(\text{ZnSO}_4) = 0.1 \text{ mol/L}$  past the first equivalence point.

### Calculation

$$\beta_{\text{Hg}} = \frac{(V_{\text{EDTA}} \times f_1 \times c_{\text{EDTA}} - V_{\text{EP1}} \times f_2 \times c_{\text{ZnSO}_4}) \times M_{\text{Hg}}}{V_{\text{S}}}$$

$\beta_{\text{Hg}}$ : Mercury content of the sample in g/L

$V_{\text{EDTA}}$ : Added volume of EDTA standard solution in mL

$f_1$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L


$V_{\text{EP1}}$ : Volume of  $\text{ZnSO}_4$  standard solution consumed up to the first equivalence point in mL

$f_2$ : Titer of the  $\text{ZnSO}_4$  standard solution

$c_{\text{ZnSO}_4}$ : Concentration of the  $\text{ZnSO}_4$  standard solution in mol/L

$M_{\text{Hg}}$ : Molar mass of mercury; 200.590 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 61** Photometric titration curve of mercury with the Optrode at 502 nm

### Further literature on this topic is given in

- *Titration Application Note AN-T-144 – Mercury analysis using automatic photometric titration*

# Mg – Magnesium

68

## General

Along with calcium, magnesium was one of the first metals to be titrated complexometrically. The titration is generally carried out in a buffered ammoniacal solution at pH = 10 and can be performed directly. Since  $\text{Mg}^{2+}$  is very often accompanied by  $\text{Ca}^{2+}$ , please refer to the separate determination of the two metals in the section on "Water hardness".

Interfering metals, such as  $\text{Cu}^{2+}$ ,  $\text{Ni}^{2+}$ ,  $\text{Co}^{2+}$ ,  $\text{Fe}^{3+}$ , and  $\text{Mn}^{2+}$  can be eliminated by precipitating them with  $\text{Na}_2\text{S}$  solution. The sample should then be filtered and the excess sulfides destroyed by oxidizing them to sulfate (e.g., with  $\text{H}_2\text{O}_2$ ).

## Table of equivalent masses

1 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L} = x \text{ mg equivalent}$

Formula	Molar mass in g/mol	Equivalent in mg
Mg	24.305	2.431
$\text{Mg}(\text{CH}_3\text{COO})_2 \cdot 4\text{H}_2\text{O}$	214.450	21.445
$\text{MgCO}_3$	84.318	8.432
$\text{MgCl}_2$	95.211	9.521
$\text{MgCl}_2 \cdot 6\text{H}_2\text{O}$	203.302	20.330
$\text{Mg}(\text{NO}_3)_2$	148.315	14.832
$\text{Mg}(\text{NO}_3)_2 \cdot 6\text{H}_2\text{O}$	256.406	25.641
MgO	40.304	4.030
$\text{Mg}(\text{OH})_2$	58.320	5.832
$\text{MgSO}_4$	120.360	12.036
$\text{MgSO}_4 \cdot \text{H}_2\text{O}$	138.380	13.838
$\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$	246.470	24.647

## Potentiometric with the Cu-ISE

### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Cu complex solution:  $c[\text{Cu}(\text{NH}_4)_2\text{EDTA}] = 0.1 \text{ mol/L}$  (e.g., Merck No. 105217)
- Ammonia buffer, pH = 10: Dissolve 54 g  $\text{NH}_4\text{Cl}$  in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 350 mL  $w/w(\text{NH}_3) = 25\%$ . Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Ammonia:  $w/w(\text{NH}_3) = 25\%$

### Analysis

If necessary, dilute acidic sample solutions with dist.  $\text{H}_2\text{O}$  to approx. 50 mL and pre-neutralize to approx. pH 4–5 with  $w/w(\text{NH}_3) = 25\%$ . Add 0.5 mL Cu complex solution and 5 mL ammonia buffer, wait for 10–20 s while stirring, and titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.

## Calculation

$$\beta_{\text{Mg}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Mg}}}{V_{\text{S}}}$$

$\beta_{\text{Mg}}$ : Magnesium content of the sample in g/L


$V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL

$f$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L

$M_{\text{Mg}}$ : Molar mass of magnesium;  
24.305 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 62** Example of a potentiometric titration curve of magnesium with the Cu-ISE

## Further literature on this topic is given in

- *Application Bulletin AB-101 – Complexometric titrations with the Cu-ISE*
- *Titration Application Note AN-T-128 – Automatic determination of magnesium in aqueous solution with the Cu-ISE*

## Photometric with the Optrode at 610 nm

### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Ammonia buffer,  $\text{pH} = 10$ : Dissolve 54 g  $\text{NH}_4\text{Cl}$  in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 350 mL  $w/w(\text{NH}_3) = 25\%$ . Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Ammonia:  $w/w(\text{NH}_3) = 25\%$
- Indicator solution: Eriochrome Black T; dissolve 100 mg Eriochrome Black T in 100 mL ethanol.

### Analysis

If necessary, dilute acidic sample solutions with dist.  $\text{H}_2\text{O}$  to approx. 100 mL and preneutralize to approx.  $\text{pH} 4\text{--}5$  with  $w/w(\text{NH}_3) = 25\%$ . Add 5 mL ammonia buffer and 0.5 mL indicator solution, and titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.

## Calculation

$$\beta_{\text{Mg}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Mg}}}{V_{\text{S}}}$$

$\beta_{\text{Mg}}$ : Magnesium content of the sample in g/L


$V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL

$f$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L

$M_{\text{Mg}}$ : Molar mass of magnesium; 24.305 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 63** Photometric titration curve of magnesium with the Optrode at 610 nm

### Further literature on this topic is given in

- *Application Bulletin AB-125 – Simultaneous determination of calcium, magnesium, and alkalinity by complexometric titration with potentiometric or photometric indication in water and beverage samples*

## Mn – Manganese

### General

Manganese behaves very similarly to  $\text{Mg}^{2+}$  in complexometric titrations. However, account must be taken of the fact that only  $\text{Mn}^{2+}$  can be titrated. In alkaline solution,  $\text{Mn}(\text{OH})_2$  forms first and then, by reaction with atmospheric oxygen,  $\text{Mn}(\text{OH})_3$ . Therefore, some precautions must be taken when performing a photometric titration. Ascorbic acid is added to the sample solution to reduce any  $\text{Mn}^{3+}$ , as well as triethanolamine hydrochloride

to prevent precipitation of  $\text{Mn}(\text{OH})_2$ . This also masks  $\text{Fe}^{3+}$  and  $\text{Al}^{3+}$  ions.  $\text{Co}^{2+}$ ,  $\text{Ni}^{2+}$ ,  $\text{Cu}^{2+}$ ,  $\text{Zn}^{2+}$ ,  $\text{Cd}^{2+}$ ,  $\text{Hg}^{2+}$ , and the Pt-group metals can be masked with KCN in case of photometric titration.  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  are also determined. If they are present, manganese must be separated beforehand, for example, as  $\text{MnO}_2$ . ( $\text{MnO}_2$  can be reduced by oxalic acid).

**Table of equivalent masses**1 mL c(Na<sub>2</sub>EDTA) = 0.1 mol/L = x mg equivalent

Formula	Molar mass in g/mol	Equivalent in mg
Mn	54.938	5.494
MnCO <sub>3</sub>	114.974	11.497
MnCl <sub>2</sub>	125.844	12.584
MnCl <sub>2</sub> ·4H <sub>2</sub> O	197.905	19.791
Mn(NO <sub>3</sub> ) <sub>2</sub> ·4H <sub>2</sub> O	251.009	25.101
MnO	70.937	7.094
MnO <sub>2</sub>	86.937	8.694
Mn(OH) <sub>2</sub>	88.953	8.895
MnSO <sub>4</sub>	151.000	15.100
MnSO <sub>4</sub> ·H <sub>2</sub> O	169.015	16.902
MnSO <sub>4</sub> ·4H <sub>2</sub> O	223.060	22.306

**Potentiometric with the Cu ISE****Reagents**

- Titrant 1: c(Na<sub>2</sub>EDTA) = 0.1 mol/L
- Titrant 2: c(CuSO<sub>4</sub>) = 0.1 mol/L; dissolve 24.97 g CuSO<sub>4</sub>·5H<sub>2</sub>O in dist. H<sub>2</sub>O, mix, and make up to 1 L with dist. H<sub>2</sub>O.
- Ammonia buffer: Dissolve 54 g NH<sub>4</sub>Cl in approx. 200 mL dist. H<sub>2</sub>O and add 350 mL w/w(NH<sub>3</sub>) = 25%. Mix and make up to 1 L with dist. H<sub>2</sub>O.
- Ammonia: w/w(NH<sub>3</sub>) = 25%

**Analysis**

If necessary, dilute the acidic sample solution, which should not contain more than 45 mg Mn<sup>2+</sup>, with dist. H<sub>2</sub>O to approx. 50 mL and then, if necessary, preneutralize to pH 3–4 with w/w(NH<sub>3</sub>) = 25%. Add 10 mL c(Na<sub>2</sub>EDTA) = 0.1 mol/L and 5 mL ammonia buffer, and titrate with c(CuSO<sub>4</sub>) = 0.1 mol/L past the first equivalence point.

**Calculation**

$$\beta_{\text{Mn}} = \frac{(V_{\text{EDTA}} \times f_1 \times c_{\text{EDTA}} - V_{\text{EP1}} \times f_2 \times c_{\text{CuSO}_4}) \times M_{\text{Mn}}}{V_{\text{S}}}$$

$\beta_{\text{Mn}}$ : Manganese content of the sample in g/L

$V_{\text{EDTA}}$ : Added volume of EDTA standard solution in mL

$f_1$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L


$V_{\text{EP1}}$ : Volume of CuSO<sub>4</sub> standard solution consumed up to the first equivalence point in mL

$f_2$ : Titer of the CuSO<sub>4</sub> standard solution

$c_{\text{CuSO}_4}$ : Concentration of the CuSO<sub>4</sub> standard solution in mol/L

$M_{\text{Mn}}$ : Molar mass of manganese; 54.938 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 64** Potentiometric titration curve of manganese with the Cu-ISE

### Photometric with the Optrode at 610 nm

#### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Triethanolamine hydrochloride:  $w/w(\text{C}_6\text{H}_{15}\text{NO}_3 \cdot \text{HCl}) = 20\%$ ; dissolve 100 g  $\text{C}_6\text{H}_{15}\text{NO}_3 \cdot \text{HCl}$  in dist.  $\text{H}_2\text{O}$ , mix, and make up to 500 mL with dist.  $\text{H}_2\text{O}$ .
- Ascorbic acid: (vitamin C), p.a.
- Ammonia buffer: Dissolve 54 g  $\text{NH}_4\text{Cl}$  in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 350 mL  $w/w(\text{NH}_3) = 25\%$ . Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Sodium hydroxide:  $c(\text{NaOH}) = 1 \text{ mol/L}$
- Indicator solution: Eriochrome Black T; dissolve 100 mg Eriochrome Black T in 100 mL ethanol.

#### Further literature on this topic is given in

- *Application Bulletin AB-101 – Complexometric titrations with the Cu-ISE*
- *Titration Application Note AN-T-120 – Automatic manganese determination in aqueous solution using the ion-selective copper electrode (Cu-ISE)*

#### Analysis

To the acidic sample solution, add 10 mL  $w/w(\text{C}_6\text{H}_{15}\text{NO}_3 \cdot \text{HCl}) = 20\%$  and a spatula tip of ascorbic acid. If necessary, dilute to approx. 75 mL with dist.  $\text{H}_2\text{O}$ . Preneutralize with  $c(\text{NaOH}) = 1 \text{ mol/L}$  to approx. pH 7. Add 5 mL ammonia buffer and 0.25 mL indicator solution, and titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.



## Calculation

$$\beta_{\text{Mn}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Mn}}}{V_{\text{S}}}$$

$\beta_{\text{Mn}}$ : Manganese content of the sample in g/L


$V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL

$f$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L

$M_{\text{Mn}}$ : Molar mass of manganese; 54.938 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 65** Example of a photometric titration curve of manganese with the Optrode at 610 nm

## Further literature on this topic is given in

- *Titration Application Note AN-T-141 – Automatic photometric determination of manganese using the Optrode*

## Ni – Nickel

### General

Nickel can be titrated very successfully in ammoniacal solutions. However, alkaline earth ions are also determined if they are present in large quantities. High concentrations of  $\text{Al}^{3+}$ ,  $\text{Fe}^{3+}$ , and  $\text{Mn}^{2+}$  can be masked by adding triethanolamine.

$\text{Ni}^{2+}$  can be separated from interfering metal ions (exception  $\text{Pd}^{2+}$ ) by precipitating it with dimethylglyoxime. The precipitate is dissolved in conc. HCl, treated with excess  $\text{Na}_2\text{EDTA}$ , and the pH is adjusted to 10 with  $\text{NH}_3$  before being titrated with, for example, a  $\text{Zn}^{2+}$  standard solution and Eriochrome Black T indicator.

**Table of equivalent masses**1 mL c(Na<sub>2</sub>EDTA) = 0.1 mol/L = x mg equivalent

Formula	Molar mass in g/mol	Equivalent in mg
Ni	58.693	5.869
NiCO <sub>3</sub>	118.702	11.870
NiCl <sub>2</sub>	129.599	12.960
NiCl <sub>2</sub> ·6H <sub>2</sub> O	237.691	23.769
Ni(NO <sub>3</sub> ) <sub>2</sub>	182.703	18.270
Ni(NO <sub>3</sub> ) <sub>2</sub> ·6H <sub>2</sub> O	290.795	29.080
NiO	74.693	7.469
NiSO <sub>4</sub>	154.756	15.476
NiSO <sub>4</sub> ·7H <sub>2</sub> O	280.863	28.086

**Potentiometric with the Cu-ISE****Reagents**

- Titrant: c(Na<sub>2</sub>EDTA) = 0.1 mol/L
- Cu complex solution:  
c[Cu(NH<sub>4</sub>)<sub>2</sub>EDTA] = 0.1 mol/L  
(e.g., Merck No. 105217)
- Ammonia buffer: Dissolve 54 g NH<sub>4</sub>Cl in approx. 200 mL dist. H<sub>2</sub>O and add 350 mL w/w(NH<sub>3</sub>) = 25%. Mix and make up to 1 L with dist. H<sub>2</sub>O.
- Sodium hydroxide: c(NaOH) = 1 mol/L


**Analysis**

If necessary, dilute the sample solution with dist. H<sub>2</sub>O to approx. 50 mL and pre-neutralize, if necessary, to pH 4–5 with c(NaOH) = 1 mol/L. Add 5 mL ammonia buffer and 0.5 mL Cu complex solution, wait for 10–20 s while stirring, and titrate with c(Na<sub>2</sub>EDTA) = 0.1 mol/L past the first equivalence point.

**Calculation**

$$\beta_{\text{Ni}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Ni}}}{V_{\text{S}}}$$

- $\beta_{\text{Ni}}$ : Nickel content of the sample in g/L
- $V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL
- $f$ : Titer of the EDTA standard solution
- $c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L
- $M_{\text{Ni}}$ : Molar mass of nickel; 58.693 g/mol
- $V_{\text{S}}$ : Sample size in mL



**Figure 66** Potentiometric titration curve of nickel with the Cu-ISE

### Photometric with the Optrode at 574 nm

#### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Ammonia buffer: dissolve 54 g  $\text{NH}_4\text{Cl}$  and 350 mL w/w( $\text{NH}_3$ ) = 25% in dist.  $\text{H}_2\text{O}$ , mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Sodium hydroxide:  $c(\text{NaOH}) = 1 \text{ mol/L}$
- Indicator solution: murexide; dissolve 0.2 g murexide (mixed 1:100 with NaCl) in 50 mL dist.  $\text{H}_2\text{O}$ .

#### Analysis

If necessary, dilute the sample solution with dist.  $\text{H}_2\text{O}$  to approx. 50 mL and then preneutralize, if necessary, to pH 4–5 with  $c(\text{NaOH}) = 1 \text{ mol/L}$ . Add 5 mL ammonia buffer and 5 mL indicator solution, and then titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point. Titrate slowly when close to the color change.


#### Further literature on this topic is given in

- *Application Bulletin AB-101 – Complexometric titrations with the Cu-ISE*
- *Titration Application Note AN-T-129 – Automatic determination of nickel in aqueous solution with the Cu-ISE*

#### Calculation

$$\beta_{\text{Ni}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Ni}}}{V_{\text{S}}}$$

- $\beta_{\text{Ni}}$ : Nickel content of the sample in g/L
- $V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL
- f: Titer of the EDTA standard solution
- $c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L
- $M_{\text{Ni}}$ : Molar mass of nickel; 58.693 g/mol
- $V_{\text{S}}$ : Sample size in mL



**Figure 67** Example of a photometric titration curve of nickel with the Optrode at 574 nm

### Further literature on this topic is given in

- *Titration Application Note AN-T-150 – Determination of nickel using automatic photometric titration*

## Pb – Lead

### General

Complexometric titration of  $\text{Pb}^{2+}$  can be carried out directly, via displacement reactions (e.g., with  $\text{MgEDTA}$ ), or by back titration in alkaline solution. Sulfates should not be present, if possible, because  $\text{Pb}^{2+}$  precipitates as  $\text{PbSO}_4$  in their presence. The alkaline solution should be adjusted to pH 10 with  $\text{NaOH}$ , and tartrate should be used as an auxiliary complexing agent. This prevents precipitation of  $\text{Pb}^{2+}$  as  $\text{Pb}(\text{OH})_2$ . However, the tartrate excess must not be too large otherwise

the effective stability constant of the indicator (Eriochrome Black T) or of the  $\text{PbEDTA}$  complex would be lowered too much.

Furthermore, photometric titration in slightly acidic medium enables determination of Pb in the presence of phosphate and alkaline earth ions and is also suitable for determinations in low concentration ranges.

**Table of equivalent masses**1 mL c(Na<sub>2</sub>EDTA) = 0.1 mol/L = x mg equivalent

Formula	Molar mass in g/mol	Equivalent in mg
Pb	207.200	20.720
Pb(CH <sub>3</sub> COO) <sub>2</sub>	325.288	32.529
Pb(CH <sub>3</sub> COO) <sub>2</sub> ·3H <sub>2</sub> O	379.334	37.933
Pb(CH <sub>2</sub> CH <sub>2</sub> ) <sub>4</sub>	323.444	32.344
PbCO <sub>3</sub>	267.209	26.721
PbCl <sub>2</sub>	278.106	27.811
PbO	223.199	22.320
Pb(OH) <sub>2</sub>	241.215	24.122
Pb <sub>3</sub> (PO <sub>4</sub> ) <sub>2</sub>	811.543	27.051
PbSO <sub>4</sub>	303.263	30.326
Pb(NO <sub>3</sub> ) <sub>2</sub>	331.210	33.121

**Potentiometric with the Cu-ISE****Reagents**

- Titrant: c(Na<sub>2</sub>EDTA) = 0.1 mol/L
- Cu complex solution: c[Cu(NH<sub>4</sub>)<sub>2</sub>EDTA] = 0.1 mol/L (e.g., Merck No. 105217)
- Ammonia buffer: Dissolve 54 g NH<sub>4</sub>Cl in approx. 200 mL dist. H<sub>2</sub>O and add 350 mL w/w(NH<sub>3</sub>) = 25%. Mix and make up to 1 L with dist. H<sub>2</sub>O.
- Sodium hydroxide: c(NaOH) = 2 mol/L
- Ammonium tartrate solution: β((NH<sub>4</sub>)<sub>2</sub>C<sub>4</sub>H<sub>4</sub>O<sub>6</sub>) = 250 g/L


**Analysis**

If necessary, preneutralize acidic sample solutions to pH 4–5 with c(NaOH) = 2 mol/L. Add 50 mL ammonium tartrate, 1 mL ammonia buffer, 0.5 mL Cu complex solution, and 5 mL c(NaOH) = 2 mol/L, wait for 30 s while stirring, and then titrate with c(Na<sub>2</sub>EDTA) = 0.1 mol/L past the first equivalence point.

**Calculation**

$$\beta_{\text{Pb}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Pb}}}{V_{\text{S}}}$$

- β<sub>Pb</sub>: Lead content of the sample in g/L  
V<sub>EP1</sub>: Volume of EDTA standard solution consumed up to the first equivalence point in mL  
f: Titer of the EDTA standard solution  
c<sub>EDTA</sub>: Concentration of the EDTA standard solution in mol/L  
M<sub>Pb</sub>: Molar mass of lead; 207.200 g/mol  
V<sub>S</sub>: Sample size in mL



**Figure 68** Potentiometric titration curve of lead with the Cu-ISE

## Photometric with the Optrode at 574 nm

### Reagents

- Titrant 1:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Titrant 2:  $c(\text{ZnSO}_4) = 0.1 \text{ mol/L}$ ; dissolve 28.8 g  $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$  in dist.  $\text{H}_2\text{O}$ , mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Acetate buffer: Dissolve 123 g sodium acetate in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 50 mL glacial acetic acid. Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Sodium hydroxide:  $c(\text{NaOH}) = 1 \text{ mol/L}$
- Indicator solution: Xylenol orange; dissolve 50 mg xylenol orange disodium salt in 50 mL dist.  $\text{H}_2\text{O}$ .

### Analysis

If necessary, dilute the sample solution with dist.  $\text{H}_2\text{O}$  to approx. 50 mL and adjust, if necessary, to pH 5–6 with  $c(\text{NaOH}) = 1 \text{ mol/L}$ . Add 10 mL acetate buffer, 5 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$ , and 0.5 mL indicator solution. Allow to react for 1 min, and titrate with  $c(\text{ZnSO}_4) = 0.1 \text{ mol/L}$  past the first equivalence point.

### Further literature on this topic is given in

- *Application Bulletin AB-101 – Complexometric titrations with the Cu-ISE*
- *Titration Application Note AN-T-130 – Automatic determination of lead in aqueous solution with the Cu-ISE*

### Calculation

$$\beta_{\text{Pb}} = \frac{(V_{\text{EDTA}} \times f_1 \times c_{\text{EDTA}} - V_{\text{EP1}} \times f_2 \times c_{\text{ZnSO}_4}) \times M_{\text{Pb}}}{V_{\text{S}}}$$

$\beta_{\text{Pb}}$ : Lead content of the sample in g/L

$V_{\text{EDTA}}$ : Added volume of EDTA standard solution in mL

$f_1$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L


$V_{\text{EP1}}$ : Volume of  $\text{ZnSO}_4$  standard solution consumed up to the first equivalence point in mL

$f_2$ : Titer of the  $\text{ZnSO}_4$  standard solution

$c_{\text{ZnSO}_4}$ : Concentration of the  $\text{ZnSO}_4$  standard solution in mol/L

$M_{\text{Pb}}$ : Molar mass of lead; 207.200 g/mol

$V_{\text{S}}$ : Sample size in mL



Further literature on this topic is given in

- Titration Application Note AN-T-140  
– Automated photometric determination of lead using the Optrode

**Figure 69** Example of a photometric titration curve of lead with the Optrode at 574 nm

## Pd – Palladium

### General

There is relatively little data in the literature on complexometric titration of the platinum group metals.  $\text{Pd}^{2+}$  is usually determined by back titration of excess EDTA. This excess is titrated with a suitable standard solution, for example, with  $\text{Zn}^{2+}$  in alkaline solution and Eriochrome Black T indicator. (In this case, the solution must be made alkaline with NaOH because  $\text{Pd}^{2+}$  reacts with ammonia to form very stable, non-titratable Pd-amine complexes). The excess EDTA is then

back-titrated with xylenol orange indicator and  $\text{Th}^{4+}$  at pH = 3.0 or  $\text{Ti}^{3+}$  at pH 4 to 5. In our example, we titrate at pH 4–5 with xylenol orange indicator using the less toxic  $\text{Zn}^{2+}$  as the titrant.

Since Pt forms very stable complexes with chloride ions,  $\text{Pd}^{2+}$  can be determined alongside Pt without interference in the presence of chloride ions.

### Table of equivalent masses

1 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L} = x \text{ mg equivalent}$

Formula	Molar mass in g/mol	Equivalent in mg
Pd	106.420	10.642
$\text{PdCl}_2$	177.326	17.733
$\text{PdCl}_2 \cdot 2\text{H}_2\text{O}$	213.357	21.336
$\text{PdCl}_4$	248.232	24.823
$\text{PdCl}_6$	319.138	31.914
$\text{Pd}(\text{NO}_3)_2$	230.430	23.043
$\text{PdO}$	122.419	12.242
$\text{PdSO}_4$	202.483	20.248
$\text{PdSO}_4 \cdot 2\text{H}_2\text{O}$	238.513	23.851

## Photometric with the Optrode at 610 nm

### Reagents

- Titrant 1:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Titrant 2:  $c(\text{ZnSO}_4) = 0.1 \text{ mol/L}$ ;  
dissolve 28.8 g  $\text{ZnSO}_4 \cdot 7\text{H}_2\text{O}$  in approx. 500 mL dist.  $\text{H}_2\text{O}$ . Add 0.5 mL  $w/w(\text{H}_2\text{SO}_4) = 96\%$ , mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ . Mix 100 mL of this solution, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Acetate buffer,  $\text{pH} = 4.7$ : Dissolve 123 g sodium acetate in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 50 mL glacial acetic acid. Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Indicator solution: Xylenol orange; dissolve 100 mg xylenol orange disodium salt in 100 mL dist.  $\text{H}_2\text{O}$ .

### Analysis

If necessary, dilute the acidic sample solution, which should not contain more than 20 mg  $\text{Pd}^{2+}$ , with dist.  $\text{H}_2\text{O}$  to approx. 100 mL and degas. Add 10 mL acetate buffer, 5 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$ , and 10 mL indicator solution, and titrate with  $c(\text{ZnSO}_4) = 0.1 \text{ mol/L}$  past the first equivalence point.

### Calculation

$$\beta_{\text{Pd}} = \frac{(V_{\text{EDTA}} \times f_1 \times c_{\text{EDTA}} - V_{\text{EP1}} \times f_2 \times c_{\text{ZnSO}_4}) \times M_{\text{Pd}}}{V_{\text{S}}}$$

$\beta_{\text{Pd}}$ : Palladium content of the sample in g/L

$V_{\text{EDTA}}$ : Added volume of EDTA standard solution in mL

$f_1$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L


$V_{\text{EP1}}$ : Volume of  $\text{ZnSO}_4$  standard solution consumed up to the first equivalence point in mL

$f_2$ : Titer of the  $\text{ZnSO}_4$  standard solution

$c_{\text{ZnSO}_4}$ : Concentration of the  $\text{ZnSO}_4$  standard solution in mol/L

$M_{\text{Pd}}$ : Molar mass of palladium; 106.420 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 70** Photometric titration curve of palladium with the Optrode at 610 nm

### Further literature on this topic is given in

- *Titration Application Note AN-T-145 – Determination of palladium using automatic photometric titration*



# Sn – Tin

## General

Both valence states of tin (II and IV) form stable complexes, for example with EDTA, which can be titrated in acidic solutions owing to their high complex formation constants ( $\log K_f = \text{approx. } 22$ ).

Furthermore, acidic solutions should also be used because tin ions form stable hydroxo complexes that are difficult or

impossible to titrate. Preneutralization must be carried out carefully with vigorous stirring.

Interfering  $\text{Pb}^{2+}$  ions are usually precipitated as  $\text{PbSO}_4$  and then removed from the solution by filtration.

## Table of equivalent masses

1 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L} = x \text{ mg equivalent}$

Formula	Molar mass in g/mol	Equivalent in mg
Sn	118.710	11.871
$\text{SnCl}_2$	189.616	18.962
$\text{SnCl}_2 \cdot 2\text{H}_2\text{O}$	225.647	22.565
$\text{SnCl}_4$	260.522	26.052
$\text{SnCl}_4 \cdot 5\text{H}_2\text{O}$	350.598	35.060
$\text{SnF}_2$	156.707	15.671
$\text{SnF}_4$	194.704	19.470
SnO	134.709	13.471
$\text{SnO}_2$	150.709	15.071
$\text{Sn}(\text{OH})_2$	152.725	15.273
$\text{Sn}(\text{OH})_4$	186.739	18.674
$\text{SnSO}_4$	214.773	21.477

## Photometric with the Optrode at 574 nm

### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Buffer:  $\text{pH} = 2.1$ ; dissolve 5.88 g citric acid and 3.58 g NaCl in 700 mL  $c(\text{HCl}) = 0.1 \text{ mol/L}$ , and then adjust to  $\text{pH } 2.1$  with  $c(\text{NaOH}) = 2 \text{ mol/L}$ . Transfer the solution into a 1 L volu-

metric flask and make up to the mark with dist.  $\text{H}_2\text{O}$ .

- Indicator solution: Xylenol orange; dissolve 25 mg xylenol orange disodium salt in 100 mL dist.  $\text{H}_2\text{O}$ .

## Analysis

To the very acidic sample solution, which should not contain more than 100 mg Sn<sup>2+</sup> or Sn<sup>4+</sup>, add 10 mL of buffer pH 2.1 and 80 mL deion. H<sub>2</sub>O. Next, add 0.5 mL indicator solution, and titrate with c(Na<sub>2</sub>EDTA) = 0.1 mol/L past the first equivalence point.

## Calculation

$$\beta_{\text{Sn}} = \frac{V_{\text{EP1}} \times f \times C_{\text{EDTA}} \times M_{\text{Sn}}}{V_{\text{S}}}$$

$\beta_{\text{Sn}}$ : Tin content of the sample in g/L


$V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL

$f$ : Titer of the EDTA standard solution

$C_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L

$M_{\text{Sn}}$ : Molar mass of tin;  
118.710 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 71** Example of a photometric titration curve of tin with the Optrode at 574 nm

## Further literature on this topic is given in

- *Titration Application Note AN-T-146 – Automatic determination of tin using photometric titration*

# Th – Thorium

## General

Thorium reacts with complexing agents to form very strong complexes with high complex formation constants. Thus it can even be titrated at pH < 1, for example, with Na<sub>2</sub>EDTA and xylenol orange indicator. Divalent cations (except Hg<sup>2+</sup>) do not interfere at such low pH values. Some

authors (Pribil, R., Burger, K., Milner, G. W. C., Edwards, J. W.) have pointed out that thorium can form even stronger complexes with sulfates so that EDTA can be liberated again at pH > 1. They exploited this situation to determine Th<sup>4+</sup> in the presence of Zr<sup>4+</sup>. The sum of Th<sup>4+</sup> and Zr<sup>4+</sup>

is determined in the first step by back-titrating excess EDTA at pH 2.5–2.7 with a  $\text{Bi}^{3+}$  standard solution and xylenol orange indicator. Next, the fully titrated solution is acidified with  $\text{HNO}_3$  to pH 1.2–1.3, and

then approx. 2 g ammonium sulfate is added. The thus liberated EDTA molecules are subsequently titrated with  $\text{Bi}^{3+}$  standard solution and the same indicator.

### Table of equivalent masses

1 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L} = x \text{ mg equivalent}$

Formula	Molar mass in g/mol	Equivalent in mg
Th	232.038	23.204
$\text{Th}(\text{NO}_3)_4 \cdot 5\text{H}_2\text{O}$	570.134	57.013
$\text{ThO}_2$	264.037	26.404
$\text{Th}(\text{OH})_4$	300.067	30.007
$\text{Th}(\text{SO}_4)_2$	424.150	42.415
$\text{Th}(\text{SO}_4)_2 \cdot 9\text{H}_2\text{O}$	586.290	58.629

### Photometric with the Optrode at 574 nm

#### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Acetate buffer: add 123 g sodium acetate and 50 mL  $w/w(\text{NH}_3) = 25\%$  to a 1 L volumetric flask, mix, and make up to the mark with dist.  $\text{H}_2\text{O}$ .
- Indicator solution: Xylenol orange; dissolve 50 mg xylenol orange disodium salt in 50 mL dist.  $\text{H}_2\text{O}$ .


#### Analysis

If necessary, dilute the sample solution with dist.  $\text{H}_2\text{O}$  to approx. 100 mL. Degas the mixture, and add 0.5 mL indicator solution and 10 mL acetate buffer. Titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.

#### Calculation

$$\beta_{\text{Th}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Th}}}{V_{\text{S}}}$$

- $\beta_{\text{Th}}$ : Thorium content of the sample in g/L
- $V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL
- $f$ : Titer of the EDTA standard solution
- $c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L
- $M_{\text{Th}}$ : Molar mass of thorium; 232.038 g/mol
- $V_{\text{S}}$ : Sample size in mL



**Figure 72** Photometric titration curve of thorium with the Optrode at 574 nm

### Further literature on this topic is given in

- *Titration Application Note AN-T-149*  
– *Determination of thorium using automatic photometric titration*

## Tl – Thallium

### General

$Tl^+$  forms very weak complexes, whereas  $Tl^{3+}$  forms strong complexes with a high complex formation constant. This metal is thus titrated in its trivalent state. The determinations are carried out in slightly acidic solution so that there is no interference, for example, from any alkaline earth ions that may be present.

Thallium is often present in its monovalent state ( $Tl^+$ ). The best method of oxidizing it to  $Tl^{3+}$  is to use aqua regia

(mixture of  $HNO_3$  and  $HCl$ ) because nitric acid alone is insufficient. The sample is treated with aqua regia and then evaporated until it is almost dry.

$Tl^+$  can be separated well from interfering metal ions by adding some potassium iodide solution to the sample solution. The thus precipitated  $TlI$  is removed by filtration and then oxidized to  $Tl^{3+}$  with aqua regia as described above.

**Table of equivalent masses**

1 mL c(Na<sub>2</sub>EDTA) = 0.1 mol/L = x mg equivalent

Formula	Molar mass in g/mol	Equivalent in mg
Tl	204.383	20.438
Tl <sub>2</sub> CO <sub>3</sub>	468.775	23.439
TlCl	239.836	23.439
TlCl <sub>3</sub>	310.742	31.074
TlI	331.287	33.129
TlNO <sub>3</sub>	266.387	26.639
Tl(NO <sub>3</sub> ) <sub>3</sub>	444.765	44.477
Tl <sub>2</sub> O	424.762	21.238
Tl <sub>2</sub> O <sub>3</sub>	456.764	22.838
Tl <sub>2</sub> SO <sub>4</sub>	504.829	25.241
Tl <sub>2</sub> (SO <sub>4</sub> ) <sub>3</sub> ·7H <sub>2</sub> O	823.061	41.153

**Potentiometric with the Cu-ISE****Reagents**

- Titrant 1: c(Na<sub>2</sub>EDTA) = 0.1 mol/L
- Titrant 2: c(CuSO<sub>4</sub>) = 0.1 mol/L  
dissolve 24.97 g CuSO<sub>4</sub>·5 H<sub>2</sub>O in approx. 500 mL dist. H<sub>2</sub>O, add 0.5 mL w/w(H<sub>2</sub>SO<sub>4</sub>) = 96%, mix, and make up to 1 L with dist. H<sub>2</sub>O.
- Acetate buffer, pH = 4.7: Dissolve 116 g ammonium acetate in approx. 200 mL dist. H<sub>2</sub>O and add 86 mL glacial acetic acid. Mix and make up to 1 L with dist. H<sub>2</sub>O.
- Ammonia: w/w(NH<sub>3</sub>) = 25%

**Analysis**

If necessary, dilute the acidic sample solution, which should not contain more than 150 mg Tl<sup>3+</sup>, with dist. H<sub>2</sub>O to approx. 50 mL. Add 5 mL acetate buffer and 10 mL c(Na<sub>2</sub>EDTA) = 0.1 mol/L, and then, if necessary, adjust to pH 4.5–4.7 with w/w(NH<sub>3</sub>) = 25%. Next, titrate with c(CuSO<sub>4</sub>) = 0.1 mol/L past the first equivalence point.

**Calculation**

$$\beta_{\text{Tl}} = \frac{(V_{\text{EDTA}} \times f_1 \times c_{\text{EDTA}} - V_{\text{EP1}} \times f_2 \times c_{\text{CuSO}_4}) \times M_{\text{Tl}}}{V_{\text{S}}}$$

$\beta_{\text{Tl}}$ : Thallium content of the sample in g/L

$V_{\text{EDTA}}$ : Added volume of EDTA standard solution in mL

$f_1$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L


$V_{\text{EP1}}$ : Volume of CuSO<sub>4</sub> standard solution consumed up to the first equivalence point in mL

$f_2$ : Titer of the CuSO<sub>4</sub> standard solution

$c_{\text{CuSO}_4}$ : Concentration of the CuSO<sub>4</sub> standard solution in mol/L

$M_{\text{Tl}}$ : Molar mass of thallium; 204.383 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 73** Potentiometric titration curve of thallium with the Cu-ISE

#### Further literature on this topic is given in

- *Application Bulletin AB-101 – Complexometric titrations with the Cu-ISE*
- *Titration Application Note AN-T-122 – Automatic thallium determination in aqueous solution using the ion-selective copper electrode (Cu-ISE)*

#### Photometric with the Optrode at 574 nm

##### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Acetate buffer,  $\text{pH} = 4.9$ : Dissolve 123 g sodium acetate in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 50 mL glacial acetic acid. Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Indicator solution: Xylenol orange; dissolve 50 mg xylenol orange disodium salt in 50 mL  $\text{H}_2\text{O}$ .

#### Analysis

To the acidic sample solution, add 5 mL acetate buffer and, if necessary, dilute with dist.  $\text{H}_2\text{O}$  to approx. 100 mL. Degas the mixture, add 0.5 mL indicator solution, and titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.

#### Calculation

$$\beta_{\text{Tl}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Tl}}}{V_{\text{S}}}$$

$\beta_{\text{Tl}}$ : Thallium content of the sample in g/L


$V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL

$f$ : Titer of the EDTA standard solution

$c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L

$M_{\text{Tl}}$ : Molar mass of thallium; 204.383 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 74** Example of a photometric titration curve of thallium with the Optrode at 610 nm

#### Metrohm Literature on this topic

- *Titration Application Note AN-T-147 – Determination of thallium using automatic photometric titration*

# Zn – Zinc

## General

The complexation behavior of  $\text{Zn}^{2+}$  is very similar to that of cadmium. It can be titrated in alkaline solution without problems. However,  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  interfere if present because they are also determined. In alkaline solution, a possible strategy is to use a first titration to determine the total sum of  $\text{Zn}^{2+}$  plus  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  with Eriochrome Black T indicator. Next, the  $\text{Zn}^{2+}$  is masked with KCN and the liberated

$\text{Na}_2\text{EDTA}$  is determined by means of a consecutive titration, for example, back titration with  $\text{MgSO}_4$  standard solution.

However, it is much simpler to carry out the (photometric) titration at approx. pH 5 because  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$  no longer interfere.

## Table of equivalent masses

1 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L} = x \text{ mg equivalent}$

Formula	Molar mass in g/mol	Equivalent in mg
Zn	65.409	6.541
$\text{Zn}(\text{CH}_3\text{COO})_2$	183.497	18.350
$\text{Zn}(\text{CH}_3\text{COO})_2 \cdot 2 \text{ H}_2\text{O}$	219.528	21.953
$\text{Zn}(\text{CN})_2$	117.444	11.744
$\text{ZnCO}_3$	125.418	12.542
$\text{ZnCl}_2$	136.315	13.632
$\text{Zn}(\text{NO}_3)_2$	189.419	18.942
$\text{Zn}(\text{NO}_3)_2 \cdot 6 \text{ H}_2\text{O}$	297.510	29.751
ZnO	81.408	8.142
$\text{Zn}(\text{OH})_2$	99.424	9.942
$\text{ZnSO}_4$	161.472	16.147
$\text{ZnSO}_4 \cdot 7 \text{ H}_2\text{O}$	287.579	28.758

## Potentiometric with the Cu-ISE

### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Ammonia buffer: Dissolve 54 g  $\text{NH}_4\text{Cl}$  in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 350 mL  $w/w(\text{NH}_3) = 25\%$ . Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Cu complex solution:  $c[\text{Cu}(\text{NH}_4)_2\text{EDTA}] = 0.1 \text{ mol/L}$  (e.g., Merck No. 105217)
- Sodium hydroxide:  $c(\text{NaOH}) = 1 \text{ mol/L}$

## Analysis

If necessary, preneutralize acidic sample solutions with  $c(\text{NaOH}) = 1 \text{ mol/L}$  to approx. pH 5–7 and, if necessary, make up to approx. 50 mL with dist.  $\text{H}_2\text{O}$ . Add 5 mL buffer solution and 1 mL Cu complex solution, wait for 20 s while stirring, and titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.

## Calculation


$$\beta_{\text{Zn}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Zn}}}{V_{\text{c}}}$$

- $\beta_{\text{Zn}}$ : Zinc content of the sample in g/L  
 $V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL  
 $f$ : Titer of the EDTA standard solution  
 $c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L  
 $M_{\text{Zn}}$ : Molar mass of zinc;  
 65.409 g/mol  
 $V_{\text{c}}$ : Sample size in mL

## Photometric with the Optrode at 610 nm

### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Ammonia buffer, pH = 10: dissolve 54 g  $\text{NH}_4\text{Cl}$  in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 350 mL w/w( $\text{NH}_3$ ) = 25%. Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Indicator solution: Eriochrome Black T; dissolve 100 mg each of Eriochrome Black T and ascorbic acid in 100 mL dist.  $\text{H}_2\text{O}$ .
- Sodium hydroxide:  $c(\text{NaOH}) = 1 \text{ mol/L}$



**Figure 75** Potentiometric titration curve of zinc with the Cu-ISE

### Further literature on this topic is given in

- *Application Bulletin AB-101 – Complexometric titrations with the Cu-ISE*
- *Titration Application Note AN-T-108 – Fully automated determination of zinc(II) in aqueous solution*

## Analysis

If necessary, preneutralize acidic sample solutions with  $c(\text{NaOH}) = 1 \text{ mol/L}$  to approx. pH 5 and dilute, if necessary, with dist.  $\text{H}_2\text{O}$  to approx. 50 mL. Add 5 mL ammonia buffer and 0.25 mL indicator solution, and titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.



### Calculation

$$\beta_{\text{Zn}} = \frac{V_{\text{EP1}} \times f \times C_{\text{EDTA}} \times M_{\text{Zn}}}{V_{\text{S}}}$$

$\beta_{\text{Zn}}$ : Zinc content of the sample in g/L


$V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL

$f$ : Titer of the EDTA standard solution

$C_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L

$M_{\text{Zn}}$ : Molar mass of zinc;  
65.409 g/mol

$V_{\text{S}}$ : Sample size in mL



**Figure 76** Example of a photometric titration curve of zinc with the Optrode at 610 nm

### Metrohm Literature on this topic

- *Titration Application Note AN-T-090 – Photometric EDTA titration of zinc sulfate according to Ph. Eur. and USP*

## Zr – Zirconium and Hf – Hafnium

### General

Similar to  $\text{Al}^{3+}$ , zirconium has a tendency to form hydroxo complexes and must therefore be titrated in acidic solution.  $\text{Fe}^{3+}$  ions interfere, but they can be reduced to non-interfering  $\text{Fe}^{2+}$  ions by adding ascorbic acid (vitamin C) (only applies to photometric titration).

Hafnium behaves exactly like  $\text{Zr}^{4+}$  and is titrated using the same procedures. Titrimetric separation of the two metals is not possible.

The EDTA complexes of both metals are relatively stable ( $\log K_f \text{ Zr} = 19.4 / \log K_f \text{ Hf} = 19.1$ ).

**Table of equivalent masses**

1 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L} = x \text{ mg equivalent}$

Formula	Molar mass in g/mol	Equivalent in mg
Hf	178.490	17.849
HfCl <sub>4</sub>	320.302	32.030
Hf(NO <sub>3</sub> ) <sub>4</sub>	426.510	42.651
HfO <sub>2</sub>	210.489	21.049
HfSO <sub>4</sub>	274.553	27.455
Zr	91.224	9.122
ZrCl <sub>4</sub>	233.036	23.304
ZrOCl <sub>2</sub> ·8H <sub>2</sub> O	322.252	32.225
Zr(NO <sub>3</sub> ) <sub>4</sub>	339.244	33.924
Zr(NO <sub>3</sub> ) <sub>4</sub> ·5H <sub>2</sub> O	429.320	42.932
ZrO <sub>2</sub>	123.223	12.322
Zr(OH) <sub>4</sub>	159.253	15.925
Zr(SO <sub>4</sub> ) <sub>2</sub>	283.349	28.335
Zr(SO <sub>4</sub> ) <sub>2</sub> ·4H <sub>2</sub> O	355.410	35.541

**Potentiometric with the Cu-ISE****Reagents**

- Titrant 1:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Titrant 2:  $c(\text{CuSO}_4) = 0.1 \text{ mol/L}$ ; dissolve 24.97 g  $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$  in approx. 500 mL dist.  $\text{H}_2\text{O}$ . Add 0.5 mL  $w/w(\text{H}_2\text{SO}_4) = 96\%$ , mix, and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Acetate buffer, pH = 4.7: Dissolve 116 g ammonium acetate in approx. 200 mL dist.  $\text{H}_2\text{O}$  and add 86 mL glacial acetic acid. Mix and make up to 1 L with dist.  $\text{H}_2\text{O}$ .
- Ammonia:  $w/w(\text{NH}_3) = 25\%$

**Analysis**

If necessary, dilute the acidic sample solution, which should not contain more than 70 mg  $\text{Zr}^{4+}$  and/or 130 mg  $\text{Hf}^{4+}$ , with dist.  $\text{H}_2\text{O}$  to approx. 50 mL. Add 5 mL acetate buffer and 10 mL  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$ , and, if necessary, adjust to pH 4.7 with  $w/w(\text{NH}_3) = 25\%$ . Allow to react for approx. 1 min while stirring, and titrate with  $c(\text{CuSO}_4) = 0.1 \text{ mol/L}$  past the first equivalence point.

## Calculation

$$\beta_{\text{Zr/Hf}} = \frac{(V_{\text{EDTA}} \times f_1 \times c_{\text{EDTA}} - V_{\text{EP1}} \times f_2 \times c_{\text{CuSO}_4}) \times M_{\text{Zr/Hf}}}{V_S}$$

- $\beta_{\text{Zr/Hf}}$ : Zirconium or hafnium content of the sample in g/L
- $V_{\text{EDTA}}$ : Added volume of EDTA standard solution in mL
- $f_1$ : Titer of the EDTA standard solution
- $c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L
- $V_{\text{EP1}}$ : Volume of  $\text{CuSO}_4$  standard solution consumed up to the first equivalence point in mL
- $f_2$ : Titer of the  $\text{CuSO}_4$  standard solution
- $c_{\text{CuSO}_4}$ : Concentration of the  $\text{CuSO}_4$  standard solution in mol/L
- $M_{\text{Zr/Hf}}$ : Molar mass of zirconium and hafnium; 91.224 g/mol and 178.490 g/mol, respectively
- $V_S$ : Sample size in mL


### Further literature on this topic is given in

- *Application Bulletin AB-101 – Complexometric titrations with the Cu-ISE*
- *Titration Application Note AN-T-123 – Automatic zirconium determination in aqueous solution with ion-selective copper electrode (Cu-ISE)*

## Photometric with the Optrode at 520 nm

### Reagents

- Titrant:  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$
- Buffer pH = 1 (glycine / hydrochloric acid): e.g., Merck 109881
- Indicator solution: Eriochrome Cyanine; dissolve 40 mg Eriochrome Cyanine in 100 mL dist.  $\text{H}_2\text{O}$ .



**Figure 77** Example of a potentiometric titration curve of zirconium with the Cu-ISE


### Analysis

To the acidic sample solution, which should not contain more than 70 mg  $\text{Zr}^{4+}$  and/or 130 mg  $\text{Hf}^{4+}$ , add 10 mL of buffer pH 1 and 1 mL indicator solution. To the solution, add 80 mL dist.  $\text{H}_2\text{O}$  and titrate with  $c(\text{Na}_2\text{EDTA}) = 0.1 \text{ mol/L}$  past the first equivalence point.

## Calculation

$$\beta_{\text{Zr/Hf}} = \frac{V_{\text{EP1}} \times f \times c_{\text{EDTA}} \times M_{\text{Zr/Hf}}}{V_{\text{S}}}$$

- $\beta_{\text{Zr/Hf}}$ : Zirconium or hafnium content of the sample in g/L  
 f: Titer of the EDTA standard solution  
 $c_{\text{EDTA}}$ : Concentration of the EDTA standard solution in mol/L  
 $V_{\text{EP1}}$ : Volume of EDTA standard solution consumed up to the first equivalence point in mL  
 $M_{\text{Zr/Hf}}$ : Molar mass of zirconium and hafnium; 91.224 g/mol and 178.490 g/mol, respectively  
 $V_{\text{S}}$ : Sample size in mL



**Figure 78** Photometric titration curve of zirconium with the Optrode at 520 nm

### Further literature on this topic is given in

- *Titration Application Note AN-T-148 – Determination of zirconium using automatic photometric titration*

# Appendix

## Index

Species	Remarks	Sensor	Pages
Al	General	Cu-ISE	40–41
Al	General	Optrode at 610 nm	41–42
Al and Mg	Alloys, antacids	Cu-ISE	29–31
Al, Ca, Fe, Mg	Cement, clinker	Optrode at 610 nm	33–37
Ba	General	Cu-ISE	44
Ba	General	Optrode at 574 nm	45
Bi	General	Cu-ISE	46–47
Bi	General	Optrode at 520 nm	47–48
Ca	General	Ca-ISE	49–50
Ca	General	Cu-ISE	50
Ca	General	Optrode at 610 nm	51
Ca, Al, Fe, Mg	Cement, clinker	Optrode at 610 nm	33–37
Cd	General	Cu-ISE	52–53
Cd	General	Optrode at 610 nm	53–54
Co	General	Cu-ISE	55
Co	General	Optrode at 574 nm	56
Cu	General	Cu-ISE	58
Cu	General	Optrode at 520 nm	59
Fe	General	Cu-ISE	61
Fe	General	Optrode at 610 nm	62
Fe, Al, Ca, Mg	Cement, clinker	Optrode at 610 nm	33–37
Ga	General	Cu-ISE	63–64
Ga	General	Optrode at 610 nm	65
Hf	General	Cu-ISE	90–91
Hf	General	Optrode at 520 nm	91–92
Hg	General	Optrode at 502 nm	67
In	General	Cu-ISE	63–64
In	General	Optrode at 610 nm	65
Mg	General	Cu-ISE	68–69
Mg	General	Optrode at 610 nm	69–70

Species	Remarks	Sensor	Pages
Mg and Al	Alloys, antacids	Cu-ISE	29–31
Mg and Zn	General	Cu-ISE	31–33
Mg, Al, Ca, Fe	Cement, clinker	Optrode at 610 nm	33–37
Mn	General	Cu-ISE	71–72
Mn	General	Optrode at 610 nm	72–73
Ni	General	Cu-ISE	74–75
Ni	General	Optrode at 574 nm	75–76
Pb	General	Cu-ISE	77–78
Pb	General	Optrode at 574 nm	78–79
Pd	General	Optrode at 610 nm	80
Sn	General	Optrode at 574 nm	81–82
Sr	General	Cu-ISE	44
Sr	General	Optrode at 574 nm	45
Sulfate	General	Ca-ISE	38–39
Th	General	Optrode at 574 nm	83–84
Tl	General	Cu-ISE	85–86
Tl	General	Optrode at 574 nm	86
Water hardness	Ca and Mg	Ca-ISE	25–26
Water hardness	Ca and Mg	Cu-ISE	26–27
Water hardness	Ca and Mg	Optrode at 610 nm	27–29
Zn	General	Cu-ISE	87–88
Zn	General	Optrode at 610 nm	88–89
Zn and Mg	General	Cu-ISE	31–33
Zr	General	Cu-ISE	90–91
Zr	General	Optrode at 520 nm	91–92

### Photometric indication of the titration equivalence point

Species	Indicator	Wavelength
Al	Xylenol orange	610 nm
Al in cement	Xylenol orange	610 nm
Ba	Phthalein purple	574 nm
Bi	Xylenol orange	520 nm
Ca	HHSNN	610 nm
Ca in cement	Murexide	610 nm
Cd	Eriochrome Black T	610 nm
Co	Murexide	574 nm
Cu	PAN	520 nm
Fe	Sulfosalicylic acid	610 nm
Fe in cement	Sulfosalicylic acid	610 nm
Ga	Xylenol orange	610 nm
Hf	Eriochrome Cyanine	520 nm
Hg	Eriochrome Black T	502 nm
In	Xylenol orange	610 nm
Mg	Eriochrome Black T	610 nm
Mg in cement	Methylthymol blue	610 nm
Mn	Eriochrome Black T	610 nm
Ni	Murexide	574 nm
Pb	Xylenol orange	574 nm
Pd	Xylenol orange	610 nm
Sn	Xylenol orange	574 nm
Sr	Phthalein purple	574 nm
Th	Xylenol orange	574 nm
Tl	Xylenol orange	574 nm
Water, calcium hardness	HHSNN	610 nm
Water, total hardness	Eriochrome Black T	610 nm
Zn	Eriochrome Black T	610 nm
Zr	Eriochrome Cyanine	520 nm

# Table of figures

96

<b>Figure 1</b>	Overlaid Vis spectra of complexed and uncomplexed pyrocatechol violet	16
<b>Figure 2</b>	Vis absorption difference spectrum of complexed and uncomplexed pyrocatechol violet	16
<b>Figure 3</b>	Overlaid Vis spectra of complexed and uncomplexed dithizone	16
<b>Figure 4</b>	Vis absorption difference spectrum of complexed and uncomplexed dithizone	16
<b>Figure 5</b>	Overlaid Vis spectra of complexed and uncomplexed Eriochrome Black T	17
<b>Figure 6</b>	Vis absorption difference spectrum of complexed and uncomplexed Eriochrome Black T	17
<b>Figure 7</b>	Overlaid Vis spectra of complexed and uncomplexed calconcarboxylic acid	17
<b>Figure 8</b>	Vis absorption difference spectrum of complexed and uncomplexed calconcarboxylic acid	17
<b>Figure 9</b>	Overlaid Vis spectra of complexed and uncomplexed hydroxynaphthol blue	18
<b>Figure 10</b>	Vis absorption difference spectrum of complexed and uncomplexed hydroxynaphthol blue	18
<b>Figure 11</b>	Overlaid Vis spectra of complexed and uncomplexed methylthymol blue	18
<b>Figure 12</b>	Vis absorption difference spectrum of complexed and uncomplexed methylthymol blue	18
<b>Figure 13</b>	Overlaid Vis spectra of complexed and uncomplexed murexide	19
<b>Figure 14</b>	Vis absorption difference spectrum of complexed and uncomplexed murexide	19
<b>Figure 15</b>	Overlaid Vis spectra of complexed and uncomplexed PAN	19
<b>Figure 16</b>	Vis absorption difference spectrum of complexed and uncomplexed PAN	19
<b>Figure 17</b>	Overlaid Vis spectra of complexed and uncomplexed phthalein purple	20
<b>Figure 18</b>	Vis absorption difference spectrum of complexed and uncomplexed phthalein purple	20
<b>Figure 19</b>	Overlaid Vis spectra of complexed and uncomplexed Tiron	20
<b>Figure 20</b>	Vis absorption difference spectrum of complexed and uncomplexed Tiron	20
<b>Figure 21</b>	Overlaid Vis spectra of complexed and uncomplexed xylenol orange	21
<b>Figure 22</b>	Vis absorption difference spectrum of complexed and uncomplexed xylenol orange	21
<b>Figure 23</b>	Overlaid Vis spectra of complexed and uncomplexed Eriochrome Cyanine R	21
<b>Figure 24</b>	Vis absorption difference spectrum of complexed and uncomplexed Eriochrome Cyanine R	21
<b>Figure 25</b>	Titration curve of the potentiometric titer determination of EDTA with the Ca-ISE	23
<b>Figure 26</b>	Titration curve of the potentiometric titer determination of EGTA with the Ca-ISE	23
<b>Figure 27</b>	Titration curve of tap water from Herisau (Switzerland), measured with the Ca-ISE	25
<b>Figure 28</b>	Titration curve of the total hardness determination	27
<b>Figure 29</b>	Titration curve of the calcium hardness determination	27



<b>Figure 30</b>	Titration curve of the total hardness determination with the Optrode at 610 nm	29
<b>Figure 31</b>	Titration curve of the calcium hardness determination with the Optrode at 610 nm	29
<b>Figure 32</b>	Titration curve of the potentiometric aluminum determination with the Cu-ISE	31
<b>Figure 33</b>	Titration curve of the potentiometric magnesium determination with the Cu-ISE	31
<b>Figure 34</b>	Titration curve of the potentiometric zinc determination with the Cu-ISE	33
<b>Figure 35</b>	Titration curve of the potentiometric magnesium determination with the Cu-ISE	33
<b>Figure 36</b>	Titration curve of the photometric determination of calcium in cement with the Optrode at 610 nm	37
<b>Figure 37</b>	Titration curve of the photometric determination of magnesium in cement with the Optrode at 610 nm	37
<b>Figure 38</b>	Titration curve of the photometric determination of iron in cement with the Optrode at 610 nm	37
<b>Figure 39</b>	Titration curve of the photometric determination of aluminum in cement with the Optrode at 610 nm	37
<b>Figure 40</b>	Titration curve of the blank value determination with the Ca-ISE	39
<b>Figure 41</b>	Titration curve of the sulfate determination with the Ca-ISE	39
<b>Figure 42</b>	Example curve of a potentiometric aluminum determination	41
<b>Figure 43</b>	Example of a photometric titration curve of aluminum at 610 nm	42
<b>Figure 44</b>	Example of a potentiometric titration curve of barium with the Cu-ISE	44
<b>Figure 45</b>	Photometric titration curve of barium with the Optrode at 574 nm	45
<b>Figure 46</b>	Potentiometric titration curve of bismuth with the Cu-ISE	47
<b>Figure 47</b>	Example of a photometric titration curve of bismuth at 520 nm	48
<b>Figure 48</b>	Example curve of a potentiometric calcium determination with the Ca-ISE	49
<b>Figure 49</b>	Potentiometric titration curve of the determination of calcium in milk with the Cu-ISE	50
<b>Figure 50</b>	Example of a photometric titration curve of the calcium determination with the Optrode at 610 nm	51
<b>Figure 51</b>	Potentiometric titration curve of cadmium with the Cu-ISE	53
<b>Figure 52</b>	Photometric titration curve of the cadmium determination with the Optrode at 610 nm	54
<b>Figure 53</b>	Example of a potentiometric titration curve of cobalt with the Cu-ISE	55
<b>Figure 54</b>	Photometric titration curve of cobalt with the Optrode at 574 nm	56
<b>Figure 55</b>	Potentiometric titration curve of copper with the Cu-ISE	58
<b>Figure 56</b>	Photometric titration curve of copper with the Optrode at 520 nm	59
<b>Figure 57</b>	Example of a potentiometric titration curve of iron with the Cu-ISE	61
<b>Figure 58</b>	Photometric titration curve of iron with the Optrode at 610 nm	62
<b>Figure 59</b>	Potentiometric titration curve of indium with the Cu-ISE	64

<b>Figure 60</b>	Example of a photometric titration curve of gallium with the Optrode at 610 nm	65
<b>Figure 61</b>	Photometric titration curve of mercury with the Optrode at 502 nm	67
<b>Figure 62</b>	Example of a potentiometric titration curve of magnesium with the Cu-ISE	69
<b>Figure 63</b>	Photometric titration curve of magnesium with the Optrode at 610 nm	70
<b>Figure 64</b>	Potentiometric titration curve of manganese with the Cu-ISE	72
<b>Figure 65</b>	Example of a photometric titration curve of manganese with the Optrode at 610 nm	73
<b>Figure 66</b>	Potentiometric titration curve of nickel with the Cu-ISE	75
<b>Figure 67</b>	Example of a photometric titration curve of nickel with the Optrode at 574 nm	76
<b>Figure 68</b>	Potentiometric titration curve of lead with the Cu-ISE	78
<b>Figure 69</b>	Example of a photometric titration curve of lead with the Optrode at 574 nm	79
<b>Figure 70</b>	Photometric titration curve of palladium with the Optrode at 610 nm	80
<b>Figure 71</b>	Example of a photometric titration curve of tin with the Optrode at 574 nm	82
<b>Figure 72</b>	Photometric titration curve of thorium with the Optrode at 574 nm	84
<b>Figure 73</b>	Potentiometric titration curve of thallium with the Cu-ISE	86
<b>Figure 74</b>	Example of a photometric titration curve of thallium with the Optrode at 610 nm	86
<b>Figure 75</b>	Potentiometric titration curve of zinc with the Cu-ISE	88
<b>Figure 76</b>	Example of a photometric titration curve of zinc with the Optrode at 610 nm	89
<b>Figure 77</b>	Example of a potentiometric titration curve of zirconium with the Cu-ISE	91
<b>Figure 78</b>	Photometric titration curve of zirconium with the Optrode at 520 nm	92



[www.metrohm.com](http://www.metrohm.com)

