

## 1 Purpose

This document describes new functions, changes, improvements and fixed bugs in the following firmware version:

- 57.1027.0009 for OS version 008

In comparison with the firmware version:

- 57.1027.0008 for OS version 007

## 2 New functions

The following new functions have been integrated:

- Modifying the dosing rate / filling rate enabled in **Prepare buret (PREP)** menu (TITR-24778).
- Automatic line feed when using the impact printer TX-900MH (TITR-25294).
- Switching off the system is now divided into "Shutting down the system" (standby) and "Switching off the system".
- Test command added to identify the new and old main board.

## 3 Improvements

The following improvements were made:

- The **Stop volume cond.** label has been added to the PDF report.
- In the PDF report, a space has been inserted between the value and the unit for the **start volume**.
- The **Save as CV** label has been removed from the PDF report as it is not required.
- The following values are now also included in the PDF report:
  - Method name
  - Determination time
  - Curve
  - Calculation statistics
  - Measuring point list
- **Monitoring** and **Time interval** have been added to the system report for the thermal printer (if **Monitoring** is activated in the parameters).

