



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

## Method parameters

Method . . . . . AB445 - ASTM D6971 & ASTM D6810 -  
Amine & Phenol  
Method saving date . . . . . 2024-03-27 09:20:45 UTC+1  
Method version . . . . . 1  
Method group . . . . . Main group  
Method status . . . . . original  
Method saved by (full name) . . . . . Metrohm International Headquarters  
Method saved by (short name) . . . . . Metrohm

## START

### Main track

#### General

##### Workplace view

Current view . . . . . on

##### Track view for live window

Live display 1 . . . . . Main track

Live display 2 . . . . . Main track

Electrode check . . . . . off

#### Application note

- Measure 1 g sand into the sample vial. - Place a stirrer bar in the sample vial. - Add 0.2 mL – 0.6 mL standard/sample into the sample vial (volume depends on sample). - Place the sample vial on the rack of the sample changer. In the «Run» window: Select sample type «Standard» to run a calibration or «Sample» to determine the remaining antioxidant in a sample. Under «Sample position» specify the position of standard/sample on the rack and enter the exact volume of standard/sample under «Sample amount».

#### Sample data variables

Name	Type	Assignment	Fixed value	Comment	Monitoring
ID1	Text	ID1		Sample identification 1	off
ID2	Text	ID2		Sample identification 2	off
ID3	Text	ID3		Sample identification 3	off
Sample type	Text	Sample type		Sample type	off
Sample position	Number	Sample position		Sample position number	off
Sample amount	Number	Sample amount		Sample amount	off
Sample amount unit	Text	Sample amount unit	mL	Sample amount unit	off

Name . . . . . ID1

Type . . . . . Text

Assignment . . . . . on. . . . . ID1

Fixed value . . . . . off. . . . .

Check at start . . . . . on



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69

2024-03-27 09:24:49 UTC+1

Comment	Sample identification 1
Name	ID2
Type	Text
Assignment	on, ID2
Fixed value	off.
Check at start	on
Comment	Sample identification 2
Name	ID3
Type	Text
Assignment	on, ID3
Fixed value	off.
Check at start	on
Comment	Sample identification 3
Name	Sample type
Type	Text
Assignment	on, Sample type
Fixed value	off.
Check at start	on
Comment	Sample type
Name	Sample amount unit
Type	Text
Assignment	on, Sample amount unit
Fixed value	on, mL
Check at start	on
Comment	Sample amount unit
Name	Sample position
Type	Number
Assignment	on, Sample position
Fixed value	off.
Check at start	on
Comment	Sample position number
Variable monitoring	off
Lower limit	
Upper limit	
Message	
Display message	on
Record message	on
Message by e-mail	off
E-mail template	
Subject	Message from viva - method 'AOx - BP 2197' - command 'Main track'
Acoustic signal	off



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

Action . . . . . off  
Cancel determination . . . . . on  
Cancel determination and series . . . . . off

Name . . . . . **Sample amount**  
Type . . . . . Number  
Assignment . . . . . on. . . . . Sample amount  
Fixed value . . . . . off. . . . .  
Check at start . . . . . on  
Comment . . . . . Sample amount  
Variable monitoring . . . . . off  
Lower limit . . . . .  
Upper limit . . . . .  
Message . . . . .  
Display message . . . . . on  
Record message . . . . . on  
Message by e-mail . . . . . off  
E-mail template . . . . .  
Subject . . . . . Message from viva - method 'AOx - BP 2197' - command 'Main track'

Acoustic signal . . . . . off  
Action . . . . . off  
Cancel determination . . . . . on  
Cancel determination and series . . . . . off

**Command comment**

ASTM D6971 & ASTM D6810 - Determination of Amine (neutral) & Phenol (alkaline) - fully automated system

**CALL Call prepare sample**

Call text	Track name	Sample type		Condition
Call prepare sample	Prepare sample	off	Sample	off

**CALL Empty measuring vessel**

Call text	Track name	Sample type		Condition
Drain	Drain	off	Sample	off

**CALL Calibration or sample**

Call text	Track name	Sample type		Condition
Call sample	Sample	on	Sample	off
Call calibration	Calibration curve	on	Standard	off

**CALL Rinse needle**



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

## CALL

Call text	Track name	Sample type		Condition
Rinse needle	Go to rinsing position	off	Sample	off

## CALL

### CALL Final rinsing

Call text	Track name	Sample type		Condition
Rinsing cycle	Rinsing cycle	off	Sample	off

## VA TRACK

### VA track neutral

Return immediately . . . . . off

## DP

### Sweep neutral

#### General/Hardware

##### Device

Device name . . . . . 884\_1  
Device type . . . . . 884 Professional VA

##### Sensors/Electrodes

Working electrode . . . . . RDE  
Sensor type . . . . . RDE/SSE  
Reference electrode . . . . . Reference electrode  
Auxiliary electrode . . . . . Auxiliary electrode  
Electrode check . . . . . off

##### Stirrer

Stirring rate . . . . . 2000 min<sup>-1</sup>  
Hydrodynamic measurement . . . . . off

#### Pretreatment

Stirring time . . . . . 5.0 s

##### Cyclovoltammetric pretreatment

Start potential . . . . . 0 V  
Vertex potential . . . . . 1 V  
Sweep rate . . . . . 1 V/s  
Cycles . . . . . 1  
Duration . . . . . 2.00 s

##### Potentiostatic pretreatment

Potential 1 . . . . . off V  
Waiting time 1 . . . . . 0.0 s  
Potential 2 . . . . . off V  
Waiting time 2 . . . . . 0.0 s  
Potential 3 . . . . . off V  
Waiting time 3 . . . . . 0.0 s  
Potential 4 . . . . . off V  
Waiting time 4 . . . . . 0.0 s  
Potential 5 . . . . . off V  
Waiting time 5 . . . . . 0.0 s



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

Equilibration time . . . . . 5.0 s

**Sweep**

Start potential . . . . . 0 V  
End potential . . . . . 1.3 V  
Potential step . . . . . 0.012 V  
Potential step time . . . . . 0.1 s  
Sweep rate . . . . . 0.120 V/s  
Pulse amplitude . . . . . 0.05 V  
Pulse time . . . . . 0.04 s  
Measuring time . . . . . 0.02 s  
Sweep duration . . . . . 10.83 s

**Post-treatment**

**Cleaning**

Cleaning potential . . . . . off V  
Cleaning time . . . . . 0.0 s

**Standby potential**

Standby potential . . . . . off V

**Potentiostat**

**Current measuring range**

Highest range . . . . . 2 mA  
Lowest range . . . . . 20  $\mu$ A  
Automatically select optimum current measuring range . . . . . on

**VA TRACK VA track alkaline**

Return immediately . . . . . off

**DP**

**Sweep alkaline**

**General/Hardware**

**Device**

Device name . . . . . 884\_1  
Device type . . . . . 884 Professional VA

**Sensors/Electrodes**

Working electrode . . . . . RDE  
Sensor type . . . . . RDE/SSE  
Reference electrode . . . . . Reference electrode  
Auxiliary electrode . . . . . Auxiliary electrode  
Electrode check . . . . . off

**Stirrer**

Stirring rate . . . . . 2000 min<sup>-1</sup>  
Hydrodynamic measurement . . . . . off

**Pretreatment**

Stirring time . . . . . 5.0 s

**Cyclovoltammetric pretreatment**

Start potential . . . . . 0 V  
Vertex potential . . . . . 1 V  
Sweep rate . . . . . 1 V/s



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

Cycles . . . . . 1  
Duration . . . . . 2.00 s  
Potentiostatic pretreatment  
Potential 1 . . . . . off V  
Waiting time 1 . . . . . 0.0 s  
Potential 2 . . . . . off V  
Waiting time 2 . . . . . 0.0 s  
Potential 3 . . . . . off V  
Waiting time 3 . . . . . 0.0 s  
Potential 4 . . . . . off V  
Waiting time 4 . . . . . 0.0 s  
Potential 5 . . . . . off V  
Waiting time 5 . . . . . 0.0 s  
Equilibration time . . . . . 5.0 s  
**Sweep**  
Start potential . . . . . -0.1 V  
End potential . . . . . 0.7 V  
Potential step . . . . . 0.012 V  
Potential step time . . . . . 0.1 s  
Sweep rate . . . . . 0.120 V/s  
Pulse amplitude . . . . . 0.05 V  
Pulse time . . . . . 0.04 s  
Measuring time . . . . . 0.02 s  
Sweep duration . . . . . 6.67 s  
**Post-treatment**  
Cleaning  
Cleaning potential . . . . . off V  
Cleaning time . . . . . 0.0 s  
Standby potential  
Standby potential . . . . . off V  
**Potentiostat**  
Current measuring range  
Highest range . . . . . 2 mA  
Lowest range . . . . . 20 µA  
Automatically select optimum current measuring range . . . . . on

**TRACK**     **Sample**  
Return immediately . . . . . off  
Delete old data . . . . . off

**CALL**     **Call dose electrolyte sam**

Call text	Track name	Sample type	Condition
Dose electrolyte	Electrolyte blank off	Sample	off

**Call blank sample**



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

**CALL BLANK** Call text . . . . . Call blank sample  
Track name . . . . . VA track neutral  
  
Condition . . . . . off

**CALL Empty vessel sample**

Call text	Track name	Sample type	Condition
Empty vessel sample	Drain	off Sample	off

**CALL Transfer sample**

Call text	Track name	Sample type	Condition
Transfer sample	Transfer from autosampler	off Sample	off

**ADD SAMPLE Add sample**

Addition  
Add manually . . . . . off  
Already added . . . . . on  
Add with dosing device . . . . . off

**STIR & PURGE Stir sample**

Device  
Device name . . . . . 884\_1  
Device type . . . . . 884 Professional VA  
  
Stir  
Stirring rate . . . . . 2000 min<sup>-1</sup>  
Switch on . . . . . off  
Switch off . . . . . off  
Duration . . . . . on  
Time . . . . . 5 s  
  
Purge  
Switch on . . . . . off  
Switch off . . . . . on  
Duration . . . . . off

**CALL COND Call cond sample**

Call text . . . . . Call cond sample  
Track name . . . . . VA track neutral  
Stop criteria  
Evaluation quantity . . . . . off  
Number of runs . . . . . on  
    Number of runs . . . . . 1  
    Action . . . . . Cancel command  
Run time . . . . . off  
Condition . . . . . off



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

**LOOP Replications sample**  
Stop criteria  
Max. run number . . . . . on  
Max. run number . . . . . 2  
Maximum run time . . . . . off  
Signal assessment for DT . . . . . off  
Condition . . . . . off

**CALL VA Sample neutral**  
Call text . . . . . Measure sample neutral  
Track name . . . . . VA track neutral  
Condition . . . . . off

**CALL Add alkaline sample**

Call text	Track name	Sample type	Condition
Add alkaline solution	Alkaline solution	off Sample	off

**CALL COND Cond sample alkaline**  
Call text . . . . . Call cond sample alkaline  
Track name . . . . . VA track alkaline  
Stop criteria  
Evaluation quantity . . . . . off  
Number of runs . . . . . on  
Number of runs . . . . . 1  
Action . . . . . Cancel command  
Run time . . . . . off  
Condition . . . . . off

**LOOP Replications sample alkal**  
Stop criteria  
Max. run number . . . . . on  
Max. run number . . . . . 2  
Maximum run time . . . . . off  
Signal assessment for DT . . . . . off  
Condition . . . . . off

**CALL VA Sample alkaline**  
Call text . . . . . Measure sample alkaline  
Track name . . . . . VA track alkaline  
Condition . . . . . off

**TRACK Calibration curve**  
Return immediately . . . . . off  
Delete old data . . . . . off

**CALL Call dose electrolyte std**





License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

Call text	Track name	Sample type	Condition
Dose electrolyte	Electrolyte blank	off Sample	off

**CALL  
BLANK**

**Call blank calibration**

Call text . . . . . Call blank calibration  
Track name . . . . . VA track neutral  
Condition . . . . . off

**CALL**

**Empty vessel standard**

Call text	Track name	Sample type	Condition
Empty vessel standard Drain		off Sample	off

**CALL**

**Transfer standard**

Call text	Track name	Sample type	Condition
Transfer sample	Transfer from autosampler	off Sample	off

**ADD STD**

**Add standard**

Standard

Solution . . . . . Standard

Addition increments

Number . . . . . 1

Addition volume 1 . . . . . = 'SD.Sample amount' \* 1000 mL

Addition

Add manually . . . . . off

Already added . . . . . on

Add with dosing device . . . . . off

**STIR &  
PURGE**

**Stir calibration**

Device

Device name . . . . . 884\_1

Device type . . . . . 884 Professional VA

Stir

Stirring rate . . . . . 2000 min<sup>-1</sup>

Switch on . . . . . off

Switch off . . . . . off

Duration . . . . . on

Time . . . . . 5 s

Purge

Switch on . . . . . off

Switch off . . . . . on

Duration . . . . . off



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

**CALL  
COND**

**Call cond calibration**

Call text . . . . . Call cond calibration  
Track name . . . . . VA track neutral  
Stop criteria  
Evaluation quantity . . . . . off  
Number of runs . . . . . on  
Number of runs . . . . . 1  
Action . . . . . Cancel command  
Run time . . . . . off  
Condition . . . . . off

**LOOP**

**Replications calibration**

Stop criteria  
Max. run number . . . . . on  
Max. run number . . . . . 2  
Maximum run time . . . . . off  
Signal assessment for DT . . . . . off  
Condition . . . . . off

**CALL VA**

**Calibration neutral**

Call text . . . . . Measure calibration neutral  
Track name . . . . . VA track neutral  
Condition . . . . . off

**CALL**

**Add alkaline calibration**

Call text	Track name	Sample type	Condition
Add alkaline solution	Alkaline solution	off Sample	off

**CALL  
COND**

**Cond calibration alkaline**

Call text . . . . . Call cond calibration alkaline  
Track name . . . . . VA track alkaline  
Stop criteria  
Evaluation quantity . . . . . off  
Number of runs . . . . . on  
Number of runs . . . . . 1  
Action . . . . . Cancel command  
Run time . . . . . off  
Condition . . . . . off

**LOOP**

**Replications calibr alkal**

Stop criteria  
Max. run number . . . . . on  
Max. run number . . . . . 2  
Maximum run time . . . . . off  
Signal assessment for DT . . . . . off  
Condition . . . . . off



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

**CALL VA Calibration alkaline**

Call text . . . . . Measure calibration alkaline  
Track name . . . . . VA track alkaline  
Condition . . . . . off

**TRACK Electrolyte blank**

Return immediately . . . . . off  
Delete old data . . . . . off

**LQH Dose electrolyte blank**

**General/Hardware**

**Dosing device**

Dosing unit . . . . . Electrolyte neutral

**Parameters**

Function . . . . . Dosing  
Port . . . . . 3  
Volume . . . . . 5 mL  
Rate . . . . . maximum mL/min

**STIR & PURGE Stir electrolyte**

**Device**

Device name . . . . . 884\_1  
Device type . . . . . 884 Professional VA

**Stir**

Stirring rate . . . . . 2000 min<sup>-1</sup>  
Switch on . . . . . off  
Switch off . . . . . off  
Duration . . . . . on  
Time . . . . . 10 s

**Purge**

Switch on . . . . . off  
Switch off . . . . . on  
Duration . . . . . off

**TRACK Alkaline solution**

Return immediately . . . . . off  
Delete old data . . . . . off

**STIR & PURGE Stir ON**

**Device**

Device name . . . . . 884\_1  
Device type . . . . . 884 Professional VA

**Stir**

Stirring rate . . . . . 2000 min<sup>-1</sup>  
Switch on . . . . . on  
Switch off . . . . . off  
Duration . . . . . off

**Purge**



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

Switch on . . . . . off  
Switch off . . . . . on  
Duration . . . . . off

**LQH Add alkaline solution**  
**General/Hardware**

Dosing device  
Dosing unit . . . . . Electrolyte alkaline

**Parameters**  
Function . . . . . Dosing  
Port . . . . . 3  
Volume . . . . . 2 mL  
Rate . . . . . maximum mL/min

**ADD AUX Volume alkaline solution**

Auxiliary solution  
Solution . . . . . Electrolyte alkaline  
Volume . . . . . = 'Add alkaline solution.VOL' \* 1000 mL  
Include volume in calculation . . . . . on  
Addition  
Add manually . . . . . off  
Already added . . . . . on  
Add with dosing device . . . . . off

**STIR & PURGE Stir alkaline**

Device  
Device name . . . . . 884\_1  
Device type . . . . . 884 Professional VA  
Stir  
Stirring rate . . . . . 2000 min<sup>-1</sup>  
Switch on . . . . . off  
Switch off . . . . . off  
Duration . . . . . on  
Time . . . . . 30 s  
Purge  
Switch on . . . . . off  
Switch off . . . . . on  
Duration . . . . . off

**TRACK Prepare sample**

Return immediately . . . . . off  
Delete old data . . . . . off

**SWING Swing waste pos**

Device  
Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor  
Target



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

Tower . . . . . 1  
Swing . . . . . External position  
Number . . . . . 2  
Parameters  
Swing rate . . . . . 55 °/s

**LIFT**

**Lift waste pos**

Device  
Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor  
Target  
Tower . . . . . 1  
Lift position . . . . . Work position mm  
Parameters  
Lift rate . . . . . 25 mm/s

**ADD AUX**

**Fill needle electrolyte**

Auxiliary solution  
Solution . . . . . Electrolyte neutral  
Volume . . . . . 1 mL  
Include volume in calculation . . . . . off  
Addition  
Add manually . . . . . off  
Already added . . . . . off  
Add with dosing device . . . . . on  
Dosing rate . . . . . 10 mL/min  
Filling rate . . . . . maximum mL/min

**MOVE**

**MOVE to sample**

Device  
Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor  
Target  
Tower . . . . . 1  
Move . . . . . Sample position  
Beaker test  
Display message . . . . . on  
Cancel determination . . . . . off  
Cancel determination and series . . . . . off  
Parameters  
Shift rate . . . . . 20 °/s  
Shift direction . . . . . auto  
Swing rate . . . . . 55 °/s

**LIFT**

**Needle → special position**

Device  
Device name . . . . . 858 Sample Processor



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

Device type . . . . . 858.0020 Professional Sample Processor  
Target  
Tower . . . . . 1  
Lift position . . . . . Special position mm  
Parameters  
Lift rate . . . . . 25 mm/s

**STIR**

**Stir sample ON**

Device  
Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor  
Stirrer  
Stirrer . . . . . 1  
Stirrer type . . . . . unknown  
Stirring rate . . . . . 11  
Action  
Switch on . . . . . on  
Switch off . . . . . off  
Duration . . . . . off  
Time . . . . . 1 min

**ADD AUX**

**Add electrolyte**

Auxiliary solution  
Solution . . . . . Electrolyte neutral  
Volume . . . . . 10 mL  
Include volume in calculation . . . . . on  
Addition  
Add manually . . . . . off  
Already added . . . . . off  
Add with dosing device . . . . . on  
Dosing rate . . . . . 10 mL/min  
Filling rate . . . . . maximum mL/min

**STIR**

**Stir sample OFF**

Device  
Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor  
Stirrer  
Stirrer . . . . . 1  
Stirrer type . . . . . unknown  
Stirring rate . . . . . 11  
Action  
Switch on . . . . . off  
Switch off . . . . . off  
Duration . . . . . on  
Time . . . . . 20 s



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

**TRACK Transfer from autosampler**

Return immediately . . . . . off  
Delete old data . . . . . off

**LIFT Needle → Work position**

**Device**

Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor

**Target**

Tower . . . . . 1  
Lift position . . . . . Work position mm

**Parameters**

Lift rate . . . . . 25 mm/s

**CALL Transfer solution**

Call text	Track name	Sample type	Condition
Peristaltic pump ON	Peristaltic pump	off Sample	off

**TRACK Peristaltic pump**

Return immediately . . . . . off  
Delete old data . . . . . off

**PUMP Peristaltic pump ON**

**Device**

Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor

**Pumps**

Tower . . . . . 1  
Pump(s) . . . . . Peristaltic  
Rate . . . . . 5

**Action**

Switch on . . . . . off  
Switch off . . . . . off  
Duration . . . . . on  
Time . . . . . 150 s

**TRACK Go to rinsing position**

Return immediately . . . . . off  
Delete old data . . . . . off

**CALL Empty vessel**

Call text	Track name	Sample type	Condition
Drain	Drain	off Sample	off

**Swing rinse needle**



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

## SWING

### Device

Device name . . . . . 858 Sample Processor

Device type . . . . . 858.0020 Professional Sample Processor

### Target

Tower . . . . . 1

Swing . . . . . External position

Number . . . . . 1

### Parameters

Swing rate . . . . . 55 °/s

## LIFT

### Needle → Rinse position

### Device

Device name . . . . . 858 Sample Processor

Device type . . . . . 858.0020 Professional Sample Processor

### Target

Tower . . . . . 1

Lift position . . . . . Work position mm

### Parameters

Lift rate . . . . . 25 mm/s

## LQH

### Rinse needle

### General/Hardware

### Dosing device

Dosing unit . . . . . WashStation

### Parameters

Function . . . . . Dosing

Port . . . . . 1

Volume . . . . . 7 mL

Rate . . . . . maximum mL/min

## CALL

### CALL WashStation

Call text	Track name	Sample type	Condition
Refill WashStation	Refill WashStation	off Sample	off

## STIR & PURGE

### Stirrer ON

### Device

Device name . . . . . 884\_1

Device type . . . . . 884 Professional VA

### Stir

Stirring rate . . . . . 3000 min<sup>-1</sup>

Switch on . . . . . on

Switch off . . . . . off

Duration . . . . . off

### Purge





License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

Switch on . . . . . off  
Switch off . . . . . on  
Duration . . . . . off

**CALL Transfer rinsing solution**

Call text	Track name	Sample type	Condition
Peristaltic pump ON	Peristaltic pump off	Sample	off

**LIFT Lift needle**

Device  
Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor  
Target  
Tower . . . . . 1  
Lift position . . . . . Home position mm  
Parameters  
Lift rate . . . . . 25 mm/s

**PUMP Dose rinse EtOH**

Device  
Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor  
Pumps  
Tower . . . . . 1  
Pump(s) . . . . . 2  
Action  
Switch on . . . . . off  
Switch off . . . . . off  
Duration . . . . . on  
Time . . . . . 3 s

**PUMP Empty tube**

Device  
Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor  
Pumps  
Tower . . . . . 1  
Pump(s) . . . . . Peristaltic  
Rate . . . . . 5  
Action  
Switch on . . . . . off  
Switch off . . . . . off  
Duration . . . . . on  
Time . . . . . 30 s

**Refill WashStation**



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

**TRACK** Return immediately . . . . . on

Delete old data . . . . . off

**WAIT WAIT 40 s**

Wait

Stop track and waiting for [Continue] . . . . . off

Stop all tracks and waiting for [Continue] . . . . . off

Waiting time . . . . . on

Time . . . . . 40

Unit . . . . . s

Message

Record message . . . . . off

Message by e-mail . . . . . off

Acoustic signal . . . . . off

**LQH Refill**

**General/Hardware**

Dosing device

Dosing unit . . . . . WashStation

**Parameters**

Function . . . . . Dosing

Port . . . . . 1

Volume . . . . . 3 mL

Rate . . . . . 4 mL/min

**TRACK Rinsing cycle**

Return immediately . . . . . off

Delete old data . . . . . off

**LOOP LOOP Rinsing**

Stop criteria

Max. run number . . . . . on

Max. run number . . . . . 2

Maximum run time . . . . . off

Signal assessment for DT . . . . . off

Condition . . . . . off

**CALL CALL Drain**

Call text	Track name	Sample type	Condition
Drain	Drain	off Sample	off

**CALL CALL Rinse**



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

Call text	Track name	Sample type	Condition
Rinse	Rinse	off Sample	off

**TRACK****Drain**

Return immediately . . . . . off  
Delete old data . . . . . off

**PUMP****PUMP Drain****Device**

Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor

**Pumps**

Tower . . . . . 1  
Pump(s) . . . . . 1

**Action**

Switch on . . . . . off  
Switch off . . . . . off  
Duration . . . . . on  
Time . . . . . 15 s

**TRACK****Rinse**

Return immediately . . . . . off  
Delete old data . . . . . off

**PUMP****PUMP Rinse****Device**

Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor

**Pumps**

Tower . . . . . 1  
Pump(s) . . . . . 2

**Action**

Switch on . . . . . off  
Switch off . . . . . off  
Duration . . . . . on  
Time . . . . . 5 s

**TRACK****Fill DU neutral**

Return immediately . . . . . off  
Delete old data . . . . . off

**LQH****DU Electrolyte neutral****General/Hardware****Dosing device**

Dosing unit . . . . . Electrolyte neutral

**Parameters**



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

Function . . . . . Fill  
Port . . . . . 2  
Rate . . . . . maximum mL/min

**TRACK Fill DU alkaline**  
Return immediately . . . . . off  
Delete old data . . . . . off

**LQH DU Electrolyte alkaline**  
**General/Hardware**  
Dosing device  
Dosing unit . . . . . Electrolyte alkaline  
**Parameters**  
Function . . . . . Fill  
Port . . . . . 2  
Rate . . . . . maximum mL/min

**TRACK Fill DU WashStation**  
Return immediately . . . . . on  
Delete old data . . . . . off

**LQH DU WashStation**  
**General/Hardware**  
Dosing device  
Dosing unit . . . . . WashStation  
**Parameters**  
Function . . . . . Fill  
Port . . . . . 2  
Rate . . . . . maximum mL/min

**TRACK Shut off**  
Return immediately . . . . . off  
Delete old data . . . . . off

**PUMP PUMPS Rinse OFF**  
Device  
Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor  
Pumps  
Tower . . . . . 1  
Pump(s) . . . . . 1+2  
Action  
Switch on . . . . . off  
Switch off . . . . . on  
Duration . . . . . off  
Time . . . . . 8.0 s

**PUMP PUMP Peri OFF**  
Device



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:24:49 UTC+1

Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor  
Pumps  
Tower . . . . . 1  
Pump(s) . . . . . Peristaltic  
Rate . . . . . 10  
Action  
Switch on . . . . . off  
Switch off . . . . . on  
Duration . . . . . off  
Time . . . . . 300.0 s

**STIR**

**STIR OFF**

Device  
Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor  
Stirrer  
Stirrer . . . . . 1  
Stirrer type . . . . . unknown  
Stirring rate . . . . . 11  
Action  
Switch on . . . . . off  
Switch off . . . . . on  
Duration . . . . . off  
Time . . . . . 10.0 s

**STIR &  
PURGE**

**STIR & PURGE OFF**

Device  
Device name . . . . . 884\_1  
Device type . . . . . 884 Professional VA  
Stir  
Stirring rate . . . . . 2000 min<sup>-1</sup>  
Switch on . . . . . off  
Switch off . . . . . on  
Duration . . . . . off  
Purge  
Switch on . . . . . off  
Switch off . . . . . on  
Duration . . . . . off

**EXIT**

**Exit track**

**CALL**

**Exit shut off**



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69

2024-03-27 09:24:49 UTC+1

Call text	Track name	Sample type	Condition
Exit shut off	Shut off	off Sample	off

#### CALL

##### Call fill DU

Call text	Track name	Sample type	Condition
Fill DU neutral	Fill DU neutral	off Sample	off
Fill DU alkaline	Fill DU alkaline	off Sample	off
Fill DU WashStation	Fill DU WashStation	off Sample	off

#### LIFT

##### Needle → Home position

###### Device

Device name . . . . . 858 Sample Processor

Device type . . . . . 858.0020 Professional Sample Processor

###### Target

Tower . . . . . 1

Lift position . . . . . Home position mm

###### Parameters

Lift rate . . . . . 25 mm/s

#### ERROR

##### Error track

#### CALL

##### Error shut off

Call text	Track name	Sample type	Condition
Error shut off	Shut off	off Sample	off