



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

Method parameters

Method AB445 – ASTM D6971 & D6810 - Amine & Phenol
Method saving date 2024-06-26 13:42:58 UTC+2
Method version 1
Method group Main group
Method status original
Method saved by (full name) Metrohm International Headquarters
Method saved by (short name) Metrohm

START

Main track

General

Workplace view

Current view on

Track view for live window

Live display 1 Main track

Live display 2 Main track

Electrode check off

Application note

- Measure 1 g sand into the sample vial. - Place a stirrer bar in the sample vial. - Add 0.05 mL – 0.6 mL reference/in-service oil into the sample vial (volume depends on sample). - Place the sample vial on the rack of the sample changer. In the «Run» window: Select sample type «Standard» to run a calibration with reference oil, or «Sample» to determine the remaining antioxidant in an in-service oil sample. Under «Sample position» specify the position of reference/in-service oil on the rack and enter the exact volume of reference/in-service oil under «Sample amount».

Sample data variables

Name	Type	Assignment	Fixed value	Comment	Monitoring
ID1	Text	ID1		Sample identification 1	off
ID2	Text	ID2		Sample identification 2	off
ID3	Text	ID3		Sample identification 3	off
Sample type	Text	Sample type		Sample type	off
Sample position	Number	Sample position		Sample position number	off
Sample amount	Number	Sample amount		Sample amount	off
Sample amount unit	Text	Sample amount unit	mL	Sample amount unit	off

Name ID1
Type Text
Assignment on ID1
Fixed value off
Check at start on



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

Comment	Sample identification 1
<hr/>	
Name	ID2
Type	Text
Assignment	on, ID2
Fixed value	off.
Check at start	on
Comment	Sample identification 2
<hr/>	
Name	ID3
Type	Text
Assignment	on, ID3
Fixed value	off.
Check at start	on
Comment	Sample identification 3
<hr/>	
Name	Sample type
Type	Text
Assignment	on, Sample type
Fixed value	off.
Check at start	on
Comment	Sample type
<hr/>	
Name	Sample amount unit
Type	Text
Assignment	on, Sample amount unit
Fixed value	on, mL
Check at start	on
Comment	Sample amount unit
<hr/>	
Name	Sample position
Type	Number
Assignment	on, Sample position
Fixed value	off.
Check at start	on
Comment	Sample position number
Variable monitoring	off
Lower limit	
Upper limit	
Message	
Display message	on
Record message	on
Message by e-mail	off
E-mail template	
Subject	Message from viva - method 'AOx - BP 2197' - command 'Main track'
Acoustic signal	off



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

Action off
Cancel determination on
Cancel determination and series off

Name **Sample amount**
Type Number
Assignment on. Sample amount
Fixed value off.
Check at start on
Comment Sample amount
Variable monitoring off
Lower limit
Upper limit
Message
Display message on
Record message on
Message by e-mail off
E-mail template
Subject Message from viva - method 'AOx - BP 2197' - command 'Main track'

Acoustic signal off
Action off
Cancel determination on
Cancel determination and series off

Command comment

AB-445/1: MVA-25 - ASTM D6971 & ASTM D6810 - Determination of Amine (neutral) & Phenol (alkaline)

CALL Call prepare sample

Call text	Track name	Sample type		Condition
Call prepare sample	Prepare sample	off	Sample	off

CALL Empty measuring vessel

Call text	Track name	Sample type		Condition
Drain	Drain	off	Sample	off

CALL Calibration or sample

Call text	Track name	Sample type		Condition
Call sample	Sample	on	Sample	off
Call calibration	Calibration curve	on	Standard	off

CALL Rinse needle



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

CALL

Call text	Track name	Sample type		Condition
Rinse needle	Go to rinsing position	off	Sample	off

CALL

CALL Final rinsing

Call text	Track name	Sample type		Condition
Rinsing cycle	Rinsing cycle	off	Sample	off

VA TRACK

VA track neutral

Return immediately off

DP

Sweep neutral

General/Hardware

Device

Device name 884_1

Device type 884 Professional VA

Sensors/Electrodes

Working electrode RDE

Sensor type RDE/SSE

Reference electrode Reference electrode

Auxiliary electrode Auxiliary electrode

Electrode check off

Stirrer

Stirring rate 2000 min⁻¹

Hydrodynamic measurement off

Pretreatment

Stirring time 5.0 s

Cyclovoltammetric pretreatment

Start potential 0 V

Vertex potential 1 V

Sweep rate 1 V/s

Cycles 1

Duration 2.00 s

Potentiostatic pretreatment

Potential 1 off V

Waiting time 1 0.0 s

Potential 2 off V

Waiting time 2 0.0 s

Potential 3 off V

Waiting time 3 0.0 s

Potential 4 off V

Waiting time 4 0.0 s

Potential 5 off V

Waiting time 5 0.0 s



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

Equilibration time 5.0 s

Sweep

Start potential 0 V
End potential 1.3 V
Potential step 0.012 V
Potential step time 0.1 s
Sweep rate 0.120 V/s
Pulse amplitude 0.05 V
Pulse time 0.04 s
Measuring time 0.02 s
Sweep duration 10.83 s

Post-treatment

Cleaning

Cleaning potential off V
Cleaning time 0.0 s

Standby potential

Standby potential off V

Potentiostat

Current measuring range

Highest range 2 mA
Lowest range 20 μ A
Automatically select optimum current measuring range on

VA TRACK VA track alkaline

Return immediately off

DP

Sweep alkaline

General/Hardware

Device

Device name 884_1
Device type 884 Professional VA

Sensors/Electrodes

Working electrode RDE
Sensor type RDE/SSE
Reference electrode Reference electrode
Auxiliary electrode Auxiliary electrode
Electrode check off

Stirrer

Stirring rate 2000 min⁻¹
Hydrodynamic measurement off

Pretreatment

Stirring time 5.0 s

Cyclovoltammetric pretreatment

Start potential 0 V
Vertex potential 1 V
Sweep rate 1 V/s



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

Cycles 1
Duration 2.00 s
Potentiostatic pretreatment
Potential 1 off V
Waiting time 1 0.0 s
Potential 2 off V
Waiting time 2 0.0 s
Potential 3 off V
Waiting time 3 0.0 s
Potential 4 off V
Waiting time 4 0.0 s
Potential 5 off V
Waiting time 5 0.0 s
Equilibration time 5.0 s
Sweep
Start potential -0.1 V
End potential 0.7 V
Potential step 0.012 V
Potential step time 0.1 s
Sweep rate 0.120 V/s
Pulse amplitude 0.05 V
Pulse time 0.04 s
Measuring time 0.02 s
Sweep duration 6.67 s
Post-treatment
Cleaning
Cleaning potential off V
Cleaning time 0.0 s
Standby potential
Standby potential off V
Potentiostat
Current measuring range
Highest range 2 mA
Lowest range 20 µA
Automatically select optimum current measuring range on

TRACK **Sample**
Return immediately off
Delete old data off

CALL **Call dose electrolyte sam**

Call text	Track name	Sample type	Condition
Dose electrolyte	Electrolyte blank off	Sample	off

Call blank sample



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

CALL BLANK Call text Call blank sample
Track name VA track neutral

Condition off

CALL Empty vessel sample

Call text	Track name	Sample type	Condition
Empty vessel sample	Drain	off Sample	off

CALL Transfer sample

Call text	Track name	Sample type	Condition
Transfer sample	Transfer from autosampler	off Sample	off

ADD SAMPLE Add sample

Addition
Add manually off
Already added on
Add with dosing device off

STIR & PURGE Stir sample

Device
Device name 884_1
Device type 884 Professional VA

Stir
Stirring rate 2000 min⁻¹
Switch on off
Switch off off
Duration on
Time 5 s

Purge
Switch on off
Switch off on
Duration off

CALL COND Call cond sample

Call text Call cond sample
Track name VA track neutral
Stop criteria
Evaluation quantity off
Number of runs on
 Number of runs 1
 Action Cancel command
Run time off
Condition off



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

LOOP Replications sample
Stop criteria
Max. run number on
Max. run number 2
Maximum run time off
Signal assessment for DT off
Condition off

CALL VA Sample neutral
Call text Measure sample neutral
Track name VA track neutral
Condition off

CALL Add alkaline sample

Call text	Track name	Sample type	Condition
Add alkaline solution	Alkaline solution	off Sample	off

CALL COND Cond sample alkaline
Call text Call cond sample alkaline
Track name VA track alkaline
Stop criteria
Evaluation quantity off
Number of runs on
Number of runs 1
Action Cancel command
Run time off
Condition off

LOOP Replications sample alkal
Stop criteria
Max. run number on
Max. run number 2
Maximum run time off
Signal assessment for DT off
Condition off

CALL VA Sample alkaline
Call text Measure sample alkaline
Track name VA track alkaline
Condition off

TRACK Calibration curve
Return immediately off
Delete old data off

CALL Call dose electrolyte std



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

Call text	Track name	Sample type	Condition
Dose electrolyte	Electrolyte blank	off Sample	off

**CALL
BLANK**

Call blank calibration

Call text Call blank calibration
Track name VA track neutral
Condition off

CALL

Empty vessel standard

Call text	Track name	Sample type	Condition
Empty vessel standard Drain		off Sample	off

CALL

Transfer standard

Call text	Track name	Sample type	Condition
Transfer sample	Transfer from autosampler	off Sample	off

ADD STD

Add standard

Standard
Solution Standard
Addition increments
Number 1
Addition volume 1 = 'SD.Sample amount' * 1000 mL
Addition
Add manually off
Already added on
Add with dosing device off

**STIR &
PURGE**

Stir calibration

Device
Device name 884_1
Device type 884 Professional VA
Stir
Stirring rate 2000 min⁻¹
Switch on off
Switch off off
Duration on
Time 5 s
Purge
Switch on off
Switch off on
Duration off



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

**CALL
COND**

Call cond calibration

Call text Call cond calibration
Track name VA track neutral
Stop criteria
Evaluation quantity off
Number of runs on
Number of runs 1
Action Cancel command
Run time off
Condition off

LOOP

Replications calibration

Stop criteria
Max. run number on
Max. run number 2
Maximum run time off
Signal assessment for DT off
Condition off

CALL VA

Calibration neutral

Call text Measure calibration neutral
Track name VA track neutral
Condition off

CALL

Add alkaline calibration

Call text	Track name	Sample type	Condition
Add alkaline solution	Alkaline solution	off Sample	off

**CALL
COND**

Cond calibration alkaline

Call text Call cond calibration alkaline
Track name VA track alkaline
Stop criteria
Evaluation quantity off
Number of runs on
Number of runs 1
Action Cancel command
Run time off
Condition off

LOOP

Replications calibr alkal

Stop criteria
Max. run number on
Max. run number 2
Maximum run time off
Signal assessment for DT off
Condition off



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

CALL VA Calibration alkaline

Call text Measure calibration alkaline
Track name VA track alkaline
Condition off

TRACK Electrolyte blank

Return immediately off
Delete old data off

LQH Dose blank neutral

General/Hardware

Dosing device

Dosing unit Electrolyte neutral

Parameters

Function Dosing
Port 3
Volume 5 mL
Rate maximum mL/min

STIR & PURGE Stir electrolyte

Device

Device name 884_1
Device type 884 Professional VA

Stir

Stirring rate 2000 min⁻¹
Switch on off
Switch off off
Duration on
Time 10 s

Purge

Switch on off
Switch off on
Duration off

TRACK Prepare sample

Return immediately off
Delete old data off

SWING Swing waste pos

Device

Device name 858 Sample Processor
Device type 858.0020 Professional Sample Processor

Target

Tower 1
Swing External position
Number 2

Parameters

Swing rate 55 °/s



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

LIFT

Lift waste pos

Device

Device name 858 Sample Processor
Device type 858.0020 Professional Sample Processor

Target

Tower 1
Lift position Work position mm

Parameters

Lift rate 25 mm/s

LQH

Fill needle neutral

General/Hardware

Dosing device

Dosing unit Electrolyte neutral

Parameters

Function Dosing
Port 1
Volume 1 mL
Rate 10 mL/min

MOVE

MOVE to sample

Device

Device name 858 Sample Processor
Device type 858.0020 Professional Sample Processor

Target

Tower 1
Move Sample position

Beaker test

Display message on
Cancel determination off
Cancel determination and series off

Parameters

Shift rate 20 °/s
Shift direction auto
Swing rate 55 °/s

LIFT

Needle → special position

Device

Device name 858 Sample Processor
Device type 858.0020 Professional Sample Processor

Target

Tower 1
Lift position Special position mm

Parameters

Lift rate 25 mm/s

Stir sample ON



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

STIR**Device**

Device name 858 Sample Processor

Device type 858.0020 Professional Sample Processor

Stirrer

Stirrer 1

Stirrer type unknown

Stirring rate 11

Action

Switch on on

Switch off off

Duration off

Time 1 min

ADD AUX**Add electrolyte neutral****Auxiliary solution**

Solution Electrolyte neutral

Volume 10 mL

Include volume in calculation on

Addition

Add manually off

Already added off

Add with dosing device on

Dosing rate 10 mL/min

Filling rate maximum mL/min

STIR**Stir sample OFF****Device**

Device name 858 Sample Processor

Device type 858.0020 Professional Sample Processor

Stirrer

Stirrer 1

Stirrer type unknown

Stirring rate 11

Action

Switch on off

Switch off off

Duration on

Time 20 s

TRACK**Alkaline solution**

Return immediately off

Delete old data off

**STIR &
PURGE****Stir ON****Device**

Device name 884_1



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

Device type 884 Professional VA

Stir

Stirring rate 2000 min⁻¹
Switch on on
Switch off off
Duration off

Purge

Switch on off
Switch off on
Duration off

LQH

Add alkaline solution

General/Hardware

Dosing device

Dosing unit Electrolyte alkaline

Parameters

Function Dosing
Port 3
Volume 2 mL
Rate maximum mL/min

ADD AUX

Volume alkaline solution

Auxiliary solution

Solution Electrolyte alkaline
Volume = 'Add alkaline solution.VOL' * 1000 mL
Include volume in calculation on

Addition

Add manually off
Already added on
Add with dosing device off

**STIR &
PURGE**

Stir alkaline

Device

Device name 884_1
Device type 884 Professional VA

Stir

Stirring rate 2000 min⁻¹
Switch on off
Switch off off
Duration on
Time 30 s

Purge

Switch on off
Switch off on
Duration off

Transfer from autosampler



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

TRACK Return immediately off

Delete old data off

LIFT Needle → Work position

Device

Device name 858 Sample Processor

Device type 858.0020 Professional Sample Processor

Target

Tower 1

Lift position Work position mm

Parameters

Lift rate 25 mm/s

CALL Transfer solution

Call text	Track name	Sample type	Condition
Peristaltic pump ON	Peristaltic pump	off Sample	off

TRACK Peristaltic pump

Return immediately off

Delete old data off

PUMP Peristaltic pump ON

Device

Device name 858 Sample Processor

Device type 858.0020 Professional Sample Processor

Pumps

Tower 1

Pump(s) Peristaltic

Rate 5

Action

Switch on off

Switch off off

Duration on

Time 150 s

TRACK Go to rinsing position

Return immediately off

Delete old data off

CALL Empty vessel

Call text	Track name	Sample type	Condition
Drain	Drain	off Sample	off

Swing rinse needle



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

SWING

Device

Device name 858 Sample Processor
Device type 858.0020 Professional Sample Processor
Target
Tower 1
Swing External position
Number 1
Parameters
Swing rate 55 °/s

LIFT

Needle → Rinse position

Device

Device name 858 Sample Processor
Device type 858.0020 Professional Sample Processor
Target
Tower 1
Lift position Work position mm
Parameters
Lift rate 25 mm/s

LQH

Fill RinsingStation

General/Hardware

Dosing device

Dosing unit WashStation

Parameters

Function Dosing
Port 1
Volume 7 mL
Rate maximum mL/min

CALL

CALL RinsingStation

Call text	Track name	Sample type	Condition
Refill RinsingStation	RinsingStation	off Sample	off

STIR & PURGE

Stirrer ON

Device

Device name 884_1
Device type 884 Professional VA

Stir

Stirring rate 3000 min⁻¹
Switch on on
Switch off off
Duration off

Purge



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

Switch on off
Switch off on
Duration off

CALL Transfer rinsing solution

Call text	Track name	Sample type	Condition
Peristaltic pump ON	Peristaltic pump off	Sample	off

LIFT Lift needle

Device

Device name 858 Sample Processor
Device type 858.0020 Professional Sample Processor

Target

Tower 1
Lift position Home position mm

Parameters

Lift rate 25 mm/s

PUMP Dose rinse EtOH

Device

Device name 858 Sample Processor
Device type 858.0020 Professional Sample Processor

Pumps

Tower 1
Pump(s) 2

Action

Switch on off
Switch off off
Duration on
Time 3 s

PUMP Empty tube

Device

Device name 858 Sample Processor
Device type 858.0020 Professional Sample Processor

Pumps

Tower 1
Pump(s) Peristaltic
Rate 5

Action

Switch on off
Switch off off
Duration on
Time 30 s

RinsingStation



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

TRACK	Return immediately	on
	Delete old data	off
WAIT	Wait 40 s	
	Wait	
	Stop track and waiting for [Continue]	off
	Stop all tracks and waiting for [Continue]	off
	Waiting time	on
	Time	40
	Unit	s
	Message	
	Record message	off
	Message by e-mail	off
	Acoustic signal	off
LQH	Refill RinsingStation	
	General/Hardware	
	Dosing device	
	Dosing unit	WashStation
	Parameters	
	Function	Dosing
	Port	1
	Volume	3 mL
	Rate	4 mL/min
TRACK	Rinsing cycle	
	Return immediately	off
	Delete old data	off
STIR & PURGE	Stir rinse ON	
	Device	
	Device name	884_1
	Device type	884 Professional VA
	Stir	
	Stirring rate	2000 min ⁻¹
	Switch on	on
	Switch off	off
	Duration	off
	Purge	
	Switch on	off
	Switch off	on
	Duration	off
LOOP	LOOP Rinsing	
	Stop criteria	
	Max. run number	on
	Max. run number	2



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

Maximum run time off
Signal assessment for DT off
Condition off

CALL CALL Drain

Call text	Track name	Sample type	Condition
Drain	Drain	off Sample	off

CALL CALL Rinse

Call text	Track name	Sample type	Condition
Rinse	Rinse	off Sample	off

TRACK Drain

Return immediately off
Delete old data off

PUMP Draining time

Device
Device name 858 Sample Processor
Device type 858.0010 Professional Sample Processor
Pumps
Tower 1
Pump(s) 1
Action
Switch on off
Switch off off
Duration on
Time 15 s

TRACK Rinse

Return immediately off
Delete old data off

PUMP Rinsing time

Device
Device name 858 Sample Processor
Device type 858.0010 Professional Sample Processor
Pumps
Tower 1
Pump(s) 2
Action
Switch on off
Switch off off
Duration on
Time 5 s



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

TRACK	Fill DU neutral Return immediately off Delete old data off
LQH	DU Electrolyte neutral General/Hardware Dosing device Dosing unit Electrolyte neutral Parameters Function Fill Port 2 Rate maximum mL/min
TRACK	Fill DU alkaline Return immediately off Delete old data off
LQH	DU Electrolyte alkaline General/Hardware Dosing device Dosing unit Electrolyte alkaline Parameters Function Fill Port 2 Rate maximum mL/min
TRACK	Fill DU WashStation Return immediately off Delete old data off
LQH	DU WashStation General/Hardware Dosing device Dosing unit WashStation Parameters Function Fill Port 2 Rate maximum mL/min
TRACK	Shut off Return immediately off Delete old data off
PUMP	PUMPS Rinse OFF Device Device name 858 Sample Processor Device type 858.0020 Professional Sample Processor Pumps Tower 1



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

Pump(s) 1+2
Action
Switch on off
Switch off on
Duration off
Time 8.0 s

PUMP

PUMP Peri OFF

Device
Device name 858 Sample Processor
Device type 858.0020 Professional Sample Processor
Pumps
Tower 1
Pump(s) Peristaltic
Rate 10
Action
Switch on off
Switch off on
Duration off
Time 300.0 s

STIR

STIR OFF

Device
Device name 858 Sample Processor
Device type 858.0020 Professional Sample Processor
Stirrer
Stirrer 1
Stirrer type unknown
Stirring rate 11
Action
Switch on off
Switch off on
Duration off
Time 10.0 s

**STIR &
PURGE**

STIR & PURGE OFF

Device
Device name 884_1
Device type 884 Professional VA
Stir
Stirring rate 2000 min⁻¹
Switch on off
Switch off on
Duration off
Purge
Switch on off
Switch off on



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

Duration off

EXIT **Exit track**

CALL **Exit shut off**

Call text	Track name	Sample type		Condition
Exit shut off	Shut off	off	Sample	off

CALL **Call fill DU**

Call text	Track name	Sample type		Condition
Fill DU neutral	Fill DU neutral	off	Sample	off
Fill DU alkaline	Fill DU alkaline	off	Sample	off
Fill DU WashStation	Fill DU WashStation	off	Sample	off

LIFT **Needle → Home position**

Device
Device name 858 Sample Processor
Device type 858.0020 Professional Sample Processor
Target
Tower 1
Lift position Home position mm
Parameters
Lift rate 25 mm/s

ERROR **Error track**

CALL **Error shut off**

Call text	Track name	Sample type		Condition
Error shut off	Shut off	off	Sample	off

Evaluation parameters

General

Sweep alkaline

Data processing

Smoothing 5

Reversed peaks off

Curve evaluation

Fixed point evaluation off

Sweep neutral



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

Data processing

Smoothing 5
Reversed peaks off

Curve evaluation

Fixed point evaluation off

Substances

Sweep alkaline

Substances - Recognition

Substance	Active	Characteristic potential	Tolerance	Min. width	Max. width	Min. measured quantity
Phenol	on	0.3 V	0.2 V	0.1 V	0.7 V	10 nA

Sweep alkaline

Substances - Baseline

Substance	Baseline type	Start base point	End base point
Phenol	Linear	Automatically	Automatically

Sweep neutral

Substances - Recognition

Substance	Active	Characteristic potential	Tolerance	Min. width	Max. width	Min. measured quantity
Amine	on	0.7 V	0.2 V	0.1 V	0.7 V	10 nA

Sweep neutral

Substances - Baseline

Substance	Baseline type	Start base point	End base point
Amine	Linear	Automatically	Automatically

Standards

Standards

Name	Standard
Amine	1000 mL/L
Phenol	1000 mL/L

Calibration

General

Calibration method External calibration
Blank value correction with evaluation off
quantity

Sweep alkaline

Calibration curves

Substance	Evaluation quantity	Curve type	Weighting
Phenol	Area	Linear regression through 0	on

Sweep neutral Calibration curves

Substance	Evaluation quantity	Curve type	Weighting
Amine	Area	Linear regression through 0	on

Results

Substance concentration in the sample

Sweep alkaline

Substance	Decimal places	Assignment
Phenol	2	none

Substance concentration in the sample

Sweep neutral

Substance	Decimal places	Assignment
Amine	2	none

Results

Additional results

Result	Places	Prefix	Unit
Peak potential	3		V
Height	2	#	A
RSD of the heights of all replications	1		%
Measured value	2	#	A
RSD of the measured values of all replications	1		%
Area	3	μ	C
RSD of the areas of all replications	1		%
Start base point	3		
End base point	3		
Standardized area	3		
Standardized height	3		
Standardized measured value	3		
Total volume	3	#	L
Zero-order coefficient	3		



License ID 1955606
Client ID 001-NBK-02363
User Metrohm

Program version viva 3.0 - 69
2024-06-26 13:44:16 UTC+2

Result	Places	Prefix	Unit
First-order coefficient	3		
Second-order coefficient	3		
Fourth-order coefficient	3		
Coefficient of determination	5		
Substance concentration in measuring vessel	2	#	
RSD of the substance concentration in measuring vessel	1		%
Amount of substance	2	#	
RSD of the substance concentration in the sample	1		%
Effective addition volume of the standard solution for the evaluation ratio	2	#	L
RSD of the effective addition volume of the standard solution for the evaluation ratio	1		%
Calibration factor DT	2	#	
RSD of the calibration factor DT	1		%
Effective addition volume of the sample solution for the evaluation ratio	2	#	L
RSD of the effective addition volume of the sample solution for the evaluation ratio	1		%

User-defined results

Result type	Result name	Formula	Unit	Decimal places	Assignment	Description
Single result	Remaining Phenol	= 'RS.Sweep alkaline.Phenol. CONC' / 'ED. Standards.Standard. Phenol.CONC' * 100	%	1	RS01	
Single result	Remaining Amine	= 'RS.Sweep neutral.% Amine.CONC' / 'ED.Standards. Standard.Amine. CONC' * 100		1	RS02	
Single result	Area Phenol mean	= 'RS.Sweep alkaline.Phenol.VAR {1}.AREA.MNV' * 1000000	µC	3	RS04	
Single result	Area Phenol RSD	= 'RS.Sweep alkaline.Phenol.VAR {1}.AREA.RSD'	%	1	RS05	
Single result	Area Amine mean	= 'RS.Sweep neutral.µC Amine.VAR{1}. AREA.MNV' * 1000000		3	RS07	
Single result	Area Amine RSD	= 'RS.Sweep neutral.% Amine.VAR{1}. AREA.RSD'		1	RS08	

Database

Name database viva