



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:26:06 UTC+1

## Method parameters

Method . . . . . AB445 - Maintenance - Cleaning measurement  
Method saving date . . . . . 2024-03-27 09:19:11 UTC+1  
Method version . . . . . 1  
Method group . . . . . Main group  
Method status . . . . . original  
Method saved by (full name) . . . . . Metrohm International Headquarters  
Method saved by (short name) . . . . . Metrohm

## START

### Main track

#### General

Workplace view

Current view . . . . . on

Track view for live window

Live display 1 . . . . . Main track

Live display 2 . . . . . Main track

Electrode check . . . . . off

#### Application note

- Specify the number of repeated cleaning measurements in the «Run» window under «Cleaning cycles». Usually 3-5 cycles are sufficient. Note! Electrolyte alkaline is used for the electrochemical cleaning of the electrodes.

#### Sample data variables

Name	Type	Assignment	Fixed value	Comment	Monitoring
ID1	Text	ID1	Cleaning	Sample identification 1	off
ID2	Text	ID2		Sample identification 2	off
Cleaning cycles	Number	ID3		Sample identification 3	off

Name . . . . . ID1

Type . . . . . Text

Assignment . . . . . on, . . . . . ID1

Fixed value . . . . . on, . . . . . Cleaning measurement

Check at start . . . . . on

Comment . . . . . Sample identification 1

Name . . . . . ID2

Type . . . . . Text

Assignment . . . . . on, . . . . . ID2

Fixed value . . . . . off, . . . . .

Check at start . . . . . on

Comment . . . . . Sample identification 2

Name . . . . . Cleaning cycles

## Measurement for electrochemical cleaning of electrodes - fully automated system

Call text	Track name	Sample type	Condition
Drain	Drain	off Sample	off

Function . . . . .	Dosing
Port . . . . .	3
Volume . . . . .	5 mL
Rate . . . . .	30 mL/min

Call text	Clean electrode
Track name	VA track
Condition	off

## Page 2 of 6



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:26:06 UTC+1

Call text	Track name	Sample type	Condition
Rinsing cycle	Rinsing cycle	off Sample	off

#### VA TRACK VA track

Return immediately . . . . . off

#### LOOP Replications

##### Stop criteria

Max. run number . . . . . on  
Max. run number . . . . . = 'SD.Cleaning cycles'  
Maximum run time . . . . . off  
Signal assessment for DT . . . . . off  
Condition . . . . . off

#### DP Sweep

##### General/Hardware

##### Device

Device name . . . . . 884\_1  
Device type . . . . . 884 Professional VA

##### Sensors/Electrodes

Working electrode . . . . . RDE  
Sensor type . . . . . RDE/SSE  
Reference electrode . . . . . Reference electrode  
Auxiliary electrode . . . . . Auxiliary electrode  
Electrode check . . . . . off

##### Stirrer

Stirring rate . . . . . 2000 min<sup>-1</sup>  
Hydrodynamic measurement . . . . . off

##### Pretreatment

Stirring time . . . . . 5.0 s

##### Cyclovoltammetric pretreatment

Start potential . . . . . -1.2 V  
Vertex potential . . . . . -0.100 V  
Sweep rate . . . . . 1 V/s  
Cycles . . . . . off  
Duration . . . . . - s

##### Potentiostatic pretreatment

Potential 1 . . . . . off V  
Waiting time 1 . . . . . 0.0 s  
Potential 2 . . . . . off V  
Waiting time 2 . . . . . 0.0 s  
Potential 3 . . . . . off V  
Waiting time 3 . . . . . 0.0 s  
Potential 4 . . . . . off V  
Waiting time 4 . . . . . 0.0 s



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:26:06 UTC+1

Potential 5 . . . . . off V  
Waiting time 5 . . . . . 0.0 s  
Equilibration time . . . . . 5.0 s

#### Sweep

Start potential . . . . . -0.1 V  
End potential . . . . . 1.0 V  
Potential step . . . . . 0.006 V  
Potential step time . . . . . 0.1 s  
Sweep rate . . . . . 0.060 V/s  
Pulse amplitude . . . . . 0.05 V  
Pulse time . . . . . 0.04 s  
Measuring time . . . . . 0.02 s  
Sweep duration . . . . . 18.33 s

#### Post-treatment

##### Cleaning

Cleaning potential . . . . . off V  
Cleaning time . . . . . 0.0 s

##### Standby potential

Standby potential . . . . . off V

#### Potentiostat

##### Current measuring range

Highest range . . . . . 2 mA  
Lowest range . . . . . 20 µA  
Automatically select optimum current measuring range . . . . . on

#### TRACK

##### Rinsing cycle

Return immediately . . . . . off  
Delete old data . . . . . off

#### LOOP

##### LOOP Rinsing

##### Stop criteria

Max. run number . . . . . on  
Max. run number . . . . . 2  
Maximum run time . . . . . off  
Signal assessment for DT . . . . . off  
Condition . . . . . off

#### CALL

##### CALL Drain

Call text	Track name	Sample type	Condition
Drain	Drain	off Sample	off

#### CALL

##### CALL Rinse



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:26:06 UTC+1

Call text	Track name	Sample type	Condition
Rinse	Rinse	off Sample	off

**TRACK**

**Drain**

Return immediately . . . . . off  
Delete old data . . . . . off

**PUMP**

**PUMP Drain**

Device

Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor

Pumps

Tower . . . . . 1  
Pump(s) . . . . . 1

Action

Switch on . . . . . off  
Switch off . . . . . off  
Duration . . . . . on  
Time . . . . . 15 s

**TRACK**

**Rinse**

Return immediately . . . . . off  
Delete old data . . . . . off

**PUMP**

**PUMP Rinse**

Device

Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor

Pumps

Tower . . . . . 1  
Pump(s) . . . . . 2

Action

Switch on . . . . . off  
Switch off . . . . . off  
Duration . . . . . on  
Time . . . . . 4 s

**TRACK**

**Shut off**

Return immediately . . . . . off  
Delete old data . . . . . off

**PUMP**

**PUMPS OFF**

Device

Device name . . . . . 858 Sample Processor  
Device type . . . . . 858.0020 Professional Sample Processor

Pumps



License ID 1955606  
Client ID 001-nbk-02849  
User Metrohm

Program version viva 3.0 - 69  
2024-03-27 09:26:06 UTC+1

Tower . . . . . 1  
Pump(s) . . . . . 1+2  
Action  
Switch on . . . . . off  
Switch off . . . . . on  
Duration . . . . . off  
Time . . . . . 8.0 s

**STIR &  
PURGE**

**STIR & PURGE OFF**

Device  
Device name . . . . . 884\_1  
Device type . . . . . 884 Professional VA  
Stir  
Stirring rate . . . . . 2000 min<sup>-1</sup>  
Switch on . . . . . off  
Switch off . . . . . on  
Duration . . . . . off  
Purge  
Switch on . . . . . off  
Switch off . . . . . on  
Duration . . . . . off

**LQH**

**Fill DU electrolyte**

**General/Hardware**

Dosing device  
Dosing unit . . . . . Electrolyte alkaline  
**Parameters**  
Function . . . . . Fill  
Port . . . . . 2  
Rate . . . . . maximum mL/min

**EXIT**

**Exit track**

**CALL**

**Exit shut off**

Call text	Track name	Sample type	Condition
Exit shut off	Shut off	off Sample	off

**ERROR**

**Error track**

**CALL**

**Error shut off**

Call text	Track name	Sample type	Condition
Error shut off	Shut off	off Sample	off