



## Compatibility of *Vision* versions with Windows™ operating systems

	Windows 2000	Windows XP (SP1, SP2, SP3)	Windows Vista	Windows 7	Windows 8	Windows 10
<b>Vision</b>						
4.0.1	✓	✓	✓	✓	✓*	–
<b>Vision Multi</b>						
4.0.1	✓	✓	✓	✓	✓*	–
<b>Vision Pharma</b>						
4.0.1	✓	✓	✓	✓	✓*	–
<b>Vision</b>						
4.0.2	✓	✓	✓	✓	✓*	–
<b>Vision Multi</b>						
4.0.2	✓	✓	✓	✓	✓*	–
<b>Vision Pharma</b>						
4.0.2	✓	✓	✓	✓	✓*	–
<b>Vision</b>						
4.0.3	✓	✓	✓	✓	✓	–
<b>Vision Multi</b>						
4.0.3	✓	✓	✓	✓	✓	–
<b>Vision Pharma</b>						
4.0.3	✓	✓	✓	✓	✓	–
<b>Vision</b>						
4.1–4.1.1.156	–	–	✓	✓	✓	✓
<b>Vision Multi</b>						
4.1–4.1.1.156	–	–	✓	✓	✓	✓
<b>Vision Pharma</b>						
4.1–4.1.1.156	–	–	✓	✓	✓	✓
<b>Vision Multi Pharma</b>						
4.1–4.1.1.156	–	–	✓	✓	✓	✓
<b>Vision</b>						
4.1.1.189 or higher	–	–	–	✓	✓	✓
<b>Vision Multi</b>						
4.1.1.189 or higher	–	–	–	✓	✓	✓
<b>Vision Pharma</b>						
4.1.1.189 or higher	–	–	–	✓	✓	✓
<b>Vision Multi Pharma</b>						
4.1.1.189 or higher	–	–	–	✓	✓	✓

\* only 64-bit operating systems are supported