



# Compatibility of *tiamo*<sup>TM</sup> versions with Metrohm devices

Device	<i>tiamo</i> <sup>TM</sup> 1.0	<i>tiamo</i> <sup>TM</sup> 1.1	<i>tiamo</i> <sup>TM</sup> 1.2	<i>tiamo</i> <sup>TM</sup> 1.2.1	<i>tiamo</i> <sup>TM</sup> 1.3	<i>tiamo</i> <sup>TM</sup> 2.0	<i>tiamo</i> <sup>TM</sup> 2.1	<i>tiamo</i> <sup>TM</sup> 2.2	<i>tiamo</i> <sup>TM</sup> 2.3	<i>tiamo</i> <sup>TM</sup> 2.4	<i>tiamo</i> <sup>TM</sup> 2.5	<i>tiamo</i> <sup>TM</sup> 3.0
<b>Titrimos</b>	702 SM Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	716 DMS Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	718 STAT Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	719 S Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	720 KFS Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	721 Net Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	736 GP Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	751 GPD Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	758 KFD Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	784 KFP Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	785 DMP Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	794 Basic Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	795 KFT Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	798 MPT Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	799 GPT Titrimo	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
	<b>Titrandos</b>	808	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
809		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
835		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
836		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
841		✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
842		-	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
851		-	-	-	-	-	-	✓	✓	✓	✓	✓
852		-	-	-	-	-	-	✓	✓	✓	✓	✓
855		-	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
857		-	-	✓	✓	✓	✓	✓	✓	✓	✓	✓
888		-	-	-	-	✓	✓	✓	✓	✓	✓	✓
890		-	-	-	-	✓	✓	✓	✓	✓	✓	✓
901		-	-	-	-	-	✓	✓	✓	✓	✓	✓
902		-	-	-	-	-	✓	✓	✓	✓	✓	✓
904		-	-	-	-	-	✓	✓	✓	✓	✓	✓
905		-	-	-	-	-	✓	✓	✓	✓	✓	✓
906		-	-	-	-	-	✓	✓	✓	✓	✓	✓
907	-	-	-	-	-	✓	✓	✓	✓	✓	✓	

Device	<i>tiamo</i> <sup>™</sup> 1.0	<i>tiamo</i> <sup>™</sup> 1.1	<i>tiamo</i> <sup>™</sup> 1.2	<i>tiamo</i> <sup>™</sup> 1.2.1	<i>tiamo</i> <sup>™</sup> 1.3	<i>tiamo</i> <sup>™</sup> 2.0	<i>tiamo</i> <sup>™</sup> 2.1	<i>tiamo</i> <sup>™</sup> 2.2	<i>tiamo</i> <sup>™</sup> 2.3	<i>tiamo</i> <sup>™</sup> 2.4	<i>tiamo</i> <sup>™</sup> 2.5	<i>tiamo</i> <sup>™</sup> 3.0
<b>Coulometers</b>												
756 KF Coulometer	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
831 KF Coulometer	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>Sample Changers</b>												
730 Sample changer	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
774 Oven Sample Changer	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
778 Sample Processor	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
789 Robotic SP XL	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
814 USB Sample Processor	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
815 Robotic USB Sample Processor	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
855 Robotic Titrosampler	–	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
864 Robotic Balance Sample Processor	–	–	–	–	–	✓	✓	✓	✓	✓	✓	✓
874 Oven Sample Processor	–	–	–	–	–	✓	✓	✓	✓	✓	✓	✓
<b>Meters</b>												
712 Conductometer	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
856 Conductivity Module	–	–	–	–	–	✓	✓	✓	✓	✓	✓	✓
<b>Miscellaneous</b>												
846 Dosing Interface	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
854 iConnect	–	–	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓
859 Titrotherm	–	–	–	–	–	–	–	✓	✓	✓	✓	✓
867 pH module	–	–	–	–	–	✓	✓	✓	✓	✓	✓	✓
Avantes spectrometer	–	–	–	–	–	✓	✓	✓	✓	✓	✓	✓
Barcode readers	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓	✓