

1 Purpose

This document describes the fixed bugs and problems as well as improvements in the following new firmware version in comparison with the previous version.

- New version: 8.0.0.91
- Old version: 0.1.1.49

2 New features

Topic	Description
New hardware	The "ID Kit" and SERS test kit were added.
Library management	Library management can be enabled and disabled.
Reprocess samples	Reprocess samples by using operating procedure parameters against enabled libraries.
Added languages	Czech, Hungarian and Turkish were added.

3 Improvements

Topic	Description
Auto-integration toggle	An auto-integration toggle when building an operating procedure was added.
Stability	Various improvements for stability.
Translations	The translations were updated.

4 Fixed bugs

Topic	Description
Mixture match (reprocess)	Mixture matching is forced when reprocessing sample data now.
ST Japan sub-library	The ST Japan sup-libraries can now be enabled independently.

