

1 Purpose

This document describes the fixed bugs and problems as well as improvements in the following new firmware version in comparison with the previous version.

- New version: 8.0.1.136
- Old version: 8.0.0.98

2 New features

Topic	Description
SERS Tip	SERS Tip can be used with ID Kit operating procedure and with other libraries for R&D purposes.
New library	New 19,000+ entries library is available: Mira Lib Bio-Rad™ Full Raman Library Collection.
Result screen	New Mixture Match result screens.
Intelligent Universal Attachment	New intelligent Universal Attachment (iUA) is implemented in firmware.
New functions	Guided user interface for iUA. iUA positions for surface, bag and bottle measurements.
Content ID operating procedure	Content ID operating procedure is available with iUA, results are separated to show the container and the content.
Hit Quality Index	HQI (Hit Quality Index) is available on the instrument.
Foreign language support	Added foreign language support to library names.
Spectral representation	Spectral representation (spectral quantification) is shown in the Mixture Matching result screen.
Renaming spectrum	Rename spectrum on the instrument.
Global Harmonization System	Global Harmonization System placards are now visible on the instrument.

