

1 Purpose

This document describes new functions, the changes, improvements and fixed bugs and problems in the following firmware version:

- 57.1008.0010
- in comparison with the firmware version:
- 57.1008.0009 (first released version)

2 New functions

The following new functions have been integrated:

- Changing the password for the **Expert** dialog type
- Factory reset
- User interface: 10 new languages added:
 - French
 - Italian
 - Japanese
 - Korean
 - Polish
 - Portuguese
 - Russian
 - Spanish
 - Thai
 - Vietnamese

3 Improvements

The following improvements were made:

- Firmware optimizations for the final instrument check in the production.

4 Fixed bugs and problems

The following bugs and problems were fixed:

- OHMTI-687 – User interface adjustment
- OHMTI-753 – User interface adjustment
- Various further user interface adjustments.
- OHMTI-763 – PDF reports: Serial numbers are added to the measuring point list.
- OHMTI-909 – Live status: In some cases, not the entire titration curve was displayed.

5 Known bugs and problems

The following bugs or problems are known:

- none