

1 Purpose

This document describes new functions, the changes, improvements and fixed bugs and problems in the following firmware version:

- 57.1008.0013 for OS version 005

in comparison with the firmware version:

- 57.1008.0011 for OS version 005

2 New functions

The following new functions have been integrated:

- none

3 Improvements

The following improvements were made:

- Optimizations for the use of special characters in method names
- General improvements

4 Fixed bugs and problems

The following bugs and problems were fixed:

- TITR-16272 Export of methods with special characters
- TITR-16601 Screen freezes and buret does not refill (in SET mode)
- TITR-16569 Incorrect calculation variable

5 Known bugs and problems

The following bugs or problems are known:

- Problems with saving or loading methods / sensors / solutions while a sample series is being processed.
- If the equivalence point criterion is set to «OFF», then the fixed point cannot be used for calculations.