

1 Purpose

This document describes new functions, changes, improvements and fixed bugs in the following firmware version:

- 57.1028.0003 for OS version 001

In comparison with the firmware version:

- 57.1028.0002 for OS version 001

2 New functions

The following new functions have been integrated:

- Switching off the system is now divided into "Shutting down the system" (standby) and "Switching off the system".

3 Improvements

The following improvements were made:

- Various write errors were fixed.
- Better recognition of devices when connected to a USB connector.
- In the method templates, the value **Off** is defined as standard for "Statistics".
- Additional columns in the sample list for method and result 1.
- During a CN firmware update, it is no longer possible for the user to enter commands via the command line.
- The time required to change the method on the start page has been reduced.
- The message that appeared when opening the system menu has been removed.
- When working with a sample list, the method of the previous determination is automatically adopted when a new determination is inserted.
- The screens are loaded at startup. This means that content is displayed more quickly when switching between screens.

4 Fixed bugs

The following bugs were fixed:

