

1 Instruments

The firmware version below can be used with the following instruments:

- 850 Professional IC
- 881 Compact IC pro
- 882 Compact IC plus
- 883 Basic IC plus
- 886 Professional Thermostat / Reactor
- 896 Professional Detector
- 920 Absorber Module
- Eco IC

2 Purpose

This document describes the new features, the improvements and the fixed bugs in the following new firmware version:

- 5.850.0114

In comparison with the firmware version:

- 5.850.0113

3 New features

Topic	Description
2.850.9020 IC Conductivity Detector MB	The 2.850.9020 IC Conductivity Detector MB is now supported.

4 Improvements

Topic	Description
None	None

5 Fixed bugs

Topic	Description
None	None

6 Known bugs

Topic	Description
None	None

7 Conformity

The firmware version 5.850.0114 does not contain any modifications that affect conformity of the instruments mentioned above regarding 21 CFR Part 11 and GAMP.