



Eco KF Titrator

2.1027.0010

The new Eco KF Titrator with integrated magnetic stirrer und touch-sensitive display is ideal for routine water content determination. Predefined methods enable a smooth and easy instrument start-up. In addition to GLP-compliant printouts on paper or as PDFs, the Eco KF Titrator also offers the option of connecting balances or of sending the determination data to a PC via PC / LIMS reports.

Below, the accessories are grouped into Scope of delivery and Optional accessories. Please keep this printout at hand for ordering replacement material. These lists may be subject to change.

Scope of delivery 2.1027.0010

Qt.	Order no.	Description
-----	-----------	-------------

1 PCS

1.1027.0010

Eco KF Titrator

The new Eco KF Titrator with integrated magnetic stirrer und touch-sensitive display is ideal for routine water content determination. Predefined methods enable a smooth and easy instrument start-up. In addition to GLP-compliant printouts on paper or as PDFs, the Eco KF Titrator also offers the option of connecting balances or of sending the determination data to a PC via PC / LIMS reports.



1 PCS

6.02712.000

Protective cover for Eco Dosimat / Eco Titrator

Protective cover for Eco Dosimat / Eco Titrator display.



1 PCS

6.1446.080

Threaded stopper / M8



1 PCS

6.1546.040

Piston tool for Eco Dosimat/Eco Titrator

Tool for removing and mounting the piston.



1 PCS

6.1602.105

Bottle attachment / GL 45 / green

For direct dosing from reagent bottles with GL 45 thread.



1 PCS

6.1619.010

Adsorber complete, for Exchange Unit



2 PCS

6.1805.090

FEP tubing / M6 / 31 cm

With light and kink protection



1 PCS

6.1805.100

FEP tubing / M6 / 40 cm

With light and kink protection



1 PCS

6.1819.020

Tube

For M6 threaded holes. For 6.1602.1XX bottle attachments



1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2016.070

Support rod / 400 mm



1 PCS

6.2023.020

Clip for SGJ 14/15

Clip for SGJ 14/16



1 PCS

6.2164.010

Power supply unit 100 - 240 V/24 V DC

Power supply unit for 915 KF Ti-Touch, 916 Ti-Touch, 917 Coulometer and all Eco Titrators / Eco Dosimat.



1 PCS

6.6081.020

USB flash drive for Eco KF Titrator

Storage medium for Eco KF Titrator



Optional accessories

Order no.	Description	
2.1029.0010	Solvent Pump Solvent Pump is used with an Eco Titrator / Eco KF Titrator / Eco Coulometer for reliably adding solvent and aspirating off titration vessel content without contact with the reagents. Device-independent use is also possible with optional 6.2164.010 power supply unit.	
2.142.0100	Custom Q3X thermal printer Compact printer with USB interface for • 900 Touch Control • 915 KF Ti-Touch • 916 Ti-Touch • 917 Coulometer • 877 / 848 Titrino plus • 865 / 876 Dosimat plus • 91X Meter (cable 6.2151.140) • Eco Dosimat / Titrator • 862 Compact Titrosampler • 870 KF Titrino plus • 899 Coulometer Paper width 60 mm (40 characters). Including 6.2151.120 USB cable.	

Compact impact printer with USB and RS232 interface for

- 900 Touch Control (+)
- 915 KF Ti-Touch (+)
- 916 Ti-Touch (+)
- 917 Coulometer (+)
- 91X Meter (+)
- Eco Dosimat/Titrator (+)
- 877/848 Titrino plus (additional 6.2151.100 required)
- 865/876 Dosimat plus (additional 6.2151.100 required)
- 862 Compact Titrosampler (additional 6.2151.100 required)
- 870 KF Titrino plus (additional 6.2151.100 required)
- 899 Coulometer (additional 6.2151.100 required)



(+) Cable included in scope of delivery

6.01406.250

50 – 150 mL titration vessel

50 – 150 mL titration vessel with screw top for use with 6014050X0 titration vessel lids.



6.02709.040

Stopper

Stopper for titration vessel lids 6014050X0.



6.02711.000 Spoon for paste

Suitable for use in Karl Fischer water content determinations of pasty substances such as honey, cosmetics, greases or various dairy products.

The spoon for paste facilitates the addition of viscous samples. It is handled like any other spatula and enables convenient, rapid and loss-free addition of pastes into the Karl Fischer titration vessel.

To dissolve viscous samples quickly and completely, we also recommend the use of a homogenizer.



6.03002.150 Cylinder unit Eco 5 mL

Cylinder unit 5 mL for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.



6.03002.210 Cylinder unit Eco 10 mL

Cylinder unit 10 mL for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.



6.03002.220 Cylinder unit Eco 20 mL

Cylinder unit 20 mL for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.



6.03002.250 Cylinder unit Eco 50 mL

Cylinder unit 50 mL for Eco Dosimat / Eco Titrator. Without dosing tubing and antidiffusion tip.



6.05610.050 KFT equipment for Eco KF Titrator

Complete accessory set for volumetric Karl Fischer titration



6.05610.060 KFT equipment for Eco KF Titrator with Polytron

Complete accessory set for volumetric Karl Fischer titration with Polytron



6.9001.050 Thermal paper for Custom Neos-N / Q3X

Thermal paper 60 mm / 25 m for Custom Neos-N / Q3X.
New: 25-year legibility with correct storage.

1 roll

