



## Transfer Tubing with holder 2 x M6 / 10 mL to Sample Processors

6.1562.130

Transfer tubing for dilution and preconcentration using a Dosino and Sample Processor.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

### Optional accessories

Order no.	Description
6.2744.200	Adapter UNF 10/32 outer / M6 inner Connecting 1/16 in. PEEK capillaries and transfer tubing (M6)

