



dPt-Titrode

6.00401.300

Electrodo de anillo de platino combinado digital para OMNIS con una membrana de vidrio de pH como electrodo de referencia.

Este electrodo no precisa mantenimiento y es apto para titulaciones redox con un valor de pH constante, por ejemplo:

- Yodometría
- Cromatometría
- Cerimetría
- Permanganometría

Este electrodo se almacena en agua destilada.
Los dTodes se pueden utilizar en OMNIS Titrator.

Partes/accesorios 6.00401.300

Qt.	Order no.	Descripción
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Accesorios opcionales

Order no.	Descripción	
6.02104.300	Electrode cable plug-in head Q / plug P, 0.55 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	
6.02104.310	Electrode cable plug-in head Q / plug P, 1.5 m Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.	