



Electrodo de anillo de platino combinado

6.0451.100

Electrodo de anillo de platino combinado con un diafragma de varilla cerámica.

Este electrodo es apto para titulaciones redox con un valor de pH variable, por ejemplo:

- Contenido de oxígeno según Winkler
- Determinación de peróxido de hidrógeno con KMnO_4
- Titulaciones de diazotación

Como electrolito de referencia y para su conservación se utiliza $c(\text{KCl}) = 3 \text{ mol/L}$.

Partes/accesorios 6.0451.100

Qt.	Order no.	Descripción
-----	-----------	-------------

1 PCS

6.1236.050

Sleeve with SGJ 14/12 mm

Sleeve with SGJ 14/12 mm, polyethylene.



1 PCS

6.2008.040

Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Accesorios opcionales

Order no.	Descripción	
6.2104.020	Electrode cable / 1 m / F For connecting electrodes with Metrohm plug-in head G to Metrohm instruments (socket F).	
6.2308.020	Electrolyte 3 mol/L KCl (250 mL) Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)	