



Tubos de calibración O₂

6.1260.000

Tubos de calibración para la calibración del Lumitrode para O₂ a una saturación de aire al 100%.

Accesorios opcionales

Order no.	Descripción
-----------	-------------

The optical sensor for measuring dissolved oxygen (DO) can be used with a 913 pH/DO meter or with a 914 pH/DO conductometer. The measuring principle of the sensor is based on luminescence quenching. The space-saving and maintenance-free sensor is suitable for DO measurement, for example in:

- Water quality control
- Wastewater industry
- Beverage production
- Fish farming

This sensor is supplied with a calibration vessel.

Where necessary, it is easy to replace the measurement cap (O₂ cap) that contains the oxygen-intensive luminophore.

