



# OMNIS Titrator Food

2.1001.2220

L'OMNIS Titrator Food vi offre il pacchetto completo per la titolazione acquosa acido/base. Il pacchetto contiene l'OMNIS Advanced Titrator con agitatore magnetico, un'unità di cilindro da 20 mL, un d-ECotrode plus per la titolazione acquosa acido/base e una licenza stand-alone del software OMNIS.

## Parti incluse 2.1001.2220

| Qt.   | Order no.   | Descrizione    |
|-------|-------------|----------------|
| 1 PCS | 6.00201.300 | dECotrode Plus |

Digital, combined pH electrode for OMNIS.  
The electrode is suitable for aqueous acid/base titrations.  
The fixed ground-joint diaphragm is insensitive to contamination.  
Reference electrolyte:  $c(\text{KCl}) = 3 \text{ mol/L}$ , storage in storage solution.  
dECotrodes can be used on OMNIS Titrators.



1 PCS

6.01600.010

OMNIS Liquid Adapter

Adapter for connecting reagent bottles with 3S technology to an OMNIS Titrator, Titration Module or Dosing Module.



---

1 PCS

6.01601.000

Bottle cap, multi-use

Reusable bottle cap with 3S technology for connecting reagent bottles to an OMNIS Titrator, Titration Module or Dosing Module.



---

1 PCS

6.02004.010

Bottle holder

Bottle holder for OMNIS Titrator, Titration Module and Dosing Module.




---

1 PCS

6.02005.000

OMNIS electrode holder

Electrode holder for one OMNIS titration system.



1 PCS

6.02005.010

Electrode Parking Station

Electrode Parking Station for electrode storage on an OMNIS titration system.



---

1 PCS

6.02100.010

Measuring module digital

Digital measurement channel for one OMNIS Titrator or Titration Module for the connection of digital "dTrodes" electrodes.



---

1 PCS

6.02104.300

Electrode cable plug-in head Q / plug P, 0.55 m

Cable for connecting a dTrobe (plug-in head Q) to the digital measuring module from OMNIS.



---

1 PCS

6.02107.000

Network cable, 1 m

Network cable (1 m) for connecting one OMNIS Sample Robot or OMNIS Titrator to an Ethernet Router, Ethernet Switch or local network



Intelligent 20 mL cylinder unit for an OMNIS Titrator, Titration Module or Dosing Module. Includes dosing tubing and antidiffusion tip.



License-based PC program for controlling complex OMNIS analysis systems.

- Sample-oriented operating procedures with numerous templates
- Graphical method editor with numerous templates
- Clearly organized and professional sample lists
- Management of instruments, work systems, sensors, solutions
- Dialog languages: German, English, Chinese



Function license "Advanced" for one OMNIS Titrator

- Liquid Handling (LQH)
- MEAS U/T/pH
- MEAS CONC (direct measurement and standard addition)
- Calibration (Cal)
- Endpoint titration without conditioning (SET)
- Monotonic equivalence point titration (MET)
- Dynamic equivalence point titration (DET)
- Titration with OMNIS Titration Module or OMNIS Dosing Module
- Parallel dosing



Enables stand-alone operation of the OMNIS software on a Windows™ computer.

Features:

- The license already includes one OMNIS instrument license.
- Must be activated via the Metrohm licensing portal.
- Not transferable to another computer.



For pistons for 807 Dosing Units or cylinder units OMNIS.



Bottle for auxiliary solutions.



Stirring Bar with magnetic core, PTFE coated, length 16 mm



1 PCS

6.1903.030

Stirring bar / 25 mm

Stirring bar with magnetic core, PTFE covering, length 25 mm.



---

1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



---

1 PCS

6.2016.050

Support rod / 300 mm



---

1 PCS

6.2307.230

Buffer solutions pH 4, 7 and 9

Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL



## Maintenance kit for pH electrodes

The kit contains:

- 50 mL cleaning solution
  - 50 mL 3M KCl solution
  - 50 mL storage solution
  - 2 Storage vessels
  - Instructions for use
- 



# Accessori opzionali

| Order no.   | Descrizione                                                                                             |                                                                                       |
|-------------|---------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------|
| 6.05609.000 | KFT equipment for OMNIS Titrator                                                                        |                                                                                       |
|             | Complete accessory set for volumetric Karl Fischer titration.                                           |    |
| 6.05609.010 | KFT equipment for OMNIS Titrator with Polytron                                                          |                                                                                       |
|             | Complete accessory set for volumetric Karl Fischer titration with Polytron.                             |    |
| 6.06001.160 | Function license KFT with conditioning                                                                  |                                                                                       |
|             | Function license for volumetric Karl Fischer titration with conditioning for an OMNIS titration system. |  |