



Ag/AgCl 12.5 cm

6.0726.100

/ 10 cm

14/15

Scope of delivery 6.0726.100

Qt.	Order no.	Description
1 PCS	6.2008.040	Storage vessel

Together with 6.2008.050 storage vessel holder. it provides a support for the electrode on 807 Dosing Units.



Optional accessories

Order no.	Description	
6.1243.030	Diaphragm for reference electrode Spare ground-joint diaphragm for 6.0726.1XX and 6.0729.1XX reference electrodes	
6.2106.020	Strand / 1 m / 2 x B Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B	
6.2106.060	Strand / 2 m / 2 x B Electrode cable with 2 x plug B (4 mm) for connecting electrodes with plug-in head B	
6.2308.020	Electrolyte 3 mol/L KCl (250 mL) Electrolyte solution $c(\text{KCl}) = 3 \text{ mol/L}$, (for Ag/AgCl reference systems)	

6.2310.010

Electrolyte KNO₃-1mol/L 250 mL

Electrolyte solution KNO₃ 1 M (reference electrolyte for combined silver electrodes)



6.2320.000

TEABr 0.4 mol/L in ethylene glycol (250 mL)

Electrolyte solution TEABr (tetraethylammonium bromide in ethylene glycol), c(TEABr) = 0.4 mol/L

