



804 Ti Stand mit Stativ

2.804.0040

Titrierstand und Controller für Propellerrührer 802 Stirrer. Ergibt zusammen mit dem optionalen 802 Stirrer eine Alternative zum Magnetrührer. Titrierstand inklusive Bodenplatte, Stativstange und Elektrodenhalter.

Lieferumfang 2.804.0040

Qt.	Order no.	Beschreibung
1 PCS	1.804.0010	804 Ti Stand without stand rod
		Titration stand and controller for 802 Stirrer. The 804 Ti Stand together with the optional 802 Stirrer provides an alternative to the magnetic stirrer. Without stand.




1 PCS

6.2001.060

Support stand without support rod

To attach an 801, 804, 803 Stirrer to a Titrando, Titrino plus, Dosimat plus or 856 and 867.



1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2016.070

Support rod / 400 mm



1 PCS

6.2021.020

Electrode holder

Electrode holder for 4 electrodes and 2 buret tips



1 PCS

6.2621.070

5 mm hex key for IC Sample Processors



1 PCS

6.2621.130

Hexagon key 2 mm

2 mm.



Optionales Zubehör

Order no.	Beschreibung	
6.2001.090	Base plate on both sides Double stand for attaching 2 x 801, 803, 804 to a Titrande	
6.2151.010	Cable MSB plug / socket Extension cable for MSB connections and for Touch Control - Titrande, Mini-DIN socket - mini-DIN plug	