



Bürettenspitze / Schlauchlänge 0.4 m

6.1541.010

Mit Antidiffusionsventil.

Lieferumfang 6.1541.010

Qt.	Order no.	Beschreibung
1 PCS	6.1446.030	Ball stopper / B-14/(15) For the movable holding of buret tips in titration vessels



1 PCS

6.1543.200

Titration tip / M6 thread

Titration tip with antidiffusion valve.



1 PCS

6.1805.100

FEP tubing / M6 / 40 cm

With light and kink protection

