



3/250 mm PTFE tubing with screw nipple

6.01805.130

Aspiration tubing with screw nipple to connect the reagent bottle to the Siphon Beaker.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Scope of delivery 6.01805.130

Qt.	Order no.	Description
1 PCS	6.1819.030	Tube For M8 threaded holes. For 6.1602.1XX bottle attachments and 6.1608.0XX bottles

