



Gold Electrode Tip

6.1204.140

Metal electrode tip, gold (Au). Electrode disk diameter 2 mm. M3 thread. Temperature range: 0 ... 40 °C.

Below, the accessories are grouped into Scope of delivery and Optional accessories.

Please keep this printout at hand for ordering replacement material.

These lists may be subject to change.

Optional accessories

Order no.	Description	
6.2103.110	Contact pin M3, 2 mm For 6.1241.060 Electrode shaft and 6.1204.1XX Electrode tips. With M3 thread. Used with VA equipment sets and for ELCD applications.	

6.2802.000

Polishing set for solid-state electrodes

Polishing set with 1 polishing cloth and about 2 g aluminum oxide powder (grain size 0.3 μm)

