




Juego de accesorios CI: Waste Collector

6.5336.000

Waste Collector para un máximo de cinco capilares de residuos.

Partes/accesorios 6.5336.000

Qt.	Order no.	Descripción
1 PCS	6.01816.010	Silicone tubing, 8/12 mm, 2 m
		Drainage tubing



1 PCS

6.1602.210

Attachment 5 x UNF 10/32 for «Waste Collector»

Attachment for holding up to 5 waste capillaries.



1 PCS

6.2057.220

Holder to Waste Collector

Holder for fastening the Waste Collector to the IC instrument.



1 PCS

6.2765.000

Container with GL 45 drainage to the "Waste Collector"

Container with drainage to the "Waste Collector".

