



843 Membrane Pump Station per sistemi Professional VA/CVS

2.843.0240

La 843 Membrane Pump Station consente il lavaggio automatico e lo svuotamento della cella di misura in sistemi Professional VA/CVS parzialmente o completamente automatici. Dispone di due pompe a membrana integrate che possono essere comandate tramite segnale remoto.

La variante per VA e CVS dispone di tutti gli accessori per il lavaggio automatico e lo svuotamento della cella di misura.

Parti incluse 2.843.0240

Qt.	Order no.	Descrizione
-----	-----------	-------------

1 PCS

1.843.0020

843 Pump Station (membrane)

The 843 Pump Station (membrane) has two built-in membrane pumps. These can be controlled directly via the interface using remote signals or manually by pressing a button.



4 PCS

6.1446.040

Threaded stopper / M6



1 PCS

6.1602.115

Bottle attachment / S 40 / green

For direct dosing from reagent bottles with S 40 thread (Merck...).



2 PCS

6.1621.000

Container 10 L

As rinsing or waste container in automated systems.



4 PCS

6.1805.530

FEP tubing / M6 / 2 m

With light and kink protection



2 PCS

6.1819.010

PTFE tube

As gas inlet in washing bottle of VA Stand



1 PCS

6.1819.110

FEP aspiration tubing to canister

Aspiration tubing to canister attachment with M6 connector. For the production of eluents.



4 PCS

6.1820.010

Adapter for 1/8" gas line/M6 thread

For connecting a 1/8" gas line to an M6 thread.



1 PCS

6.1828.020

Five-way M6 tubing connector for canister

For connecting up to 5 tubings



1 PCS

6.1829.070

Tubing set for Professional VA/CVS instruments

Tubing set for the measuring head of Professional VA/CVS instruments. Comprised of 4 FEP tubings with M6 connector, tapering tip and protected against kink (2 x 82 mm long und 2 x 110 mm long).



1 PCS

6.2053.000

Cable clip

Cable clip for fastening cables and tubes



1 PCS

6.2621.000

Adjustable wrench

Max. opening: 20 mm. For IC instruments



Accessori opzionali

Order no.	Descrizione	
6.2141.200	Cable 732/819(remote)-733/820(fill/inject) Cable IC Detector (remote) - Advanced Separation Center (fill/inject)	
6.2141.300	Remote cable 843 Pump Station to Sample Processor Remote cable for the direct connection of the 843 Pump Station to the pump connectors of the Sample Processors	
6.2148.010	Remote Box MSB Additional remote interface for the connection of devices that can be controlled via remote lines. With permanently attached cable.	