



Metrosep Amino Acids 1 - 100/4,0

6.4001.410

Colonna di separazione per la determinazione di amminoacidi tramite rilevamento UV/VIS secondo PCR con ninidrina a 120 °C. La Metrosep Amino Acids 1 - 100/4,0 è la colonna di separazione standard per amminoacidi. La colonna è basata su un materiale di polistirene /divinilbenzene solfonato. La determinazione degli amminoacidi avviene mediante rilevamento fotometrico dopo una derivatizzazione post-colonna con ninidrina.

La Metrosep Amino Acids 1 - 100/4,0 consente la separazione, nella ricerca e nella pratica, di fino a 44 amminoacidi inclusi tutti gli amminoacidi presenti in natura.

Parti incluse 6.4001.410

Qt.	Order no.	Descrizione
2 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Accessori opzionali

Order no.	Descrizione	
6.1011.030	Metrosep RP 2 Guard/3.5	
	The Metrosep RP 2 Guard/3.5 is a guard column for universal use. It reliably protects the analytical separation column against contaminations. It reliably protects the analytical separation column against contamination, removing the smallest particles, traces of iron oxide and bacteria. The Metrosep RP 2 Guard/3.5 helps to reduce costs; its filter disk can be replaced very easily.	
6.1011.040	Metrosep RP 3 Guard HC/4.0	
	The Metrosep RP 3 Guard HC/4.0 is a guard column for universal use. It reliably protects the analytical separation column against contamination, securely removing lipophilic organic contamination, the smallest particles, traces of iron oxide and bacteria. The guard column is based on a polymer material and thanks to its larger pack volume, has a significantly higher capacity than the Metrosep RP 2Guard/3.5. It can be used throughout the entire pH range.	
6.1011.130	Spare filter for RP 2 Guard/3.5, 10 pieces	