



12 x 250 mL

6.2041.310

810 Sample Processor 814 USB Sample Processor 855 Robotic Titrosampler "Light" 12
6.1432.320 6.1432.323 6.1453.220 6.1453.250

(PVC)

Optional accessories

Order no.	Description
6.1432.320	Sample beaker, clear glass, 250 mL, 1 piece
	Sample beaker, clear glass, 250 mL, suitable for sample racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820, 6.2068.020 and 6.2068.030, 1 piece



6.1432.323

Sample beaker, amber glass, 250 mL, 1 piece

Sample beaker, amber glass, 250 mL, suitable for sample racks:
6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510, 6.2041.820,
6.2068.020, and 6.2068.030, 1 piece



6.1453.220

Sample beaker, plastic (PP), 200 mL, 1 piece

Sample beaker, plastic (polypropylene, PP), 200 mL, suitable for
sample racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510,
6.2041.820, 6.2068.020, and 6.2068.030, 1 piece



6.1453.250

Sample beaker, plastic (PP), 250 mL, 1 piece

Sample beaker, plastic (polypropylene, PP), 250 mL, suitable for
sample racks: 6.2041.310, 6.2041.400, 6.2041.430, 6.2041.510,
6.2041.820, 6.2068.020, and 6.2068.030, 1 piece



6.1458.010

Titration head, 6 x SGJ 14 and 3 x SGJ 9 openings.

With 6xSGJ 14 and 3x SGJ 9 openings.

